


varian

Automation with Velocity™
Intelligent Automation for
Treatment Planning Workflows

Anthony F. Waller
Product Manager




Automation

1. DICOM Automation
2. Structure Scripting
3. ACTOR & ACTIVE dosimetry
4. Developer Mode API

varian

DICOM Automation





varian

ARIA™ Sync & Save to ARIA

- Velocity system automatically synchronizes imaging from ARIA
- One-click "Save to ARIA" saves back to ARIA
- Saves all data viewed, and any related information
- Configurable for Eclipse-only sites as well


4 VARIAN CONFIDENTIAL, PROPRIETARY, DISCLOSED SOLELY FOR REGULATORY RECIPIENT ONLY



SmartQuery PACS automation

- Add patients to PACS-specific "watchlist"
- Recent patients are added by default
- Automatically pull new imaging for patients
- Configurable and customizable for specific requirements and workflows


varian



DICOM triggers Query and Import

- Set up actions for automatic behavior based on DICOM query results and/or import
- Scriptable to include pulling related objects, creating body contours, pushing data
- Allows sending messages to users as notifications
- XML for configuration

6 VARIAN CONFIDENTIAL, PROPRIETARY, DISCLOSED SOLELY FOR REGULATORY RECIPIENT ONLY






Contour Scripting Automation for structure operations

- XML Scripting language
- Margins
- Boolean
- Smoothing
- Post-Processing (hole removal, extraction, etc)
- Set or change colors and names
- Create, delete, edit

© VARIAN CONFIDENTIAL PROPRIETARY. DISCLOSED ONLY FOR SRS/DATE RECIPIENT ONLY

varian

Edit and Validate Import, Export, Validate, and View



© VARIAN CONFIDENTIAL PROPRIETARY. DISCLOSED ONLY FOR SRS/DATE RECIPIENT ONLY


varian



Adaptive Calculation and Tracking
 Navigators for use with Eclipse and ARIA

ACTOR: Commissioning	ACTOR: Plan Generation	ACTIVE Dosimetry: Accumulation
<ul style="list-style-type: none"> Used for commissioning of the system Verifies that system reproduces original dose for original plan Aids in setup of electron density curves, dose calculation models, etc. 	<ul style="list-style-type: none"> Bulk deformable for online-imaging Creates "adaptive CTs" for dose calculation Copies RT-Plans onto adaptive CTs and pushes into Eclipse Reports structure changes over time & Deformable QA 	<ul style="list-style-type: none"> Accumulates delivered doses Compares plan vs delivered DVH metrics Creates subtraction dose volume for hot/cold analysis


varian



AdaptiveCT
 Planning CT
 Hounsfield units
 with treatment
 anatomy

varian

Scalar Magnitude of Deformation Vectors (raw data)




22 VARAN CONFIDENTIAL/PROPRIETARY. DISCLOSED SOLELY FOR IMMEDIATE RECIPIENT ONLY

varian



Min/Max of Scalar Magnitude of Deformation Vectors




23 VARAN CONFIDENTIAL/PROPRIETARY. DISCLOSED SOLELY FOR IMMEDIATE RECIPIENT ONLY

varian



Box/Whiskers of Jacobian Determinants (natural log)




24 VARAN CONFIDENTIAL/PROPRIETARY. DISCLOSED SOLELY FOR IMMEDIATE RECIPIENT ONLY

varian



Architecture API + VelocityEngine



28 VARIAN CONFIDENTIAL/PROPRIETARY. DISCLOSED SOLELY FOR IMMEDIATE RECIPIENT ONLY

Languages Supported VelocityEngine packages

- **C#**
 - .NET assembly
 - .NET 4.0
 - x64 architecture
- **Python**
 - Python 2.7
 - CPython interpreter (standard)
 - Installable as Python wheel

```

static void Main(string[] args) {
    var engine = new VelocityEngine();
    engine.LoginToGrid(USERNAME, PASS, GRID_IP, GRID_PORT, GRID_DB);
    engine.LoadPatient(PATIENT_ID);
}
  
```

```

import velocity
engine = velocity.VelocityEngine()
engine.LoginToGrid(USERNAME, DB_PASS, DB_IP, DB_PORT, DB_NAME)
engine.LoadPatient(PATIENT_ID)
  
```

varian

VARIAN CONFIDENTIAL/PROPRIETARY. DISCLOSED SOLELY FOR IMMEDIATE RECIPIENT ONLY

Developer Mode API details & Restrictions

- Works with Velocity 4
- Initial release is non-clinical
 - Requires Developer Mode license
- "Script" user type
 - Each login consumes a license, released on log out
 - Workstation is single-login
 - Grid allows multiple script logins with a single user

30 VARIAN CONFIDENTIAL/PROPRIETARY. DISCLOSED SOLELY FOR IMMEDIATE RECIPIENT ONLY

varian

How to use API Uses of Developer Mode


- **Automate** – take a manual task and automate it for bulk processing, e.g. deformable registration / quality assurance
- **Replace** – replace the Velocity client application with your own custom software
- **Leverage** – take advantage of utility functions for working with DICOM data and registration objects
- **Patient Transfer** – export DICOM data to Eclipse, or transfer patients to a clinical Velocity database

31 VARIAN CONFIDENTIAL/PROPRIETARY; DISCLOSED SOLELY FOR IMMEDIATE RECIPIENT ONLY

varian

Example Script Deformable QA structure “highlighting” areas of uncertainty

Use case:
Automate deformable and create structures for review




32 VARIAN CONFIDENTIAL/PROPRIETARY; DISCLOSED SOLELY FOR IMMEDIATE RECIPIENT ONLY

varian

Example Script Patient specific “adaptive” QA of delivered dose

Use case:
Automatically calculate dose for CBCT imaging based on treatment shifts




33 VARIAN CONFIDENTIAL/PROPRIETARY; DISCLOSED SOLELY FOR IMMEDIATE RECIPIENT ONLY

varian

VarianDeveloper GitHub

Free download



34 VARIAN CONFIDENTIAL/ PROPRIETARY. DISCLOSED SOLELY FOR IMMEDIATE RECIPIENT ONLY

varian



