

SIEMENS

CARE Right

Computed Tomography, committed to the right dose
Ivo Driesser
Austin, July 22nd 2014

Restricted © Siemens AG 2014 All rights reserved. Answers for life.

SIEMENS

Is 1 mSv the right dose for every patient?


Male, 43
normal CT scan of the abdomen

Female, 54
normal CT scan of the abdomen

Female, 99
CT scan of whole body used to detect metastases

Restricted © Siemens AG 2014 All rights reserved.
Page 2 07-20-2014

SIEMENS



In this particular case, 0.9 mSv was the right dose


In this particular case, 2.4 mSv was the right dose

In this particular case, 7.9 mSv was the right dose

Restricted © Siemens AG 2014 All rights reserved.
Page 3 07-20-2014

SIEMENS

For the ideal balance between radiation dose and image quality:



Restricted © Siemens AG 2014 All rights reserved.
Page 4 07-20-2014

Courtesy of Medical General Hospital / Hannover, Israel

CARE Right
The unique Siemens approach

SIEMENS



Restricted © Siemens AG 2014 All rights reserved.
Page 5 07-20-2014

CARE Right
Committed to the Right Dose in CT


SIEMENS



Restricted © Siemens AG 2014 All rights reserved.
Page 6 07-20-2014

Right dose technology
Proven innovation and technological leadership

SIEMENS



Restricted © Siemens AG 2014 All rights reserved.
Page 7 07-20-2014


*This product is not available for sale in the US.

CARE Dose4D
Minimal dose, maximum quality and fully automated

SIEMENS

Fully automated dose modulation in real time



- Highest dose savings without IQ compromise
- Protocols automatically adapted
- Modulation based on the topogram and online
- Fully adjustable with CARE Dose Configurator



Restricted © Siemens AG 2014 All rights reserved.
Page 8 07-20-2014

CARE Dose4D
Minimal dose, maximum quality and fully automated

SIEMENS



Different areas of the body
have different shapes

And therefore need
customized dose
(= tube current)


Restricted © Siemens AG 2014 All rights reserved.
Page 9 07-20-2014

For internal use only.

CARE Dose4D

Minimal dose, maximum quality and fully automated


SIEMENS



Restricted © Siemens AG 2014. All rights reserved.
Page 10 07-22-2014

Fully automated dose modulation in x, y and z direction in real time


- Adjustments follow the configurable curve to best fit the clinical needs
- Every organ characteristic has a configurable curve



CARE kV

First automated voltage setting

SIEMENS



Allows contrast-to-noise ratio to be optimized

- Up to 60% dose reduction can be achieved without compromise in image quality
- Zero click automated kV setting for each individual patient
- The three different modes of CARE kV
 - On:** Optimized tube voltage and tube current based on the patient's size and clinical question
 - Off:** The scan is performed with the kV setting defined in the respective scan protocol.
 - Semi:** A user specified kV setting will be used and mAs-values will be adjusted.


Page 11

CARE kV

First automated voltage setting

SIEMENS


Optimization of kV to minimize patient dose at maintained image quality



Page 12

CARE kV
First automated voltage setting


SIEMENS



Page 13

CARE kV
First automated voltage setting

SIEMENS



Page 14

Pediatric Imaging with 70 kV
Aortic angiography

SIEMENS



CARE kV
First automated voltage setting

SIEMENS



SOMATOM
Force

collimation: 192 x 0.6mm
scan time: 1 s
scan length: 740 mm
rotation time: 0.25 s
tube setting:
80 kV, 140 mA
CTDIvol: 2.09 mGy
DIAP: 154.6 mAs
2.35 mSv

80 kV
2.32 mSv
20 ml Contrast
medium


Internal © Siemens AG 2014 All rights reserved.
Page 16 07-22-2014

Courtesy of UMM, Mannheim, Germany

ADMIRE – Advanced Modeled Iterative Reconstruction
The latest generation iterative reconstruction

SIEMENS

Improved IQ at
excellent dose values



On the fly powerful dose reduction, excellent
image quality and everyday suitability

- Potential to lower radiation up to 60% in all clinical applications
- Superb details, delineation and sharpness of organ borders
- Significant artifact reduction, e.g. streak artifacts in the shoulder region
- PACS-Ready workflow even for thick slices
- Reader-Ready reconstructions


Page 17

* In clinical practice, the use of ADMIRE may reduce CT patient dose depending on the selected dose, contrast, and scan protocol. A comparison with a standard and a low-dose scan of the same patient is shown in the representative image to the right. The ADMIRE scan was reconstructed with a standard and a low-dose scan of the same patient using the ADMIRE reconstruction algorithm. Note: CT values, lung volumes, and other parameters were measured in a Siemens 256 phantom. Low-dose data reconstructed with ADMIRE showed the same image quality compared to full-dose data based on this test. Data on file.

ADMIRE – Advanced Modeled Iterative Reconstruction
The latest generation iterative reconstruction

SIEMENS

CTDIvol 0.04 mGy



WFBP

SAFIRE

ADMIRE

Handzettel 6

ADMIRE – Advanced Modeled Iterative Reconstruction
The latest generation iterative reconstruction

SIEMENS

WFBP

ADMIRE

SOMATOM
Force

collimation: 192 x 0.6 mm
scan time: 6.5s
scan length: 326 mm
rotation time: 0.25s
tube setting:
70 kV, 142 mAs
CTDIvol: 1.07 mGy
DLP: 84 mGy cm
a11 down: 1.2 mSv

Page 19

07-22-2014

Courtesy of UMM, Mannheim, Germany

HIMCR

ADMIRE – Advanced Modeled Iterative Reconstruction
The latest generation iterative reconstruction

SIEMENS

The screenshot displays the Siemens ADMIRE software interface. On the left, there are checkboxes for 'Topogram', 'Head', 'Neck', 'Pause', 'Topogram', and 'ThorAbd'. Below these are 'Load', 'Save', and 'Print' buttons. The main area shows a 'Recon job' list with columns for 'Recon job', 'SI-M', 'Strength', 'Kernel', 'FAST', 'Window', and 'Mediastern'. The 'Series description' for 'ThorAbd 7.0 BwD 3' is shown on the right, including 'Recon job type', '3D', 'Recon region', 'Narrow', 'Wide', 'CAP - Lung apex to ischium', 'Begin position', 'End position', 'Image order', 'Cranio-caudal', 'Increment', '7.0 mm', 'No. of images', '29', 'Comments', and 'Auto Tasking'.

Restricted © Siemens AG 2013 All rights reserved.

Page 20

07-22-2014

FAST Planning
Radiation reduction starts with good planning

SIEMENS

The screenshot shows the Siemens FAST Planning software interface. It features a 3D visualization of a patient's anatomy with a blue planning region. Below the visualization, there are various input fields and buttons for planning parameters, including 'Recon job', 'SI-M', 'Strength', 'Kernel', 'FAST', 'Window', and 'Mediastern'. The 'Series description' for 'ThorAbd 7.0 BwD 3' is also visible on the right.

Safe, fast and standardized

- Automated organ-based scan range planning based on anatomical landmarks
- Dose saving by perfect planning
- Auto iso-center positioning
- Optimized workflow

Restricted © Siemens AG 2013 All rights reserved.

Page 21

07-22-2014

CARE Right
Committed to the Right Dose in CT

SIEMENS




Restricted © Siemens AG 2014 All rights reserved.
Page 22 07-22-2014

Right dose levels
Understanding of absolute dose values

SIEMENS


Example: Comparison of radiation exposure due to reduction of flight level



Restricted © Siemens AG 2014 All rights reserved.
Page 23 07-22-2014


Adult Imaging with 80 kV
Aortic CTA

SIEMENS



Right dose levels
Regional reference values as a guidance

SIEMENS

		Switzerland	Germany	Europe	USA
Head		65 mGy	65 mGy	60 mGy	75 mGy
Thorax		15 mGy	12 mGy	30 mGy	n.a.
Abdomen		15 mGy	20 mGy	35 mGy	25 mGy

Average sized patient

Restricted © Siemens AG 2014 All rights reserved.
Page 25 07-22-2014

Right dose levels
Recommendations from experts

SIEMENS

SIERRA


Siemens Radiation Reduction Alliance:

- 15 key opinion leaders in radiology
- Drive research to improve radiation reduction technologies
- Create recommendations for clinical practice everywhere

Restricted © Siemens AG 2014 All rights reserved.
Page 26 07-22-2014

Right dose levels
Base for Siemens CT standard scan protocols

SIEMENS



Restricted © Siemens AG 2014 All rights reserved.
Page 27 07-22-2014

CARE Right
Committed to the Right Dose in CT

SIEMENS



Restricted © Siemens AG 2014 All rights reserved.
Page 28 07-20-2014

Right dose management
The sustainable dose management

SIEMENS

DoseMAP
Access, report and manage dose data

EduCARE
Continuous education of staff

Optimize CARE CT
Optimization of the current dose situation



Restricted © Siemens AG 2014 All rights reserved.
Page 29 07-20-2014

Right dose management
The Siemens CT Dose Management Program - DoseMAP

SIEMENS




Restricted © Siemens AG 2014 All rights reserved.
Page 30 07-20-2014

Analyze Dose Assess the dose situation

SIEMENS

Pre-examination dose checkup

- CARE Dashboard*:
Overview of used dose features per scan
- CARE Profile*:
Visualization of dose distribution prior to scan at every z-position



Restricted © Siemens AG 2014 All rights reserved.
Page 31 07-22-2014


* Availability depending on scanner

Protect Dose Manage access to scan protocols

SIEMENS

Access Protection*

- Increased security by managing user administration rights with Scan Protocol Lock



Restricted © Siemens AG 2014 All rights reserved.
Page 32 07-22-2014


* Availability depending on scanner

Protect Dose Automated monitoring of dose values

SIEMENS

Dose Notification* and Dose Alert*

- May help to protect patients from over-radiation
- Warns the operator in case set dose thresholds are exceeded




Restricted © Siemens AG 2014 All rights reserved.
Page 33 07-22-2014

* Availability depending on scanner

Right dose management EduCARE

SIEMENS



Large portfolio available to increase your skills

- Dedicated Siemens YouTube Channel
- "How to Guides"
- Webinars
- e-learnings
- Siemens Healthcare "Guide to right dose"

Restricted © Siemens AG 2014 All rights reserved.
Page 34 07-20-2014

Right dose management Optimize CARE CT – 4 steps to optimize CT Dose

SIEMENS



Restricted © Siemens AG 2014 All rights reserved.
Page 35 07-20-2014

CARE Right Committed to the Right Dose in CT

SIEMENS



Restricted © Siemens AG 2014 All rights reserved.
Page 36 07-20-2014


CARE Right
Committed to the Right Dose in CT

SIEMENS

CARE Right
Committed to the Right Dose in CT

Restricted © Siemens AG 2014 All rights reserved.
Page 37 07-22-2014

SIEMENS



CARE Right
Computed Tomography.

Restricted © Siemens AG 2014 All rights reserved. Answers for life.
