

**TOSHIBA**  
**Aquillon CT Dose Management Tools**

Erin Angel, PhD, DABR  
Clinical Sciences Manager, CT  
Clinical Collaborations

**TOSHIBA**  
Leading Innovation 360

A slide with a white and red background. The top left features the Toshiba logo and the title 'Aquillon CT Dose Management Tools'. Below the title is the name 'Erin Angel, PhD, DABR' and her title 'Clinical Sciences Manager, CT Clinical Collaborations'. The bottom left corner has the Toshiba logo and the slogan 'Leading Innovation 360'.

---

---

---

---

---

---

---


---

**Agenda**

- > Aquillon Family of CT Scanners
- > Dose Safeguards
- > Dose Management
- > Dose Reduction
- > Dose Measurement
- > Dose Reporting
- > Where to get more info...



**TOSHIBA**  
Leading Innovation 360

A slide with a red header containing the word 'Agenda'. Below the header is a list of seven items, each preceded by a right-pointing chevron. To the right of the list is a photograph of a Toshiba Aquillon CT scanner. The bottom of the slide features the Toshiba logo and the slogan 'Leading Innovation 360'.

---

---

---

---

---

---

---

---

**The Aquillon Family**



**TOSHIBA**  
Leading Innovation 360

A slide with a red header containing the text 'The Aquillon Family'. Below the header are five images of different Aquillon CT scanner models, each with its name and 'VISION EDITION' (where applicable) written below it. The bottom of the slide features the Toshiba logo and the slogan 'Leading Innovation 360'.

---

---

---

---

---

---

---

---

### SURE Exposure: Integrated Dose Reduction

- > Dose Safeguards
- > Dose Management
- > Dose Reduction

The diagram illustrates a four-step workflow: 1. Patient Registration (software interface), 2. Acquisition (X-ray beam hitting a patient), 3. Raw Data (vertical scan image), and 4. Reconstruction (axial CT scan image).

**TOSHIBA**  
Leading Innovation 3D

---

---

---

---

---

---

---

---

---

---

### Integrated Dose Management

Integrated Optimized Workflow to Minimize Dose

Scan	Start	End	Start Pos	End Pos	Scan Width	CTDI vol	DI, P, a
1	P	***	***	0.0	1000.0	0.0	0.0
2	A	***	***	1000.0	0.0	0.0	0.0
3	P	00:00.0	0.0	0.0	0.0	0.0	0.0

**TOSHIBA**  
Leading Innovation 3D

---

---

---

---

---

---

---

---

---

---

### Scan Planning

- > Real-time dose display
- > Dose displayed on console prior to exam

No.	Start	End	Start Pos	End Pos	Scan Width	CTDI vol	DI, P, a
1	P	***	***	0.0	1000.0	0.0	0.0
2	A	***	***	1000.0	0.0	0.0	0.0
3	P	00:00.0	0.0	0.0	0.0	0.0	0.0

**TOSHIBA**  
Leading Innovation 3D

---

---

---

---

---

---

---

---

---

---



### PROTECT™

> PROTECT is a partnership designed to help Aquilion customers make CT safer for all patients.



> PROTECT team will work closely with customer to **evaluate, plan,** and **implement** a strategy to lower dose, maintain image quality and meet the customers PROTECT goals.



TOSHIBA  
Leading Innovation 3D

---

---

---

---

---

---


---

---

---

---

### Dose Reduction Technologies



TOSHIBA  
Leading Innovation 3D

---

---

---

---

---

---


---

---

---

---

### Quantum Detector



TOSHIBA  
Leading Innovation 3D

---

---

---

---

---

---

---

---

---

---

### 320-Row Volume CT Scan Modes

Helical

Ultrahelical

Volume

Johnston et al. "Comparison of radiation dose estimates, image noise, and scan duration in pediatric body imaging for volumetric and helical modes on 320-detector CT and helical mode on 64-detector CT." Pediatric Radiol. 2013 Sep;43(9):1117-27

**TOSHIBA**  
Leading Innovation 3D

---

---

---

---

---

---

---

---

### Dynamic Volume CT: Lower Dose Volume Scanning

- > Why is volume acquisition lower dose than helical?
  - > No helical overranging
  - > Only one rotation of overbeam/penumbra
  - > No low pitch overlap for gated studies

**ONE Aquilion**

**TOSHIBA**  
Leading Innovation 3D

---

---

---

---

---

---

---

---

### Overbeaming

- > Overbeaming is the excess radiation exposure beyond the collimated beam in z-direction
  - > "penumbra"
- > Wider collimation decreases relative penumbra
- > Decrease # of gantry rotations inherently decreases this effect

**TOSHIBA**  
Leading Innovation 3D

---

---

---

---

---


---

---

---

### Overranging

- > Helical scanning requires raw data on either end of planned scan length for reconstruction
- > Unnecessary radiation exposure to tissue beyond desired view
- > Accounts for higher percentage of total exposure with smaller scan lengths (cardiac, neuro, pediatric)
- > Active collimation reduces
- > Does not occur with true axial volume scanning
  - > It's just a really thick axial slice



Schilham A et al. Radiographics 2010;30:1057-1067

TOSHIBA  
Leading Innovation 3D

---

---

---

---

---

---

---

---

---


---

---

---

### Active Collimation

- > Helical overranging minimized for lower dose image acquisition



TOSHIBA  
Leading Innovation 3D

---

---

---

---

---

---

---

---

---

---

---

---

### SureExposure Tube Current Modulation

- > XYZ mA modulation
- > Dual Scanogram based
- > Automatic dose savings

A screenshot of the SureExposure software interface. On the left is a control panel with buttons for 'Scan', 'Back', and 'Done'. In the center is a scanogram of a chest with three horizontal lines (yellow, green, blue) indicating different scan levels. On the right is a graph showing 'mA' modulation across the scan length, with a yellow dot at the top and a blue dot at the bottom, indicating automatic dose modulation.

TOSHIBA  
Leading Innovation 3D

---

---

---

---

---

---

---

---


---

---

---

---

### Tech Assist Lateral Slide



- > 8.4 cm's of Lateral Table Movement
- > Assist Technologist in Positioning/Moving Patient
- > Center Patient in CT Bore, Improved mA Modulation (Dose) and Image Quality

**TOSHIBA**  
Leading Innovation 360

---

---

---

---

---


---

---

---

### Auto Couch Height Compensation

- > SUREExposure uses patient size for exposure calculation
- > Incorrect positioning causes errors in patient size calculation



**TOSHIBA**  
Leading Innovation 360

---

---

---

---

---


---

---

---

### Auto Couch Height Compensation

- > SUREExposure will compensate for incorrect patient positioning to ensure accurate body size calculation and exposure dose



**TOSHIBA**  
Leading Innovation 360

---

---

---

---

---


---

---

---

### Variable Helical Pitch vHp

- > Changing table speeds within the same exam saves time, dose and cost
- > Combining ECG and non ECG gated scans into one makes excellent use of contrast media with significant dose reduction



The diagram illustrates the Variable Helical Pitch (vHp) technique. It shows a cross-section of a patient's torso with the heart and major blood vessels highlighted in purple. A vertical line represents the X-ray beam's path. The top section is labeled 'ECG Gated Helical Pitch', the middle section is 'Table Speed Change' (indicated by a red triangle), and the bottom section is 'Standard Helical Pitch'.

**TOSHIBA**  
Leading Innovation 3D

---

---

---

---


---

---

---

---

### Dose Reduction Technologies



The flowchart shows the process of dose reduction technologies. It starts with 'Patient Registration' (a software interface), followed by 'Acquisition' (a diagram of the X-ray beam), then 'Raw Data' (a grayscale image of a body part), and finally 'Reconstruction' (a color CT scan image).

**TOSHIBA**  
Leading Innovation 3D

---

---

---

---

---

---

---

---

### AIDR 3D Dose reduction



- > **Integrated**  
Preset in CT protocol
- > **Automated**  
Increase clinical workflow
- > **Adaptive**  
Personalized dose reduction

**Dose Reduction Simplified**

**TOSHIBA**  
Leading Innovation 3D

---

---

---

---

---


---

---

---

### Adaptive Iterative Dose Reduction (AIDR) 3D

AIDR 3D is the latest evolution and Toshiba's 3<sup>rd</sup> generation of iterative reconstruction technology that has been fully integrated in to the imaging chain.



The scanner model analyzes the physical properties of the CT system at the time of the acquisition.

**TOSHIBA**  
Leading Innovation 3D

---

---

---

---

---

---

---


---

---

---

### Adaptive Iterative Dose Reduction (AIDR) 3D

AIDR 3D is the latest evolution and Toshiba's 3<sup>rd</sup> generation of iterative reconstruction technology that has been fully integrated in to the imaging chain.



The statistical modeling characterizes both electronic and quantum noise patterns in projection space.

**TOSHIBA**  
Leading Innovation 3D

---

---

---

---

---

---

---


---

---

---

### Adaptive Iterative Dose Reduction (AIDR) 3D

AIDR 3D is the latest evolution and Toshiba's 3<sup>rd</sup> generation of iterative reconstruction technology that has been fully integrated in to the imaging chain.



Adaptive processing in the raw domain based on the measured scanner and statistical modeling that reduces noise and artifacts.

**TOSHIBA**  
Leading Innovation 3D

---

---

---

---

---

---

---

---

---

---

### Adaptive Iterative Dose Reduction (AIDR) 3D

AIDR 3D is the latest evolution and Toshiba's 3<sup>rd</sup> generation of iterative reconstruction technology that has been fully integrated in to the imaging chain.

**Adaptive**                      **Iterative**

A sophisticated iterative technique that optimizes reconstructions for the particular body region being scanned, dramatically reducing noise while preserving spatial resolution and image texture.

**TOSHIBA**  
Leading Innovation 3D

---

---

---

---

---

---

---

---

---

---

---

---

### Adaptive Iterative Dose Reduction (AIDR) 3D

AIDR 3D is the latest evolution and Toshiba's 3<sup>rd</sup> generation of iterative reconstruction technology that has been fully integrated in to the imaging chain.

**Adaptive**                      **Iterative**

**Dose Reduction - AIDR 3D**

**TOSHIBA**  
Leading Innovation 3D

---

---

---

---

---

---

---

---

---

---

---

---

### AIDR 3D Key Features - Integrated

- > **Integrated**
  - > When AIDR 3D is used in a protocol, <sup>SURE</sup>Exposure automatically lowers the tube current prior to imaging the patient adjusting for the expected noise reduction from AIDR 3D. Thus dose reduction is achieved without the need to manually adjust those protocols to achieve lower doses than FBP.

**TOSHIBA**  
Leading Innovation 3D

---

---

---

---

---

---

---

---

---

---

---

---

### AIDR 3D Key Features - Integrated



In clinical practice, the use of AIDR may reduce CT patient dose depending on the clinical task, patient size, anatomical location and clinical practice. A consultation with a radiologist and a physicist should be made to determine the appropriate dose to obtain diagnostic image quality for the particular clinical task.

**TOSHIBA**  
Leading Innovation 3D

---

---

---

---

---

---

---

---

---

---

---

---

### AIDR 3D Key Features - Automated

- > Automated
  - > AIDR 3D is pre-set in the protocol
  - > No need for an extra post processing step
  - > No extra steps
  - > Implementing AIDR 3D on-site is virtually invisible to the technologists with exception of lower doses and/or improved IQ. No additional training necessary.

**TOSHIBA**  
Leading Innovation 3D

---

---

---

---

---

---

---

---

---

---

---

---

### AIDR 3D Key Features - Adaptive

- > Adaptive
  - > No need to select levels
  - > AIDR 3D adapts to the body region/imaging category
  - > Several parameters within AIDR 3D are automatically adapted such as:
    - > Blending level
    - > Number of iterations,
    - > Considers edge detail and noise sensitivity
    - > Strength of signal at detectors
- > A note on adaptive
  - > Although AIDR 3D is fully adaptive, it is equipped with three strengths ("mild", "standard", and "strong") primarily for research
  - > The "standard" setting is used for majority of protocols. The "strong" setting is used for the less texture-sensitive task of body perfusion.

**TOSHIBA**  
Leading Innovation 3D

---

---

---

---

---

---

---

---

---

---

---

---









Thank You

**TOSHIBA**  
Leading Innovation 3D

---

---

---

---

---

---

---

---

**TOSHIBA**  
Leading Innovation >>>

**TOSHIBA**  
Leading Innovation 3D

---

---

---

---

---



---

---

---

**Aquilion CT, Double Slice Technology**

- > Reconstruct double the slices per rotation
- > Up to 0.5 mm x 320 detector row
- > Up to 640 slice reconstruction
- > 3D image quality improvement
- > Partial volume effects minimized
- > Data acquired as a volume

Without Double Slice Technology	With Double Slice Technology
	

No Additional Dose

**TOSHIBA**  
Leading Innovation 3D

---

---

---

---

---

---

---

---

### Aquilion CT, Double Slice Technology

**Double Slice Technology**  
coneXact™

80 x 0.5mm Detector Rows

**No Additional Dose**

**TOSHIBA**  
Leading Innovation 3D

---

---

---

---

---

---

---

---

### Noise Reduction

> Volume Mode, 40 cm water phantom, 120 kVp, LFOV, 80x0.5 mm

**Abd Protocol (0.5mm image thickness)**

**Lung Protocol (0.5mm image thickness)**

**TOSHIBA**  
Leading Innovation 3D

---

---

---

---

---

---

---

---

### Spatial Resolution

> Abdomen Protocol, 120 kVp

**Modulation Transfer Function (MTF)**

Method	Resolution (lp/cm)
FBP	~10
	~5
	~2
ADR	~15
	~10
	~5

**TOSHIBA**  
Leading Innovation 3D

---

---

---

---

---


---

---

---

## Texture Preservation

➤ Abdomen protocol, 250 mAs for FBP and 75 mAs for ADR 3D (70% dose difference)



**TOSHIBA**  
Leading Innovation 3D

---

---

---

---

---

---

---

---