

# Using DoseLab to Perform TG-142 Imaging QA

Nathan Childress, Ph.D., DABR



# Introduction

## **Task Group 142 report: Quality assurance of medical accelerators<sup>a)</sup>**

Eric E. Klein<sup>b)</sup>

*Washington University, St. Louis, Missouri*

Joseph Hanley

*Hackensack University Medical Center, Hackensack, New Jersey*

John Bayouth

*University of Iowa, Iowa City, Iowa*

TG-142 is a comprehensive QA protocol

- Covers nearly every aspect of machine and safety QA
- Recommends quantitative results
- Recommends high testing frequencies (ie, monthly imaging QA)
- Essentially requires specialized QA software to perform

# My Background

- ▶ 2001–2004
  - Developed open-source DoseLab for PhD at MD Anderson
- ▶ 2004–2010
  - Worked at Methodist Hospital as a clinical physicist
- ▶ 2008–2010
  - Commissioned a Varian 21iX at a satellite facility
  - Implemented a TG-142 program using PipsPro
- ▶ 2010–present
  - Developed DoseLab TG-142 to perform automatic QA
  - Founded Mobius Medical Systems, LP to manufacture and support DoseLab, FractionLab, Mobius3D, and MobiusFX




# My TG-142 Implementation in 2009

- ▶ I noticed many inefficiencies
- ▶ The software did not extract the maximum amount of results from each image
- ▶ Trips inside the vault to setup phantoms were not minimized
- ▶ Not all modules stored results in a database, so custom Excel sheets had to be created



# Why Develop DoseLab TG-142?


- ▶ Maximize QA efficiency
  - Fewer images & trips to vault
  - Standardized PDF reporting
  - Automatic computations
- ▶ Quantify and classify all results
  - **Pass**/**warn**/**fail** tolerances
  - Database trending
- ▶ Comprehensive tool set
  - Compatible with EPID, film, CR, etc.
  - Compatible with phantoms from SNC, Standard Imaging, etc.
  - Log file analysis for Varian / Elekta



DoseLab is used by the RPC for  
on-site TG-142 audits

# Recommended TG-142 Strategy

- ▶ Create a QA patient in your R+V system
- ▶ Add fields for all needed EPID, kV, and CBCT measurements
- ▶ Use all electronic measurements – no film
- ▶ Use your R+V system to control the linac and store images
- ▶ Export images to TG-142 software and analyze
- ▶ This makes future QA easier
  - Open QA patient
  - Deliver fields
  - Export and analyze






# Monthly MV and kV Imaging

| Daily <sup>a</sup>  |             |             |
|---|-------------|-------------|
| <b>Planar kV and MV (EPID) imaging</b>                                |             |             |
| Collision interlocks  | Functional  | Functional  |
| Positioning/repositioning   | $\leq 2$ mm | $\leq 1$ mm |
| Imaging and treatment coordinate coincidence<br>(single gantry angle) | $\leq 2$ mm | $\leq 1$ mm |

| Monthly  |                       |             |
|--|-----------------------|-------------|
| <b>Planar MV imaging (EPID)</b>  |                       |             |
| Imaging and treatment coordinate coincidence<br>(four cardinal angles) | $\leq 2$ mm           | $\leq 1$ mm |
| Scaling <sup>b</sup>   | $\leq 2$ mm           | $\leq 2$ mm |
| Spatial resolution   | Baseline <sup>c</sup> | Baseline    |
| Contrast   | Baseline              | Baseline    |
| Uniformity and noise   | Baseline              | Baseline    |
| <b>Planar kV imaging<sup>d</sup></b>                                   |                       |             |
| Imaging and treatment coordinate coincidence<br>(four cardinal angles) | $\leq 2$ mm           | $\leq 1$ mm |
| Scaling  | $\leq 2$ mm           | $\leq 1$ mm |
| Spatial resolution   | Baseline              | Baseline    |
| Contrast   | Baseline              | Baseline    |
| Uniformity and noise   | Baseline              | Baseline    |

# Supported Phantoms

- ▶ DoseLab has a unique approach to multiple phantom support
- ▶ DoseLab supports all common phantoms:
  - SNC MV and kV ImagePro phantoms
  - Leeds TOR 18FG (included with IGRT linacs)
  - Las Vegas (no spatial resolution, included with linacs)
  - Standard Imaging QC-3
  - Standard Imaging QC-kV1
  - iba DIGI-13
  - PTW EPID QC (no positioning or scaling)
  - User-customizable additions






# Phantom Selection

- ▶ Your linac comes with Leeds and Las Vegas phantoms
- ▶ Leeds has all required modules to fulfill TG-142, but is difficult to set up for automatic analysis
  - No phantom stand
  - No crosshair marks
- ▶ Las Vegas does not have spatial resolution segments
  - And still no phantom stand or crosshair marks
- ▶ Most users purchase aftermarket phantoms to overcome these difficulties

# Aftermarket Phantoms

- ▶ Easy setup – they include a stand and crosshair markings
- ▶ Have all required modules to fulfill TG-142
- ▶ Work well with automatic analysis




# MV and kV Phantom Setup




# Image Acquisition

- ▶ Step 1 – Perform manufacturer's recommended acceptance test to verify correct system performance
- ▶ Step 2 – Perform monthly tests using clinical imaging protocols to set baseline values






# DoseLab MV and kV QA Analysis



- ▶ 1 image/imaging energy (2 total)
- ▶ MV/kV phantoms
- ▶ Analyzes all TG142 parameters
- ▶ Includes scaling and positioning, without separate image

# Imaging Calculations




- ▶ Spatial resolution – uses MTF rather than subjective “How many line pairs can I see?”
- ▶ Contrast
- ▶ CNR
- ▶ Uniformity
- ▶ Scaling
- ▶ Positioning offset
- ▶ CT only:
  - Geometric distortion
  - Slice width
  - HU deviation





# Imaging Parameters: Offsets


- ▶ IGRT systems are calibrated to know the isocenter position in each image
- ▶ To calculate the offset from expected (radiation or laser isocenter), a target can be placed at isocenter and imaged
- ▶ For planar imaging, DoseLab extracts the isocenter location in DICOM tags and compares the phantom position to a baseline position
  - This occurs during image quality analysis



# CBCT Catphan Setup



# CT and CBCT QA



- ▶ 1 image set
- ▶ Catphan phantom
- ▶ Analyzes all TG142 parameters
- ▶ Also supports Gammex, CIRS, and GE phantoms

# Image Quality Tolerances

| Monthly  |                       |             |
|--|-----------------------|-------------|
| <b>Planar MV imaging (EPID)</b>  |                       |             |
| Imaging and treatment coordinate coincidence<br>(four cardinal angles) | $\leq 2$ mm           | $\leq 1$ mm |
| Scaling <sup>b</sup>   | $\leq 2$ mm           | $\leq 2$ mm |
| Spatial resolution   | Baseline <sup>c</sup> | Baseline    |
| Contrast   | Baseline              | Baseline    |
| Uniformity and noise   | Baseline              | Baseline    |
| <b>Planar kV imaging<sup>d</sup></b>                                   |                       |             |
| Imaging and treatment coordinate coincidence<br>(four cardinal angles) | $\leq 2$ mm           | $\leq 1$ mm |
| Scaling  | $\leq 2$ mm           | $\leq 1$ mm |
| Spatial resolution   | Baseline              | Baseline    |
| Contrast   | Baseline              | Baseline    |
| Uniformity and noise   | Baseline              | Baseline    |

- ▶ Results depend on many factors
  - Software (several sets of formulas exist)
  - Phantom
  - Imaging technique (kV, mAs)
  - Setup (at isocenter, on panel, through couch)
- ▶ Manufacturers / AAPM do not recommend tolerances
- ▶ Your clinic's data is used to establish your baseline and tolerance levels

# Image Quality Tolerances


| Monthly  |                       |             |
|--|-----------------------|-------------|
| <b>Planar MV imaging (EPID)</b>  |                       |             |
| Imaging and treatment coordinate coincidence<br>(four cardinal angles) | $\leq 2$ mm           | $\leq 1$ mm |
| Scaling <sup>b</sup>   | $\leq 2$ mm           | $\leq 2$ mm |
| Spatial resolution   | Baseline <sup>c</sup> | Baseline    |
| Contrast   | Baseline              | Baseline    |
| Uniformity and noise   | Baseline              | Baseline    |
| <b>Planar kV imaging<sup>d</sup></b>                                   |                       |             |
| Imaging and treatment coordinate coincidence<br>(four cardinal angles) | $\leq 2$ mm           | $\leq 1$ mm |
| Scaling  | $\leq 2$ mm           | $\leq 1$ mm |
| Spatial resolution   | Baseline              | Baseline    |
| Contrast   | Baseline              | Baseline    |
| Uniformity and noise   | Baseline              | Baseline    |

- ▶ Routine QA is used in conjunction with acceptance testing to determine that a system is operating properly and stays operating properly
- ▶ Acceptance tests are typically very different than quantitative routine QA tests



# Imaging QA and Non-Standard Formulas

- ▶ There is no industry standard set of basic formulas in diagnostic imaging or radiation therapy
- ▶ DoseLab, like nearly every other software package, uses its own set of imaging QA formulas
- ▶ DoseLab's formulas were designed to produce consistent results for automatically-placed ROIs






# Comparing Imaging QA Results

- ▶ In therapy, we are used to being able to directly compare performance data between machines and hospitals
- ▶ This is absolutely not the case with imaging QA, due to results depending on formulas, phantoms, setup, imaging techniques, etc.
- ▶ Even holding these parameters constant between different machines can lead to quantitative differences that do not indicate performance issues
- ▶ TG-142 recommends comparing results to “Baseline”

# CBCT Positioning




- ▶ 2 EPID images after CBCT positioning
- ▶ WL-QA phantom
- ▶ Many ways to accomplish this
- ▶ WL is easy (<5 min) and accurate within 0.1 mm
- ▶ Can be repeated daily for frameless SRS / SBRT

# Monthly QA Summary

- ▶ 2 photon energies with IGRT, CBCT, and VMAT
  - 15 images analyzed (1 is a CBCT set)
  - 5 phantoms
  - 7 software modules (Or 1 module – AutoQA)
  - 13 PDF reports (can be merged into one)
  - 170 numerical values written to database
  - Everything can be performed with EPID / kV imagers
- ▶ ~70 minutes, after clinic's initial setup period
  - Machine time: 45 minutes
  - Export images from R+V system: 10 minutes
  - DoseLab / FractionCHECK analysis: 15 minutes

# Result Documentation



- ▶ Documentation is essential
- ▶ PDF reports can show original image and ROIs
- ▶ Database can save and trend numeric results

# Issues Discovered During TG-142

- ▶ Scaling discrepancy  $> 3$  mm: Imager SID needs recalibration
  - May also affect CBCT positioning
- ▶ Contrast or uniformity issues: Imaging panel needs recalibration (dark field / flood field recalibration)
- ▶ Several bad pixels: Recalibrate or replace imaging panel
- ▶ Pin cushion distortion in image: Replace imaging panel
- ▶ CBCT imaging issues typically indicate need to recalibrate

# Web Resources

- ▶ “A Practical Guide to TG-142 QA”  
<http://www.medphysfiles.com/>
  - 50 page TG-142 procedure guide
  - XLS file of TG-142 tables
  - Written by Jimmy Jones
  
- ▶ “Modified Winston-Lutz Test for IGRT Setups”  
<http://www.medphysfiles.com/>
  - CBCT Winston-Lutz Guide
  - Written by Nathan Childress





# Summary

- ▶ TG-142 requires a lot of QA
- ▶ Specialized software is necessary to be TG-142 compliant
- ▶ Imaging QA is very different from therapy QA
  - Results difficult to compare
  - Nearly everything based on deviation from baseline
- ▶ Software can automate tasks:
  - Analysis of all results
  - Documentation
  - Database saving and trending