TOSHIBA Leading Innovation	57 ^m AAPM Annual Meeting, July 12-16, Anaheim, CA
	vanced Statistical Iterative
Keconst	ruction Method on Toshiba
	CT System
٠	Wenli Wang
Made for Life CONFIDENTIAL &	PROPRETARY TOSHIBA MEDICAL RESEARCH INSTITUTE USA

Outline

- History of ${f S}$ tatistical ${f I}$ terative ${f R}$ econstruction (SIR)
- · Key components of SIR
- · Potential benefits and challenges of SIR
- · Toshiba's SIR: FIRST
- Challenges & methods in evaluating SIR image quality
- Phantom and clinical image comparison of FIRST vs. FBP
- Quiz

 $\begin{tabular}{ll} \textbf{Scientific field: "Statistical Iterative Reconstruction"} &= \textbf{Marketing world} \\ "Model-based Iterative Reconstruction" \\ \end{tabular}$

Disclaimer: FIRST is currently pending FDA 510(k) clearance. Availability of this software is regionally dependent. Please check with your local Toshiba representative to determine availability.

TOSHIBA Leading Innovation

CONFIDENTIAL & PROPRIETARY

OSHIBA MEDICAL RESEARCH INSTITUTE USA

History of statistical iterative reconstruction

- Statistical recon has been routinely used in PET & SPECT
 - Shepp LA, Vardi Y. Maximum likelihood reconstruction for emission tomography. IEEE Trans Med Imag. 1982
- First statistical recon for mono-energy transmission tomography
 - Lange K, Carson R. EM reconstruction algorithms for emission and transmission tomography. J Comput Assist Tomogr 1984
- First statistical recon for x-ray CT with Bayesian estimation
 - Sauer K, Bouman C. A local update strategy for iterative reconstruction from projections. IEEE Trans Sig Proc 1993
- Why SIR from emission tomography is not directly transferable to CT?
 - X-ray source has wide energy spectrum instead of mono-energy
 - $\,-\,$ CT is energy integration detector instead of photon counting
 - CT has high photon flux/data sampling, FBP gives "satisfactory" result
 - CT requires complex data calibration and pre-processing
 - CT requires almost real time reconstruction speed

т	0	S	н	ı	В	Α
	- 40			_		


CONFIDENTIAL & PROPRIETA

SHIBA MEDICAL RESEARCH INSTITUTE USA


Key components of statistical iterative recon		
• Objective function $f^* = \underset{f}{\operatorname{argmax}} \underbrace{\left(\underset{f}{\Pr(f g)} \right)} = \underset{f}{\operatorname{argmax}} \underbrace{\left(\underset{f}{\log \Pr(g f)} + \underset{log-filedihood}{\log \Pr(f)} \right)}_{\text{Log-filedihood}}.$	-	
Posterior probability Log-likelihood Log-prior Log-likelihood/data consistency term Data domain: pre-log vs post-log for CT Data statistical modeling: Gaussian, Poisson, Compound Poisson for CT One example: post-log data domain, Gaussian distribution (Weighted-Least-Square)		
Log-prior/regularizer term (optional) Regularize the ill-posed data consistency term Constrain neighborhood pixels Vast forms: convex versus non-convex MRF, variational function, sparsity prior,	-	
nonlocal-mean patch, etc One example: smooth TV, enhance smoothness => small pixel value difference preserveedge => Large pixel value difference		
TOSHIBA CONFEDENTIAL & PROPRETARY TOSHIBA MEDICAL RESEARCH INSTITUTE USA 4	Ι.	
Potential benefits of statistical iterative recon Be able to model geometry & physics effects into system matrix	١ .	
Rexible to accommodate any scanner geometry Can model blur due to focal spot size, detector element & shift-variant anode angle Can model physical factors, such as beam-hardening & scatter Be able to design special statistical weight to suppress noise or artifacts Be able to design special regularizer to achieve uniform local impulse	-	
response function, at least in uniform regions Be able to adjust the global strength of the regularizer to achieve desired global image noise appearance	-	
TOSHIBA CONFIDENTIAL LPROPRETARY TOSHIBA MEDICAL RESEARCH INSTITUTE USA 5	١.	
Potential challenges of statistical iterative recon Objective function (especially w/ regularizer) is quite complex and thus		
needs optimization algorithm to get to the global maximum solution. — PCG, ICD (Sauer & Bouman), SPS for non-quadratic prior (U. Michigan), ASD-POCS (U. Chicago), Chambolle-Pock (U. Chicago), ADMM (U. Michigan), FISTA (Beck and Teboulle), etc.		
Due to the complexity of CT system matrix modeling, algorithm takes longer iterations to converge. To get to the right solution quickly, need		
 Algorithm acceleration GPU/CPU acceleration 	-	
	-	

Toshiba's statistical iterative recon:

FIRST) rward-projected Model-based Iterative Reconstruction SoluTion



Factors affecting statistical recon image quality



TOSHIBA

Lading Innovation >>>

CONFIDENTIAL & PROPRETARY

TOSHIBA MEDICAL RESEARCH INSTITUTE USA 8

Challenges & methods in evaluating SIR (1)

- Image noise is controlled by
 - Dose
 - Statistical weight design
 - Regularizer design and its internal parameter(s)
 - Regularizer global strength
- · Image noise may also vary with pixel's
 - Spatial location
 - Contrast level
- Image noise may even affect noise texture
 - NPS
 - Modified Q-matrix (ref: Zhang M. et. Al., Blocky artifacts assessment in PET images, MIC 2015)

TOSHIBA Leading Innovation >>>	CONFIDENTIAL & PROPRIETARY	TOSHIBA MEDICAL RESEARCH INSTITUTE USA	Ď,
Leading Innovation >>>			

Challenges & methods in evaluating SIR (2)

- · Image noise definitions
 - $-\sigma^2$ in a VOI (e.g., N pixels) in one noisy sample

$$\sigma^{2} = \frac{1}{N} \sum_{(x,y,z)} (I(x,y,z) - \overline{I})^{2}$$

 $- \ \sigma^2\!(x,\,y,\,z) \ in \ noisy \ ensembles \ (e.g.,\,M \ samples)$

$$\sigma^{2}(x, y, z) = \frac{1}{M} \sum_{m=1}^{M} (I_{m}(x, y, z) - \overline{I}(x, y, z))^{2}$$

- Noise power spectrum (NPS)

$$NPS(f_x, f_y, f_z) = \frac{\Delta x \Delta y \Delta z}{N_x N_y N_z} \frac{1}{M} \sum_{m=1}^{M} \left(DFT_{3D} \left[I_m(x, y, z) - \overline{I}(x, y, z) \right] \right)^2$$

- Reference:
 - Li K, Tang J, Chen GH, Statistical model based iterative reconstruction (MBIR) in clinical CT systems: experimental assessment of noise performance, Medical Physics, 2014

TOSHIBA

CONFIDENTIAL & PROPRIETAR

SHIRA MEDICAL RESEARCH INSTITUTE USA 1


Challenges & methods in evaluating SIR (3)

- Image resolution is controlled by
 - Data sampling
 - Forward/back projector
 - Regularizer design and its internal parameter(s)
- Image resolution may depend on contrast level, pixel location & photon counts
 - Contrast level: use Catphan or ACR phantom's disk inserts to evaluate ESF-> MTF in different materials
 - Pixel location: use large size phantom, or shift phantom so that disk insert has different distance to scanner's ISO center
 - Photon counts: different mAs
- Reference:
 - Li K, Garrett J, Ge YS, Chen GH, Statistical model based iterative reconstruction (MBIR) in clinical CT systems, Part II. Experimental assessment of spatial resolution performance, Medical Physics, 2014

TOSHIBA Leading Innovation


CONFIDENTIAL & PROPRIETA

TOSHIBA MEDICAL RESEARCH INSTITUTE USA 11



FIRST has lower NPS than FBP

- FIRST's NPS is lower than FBP at all dose levels
- At the same dose, when regularizer strength increases FIRST's NPS gets lower and its peak frequency shifts towards the low-frequency.




TOSHIBA

CONFIDENTIAL & PROPRIETAR

OSHIBA MEDICAL RESEARCH INSTITUTE USA 1

FIRST NPS's peak-freq vs dose


- FBP's NPS curve shape (peak-freq) remains the same for all dose levels
- FIRST's curve shape (peak-freq) has negligible shift with reducing dose if its regularizer strength is the same.



TOSHIBA Leading Innovation


CONFIDENTIAL & PROPRIETARY

TOSHIBA MEDICAL RESEARCH INSTITUTE USA 14



Advantage of FIRST: improved MTF

- FBP's MTF does not depend on dose or contrast
- FIRST's MTF depends more on contrast less on dose




TOSHIBA Leading Innovation

CONCIDENTIAL & DECEDIETAD

TOSHIRA MEDICAL PESEARCH INSTITUTE USA 16

Advantage of FIRST: LCD w/ model observer





SKE/BKE NPWE observer 2AFC experiment

TOSHIBA Leading Innovation

ONFIDENTIAL & PROPRIETA

OSHIBA MEDICAL RESEARCH INSTITUTE USA 1

FIRST: lung with dose reduction EBP. 188mAs WL/WW-40/400 FIRST, 94 mAs WL/WW-40/400 WL/WW-550/1600 0.5 mm dice 'Courtesy of NHLBI'

FIRST: abdomen with better resolution & noise FBP, 166 mAs FIRST, 166 mAs Stockles-35 m stockles-

Summary and conclusion

- FIRST is a fully statistical reconstruction algorithm, with system matrix modeling cone-beam geometry and x-ray optics, edge-preserving regularizer reducing noise but preserving fine anatomical structure.
- $\bullet \;\;$ FIRST can reduce dose but preserve image quality as FBP at normal dose.
- At the same dose, FIRST has higher image resolution/contrast and lower noise than FBP.

TOSHIBA
Leading Innovation

CONFIDENTIAL & PROPRIETARY

TOSHIBA MEDICAL RESEARCH INSTITUTE USA

Special thanks to

Kirsten Boedeker Jian Zhou Sachin Moghe Kurt Schultz Naruomi Akino Satoru Nakanishi

from Toshiba for providing data and plots.

TOSHIBA

CONFIDENTIAL & PROPRIETA

TOSHIBA MEDICAL RESEARCH INSTITUTE USA 21

Quiz (1)	()
The key difference between statistical reco reconstruction is A. Statistical recon is iterative	·
B. Statistical recon assumes data follows a concentration of the statistical recon reduces dose D. All of above	ertain statistical model
TOSHIBA COMFEDENTIAL & PROPRETARY TOSHI Leading innovation >>>	BA MEDICAL RESEARCH INSTITUTE USA 22
Quiz (1)	
The key difference between statistical reco reconstruction is A. Statistical recon is iterative	
B. Statistical recon assumes data follows a concentration of the statistical recon reduces dose D. All of above	ertain statistical model
Correct Answer: B	
TOSHIBA CONFOCNTIAL & PROPRETARY TOSHI	BA MEDICAL RESEARCH INSTITUTE USA 23
Quiz (2)	
 In statistical reconstruction, which factor d resolution: A. Data sampling 	oes <u>not</u> impact image
B. Statistical weight C. Forward & back projector D. Regularizer	
TOCHIDA	
TOSHIBA CONFIDENTIAL & PROPRIETARY TOSHI Leading Innovation >>>	BA MEDICAL RESEARCH INSTITUTE USA 24

Quiz (2)	
 In statistical reconstruction, resolution: A. Data sampling B. Statistical weight 	which factor does not impact image
B. Statistical Weight C. Forward & back projector D. Regularizer	
Correct Answer: B	
TOSHIBA Leading Innovation >>> CONFIDENTIAL & PR	ROPRETARY TOSHIBA MEDICAL RESEARCH INSTITUTE USA 25
Quiz (3)	100
In statistical reconstruction, A. Statistical weight B. Forward & back projector C. Regularizer	which factor(s) impact image noise:
D. A and B E. A and C F. B and C	
TOSHIBA Leading Innovation >>> CONFIDENTIAL & PR	ROPRETARY TOSHIBA MEDICAL RESEARCH INSTITUTE USA 26
Quiz (3)	
A. Statistical weightB. Forward & back projectorC. Regularizer	which factor(s) impact image noise: r
D. A and B E. A and C F. B and C	
Correct Answer: E	
TOSHIBA Leading Innovation >>>	ROPRIETARY TOSHIBA MEDICAL RESEARCH INSTITUTE USA 27