

# MRI signal in biological tissues Proton, Spin, T1, T2, T2\*

Krzysztof Gorny, Ph.D., DABR

Associate Professor of Medical Physics Mayo Clinic College of Medicine, Rochester, MN

#### **Outline**

- Nuclear magnetism
   nuclear magnetic moment
   protons spin ½ system
   Larmor Equation
   Resonance
- Resonance
   Magnetic Resonance in tissues
   Longitudinal magnetization; spin-lattice relaxation time, T1
   RF pulses
   Transverse magnetization; T2\*, T2 relaxation
   Spin Echo
   T1, T2 contrast

#### MRI scanner

Very strong (1.5T, 3.0T, 7T), uniform magnetic field – to magnetically polarize the patient
 Radio-frequency coils – to generate and detect MRI signal in tissues (at resonance)
 Gradient coils – to spatially encode signal for reconstruction (loud noise)







What is the source of MRI signal? Signal in conventional MRI comes from Hydrogen nuclei (a.k.a. "protons") • Hydrogen atom • Negative electron nucleis – positive proton Forcom



### What is the source of MRI signal? Signal in conventional MRI comes from Hydrogen nuclei (a.k.a. "protons")







## Hydrogen nuclei ("protons")







#### Spin-lattice relaxation



#### Spin-lattice relaxation





#### Spin-lattice relaxation



#### Spin-lattice relaxation



#### Spin-lattice relaxation

- When B<sub>0</sub> = 0 proton moments are oriented randomly. Tissue magnetization = 0.
   When B<sub>0</sub> > 0 proton moments rearrange exchanging energy with the lattice.
   Preferential orientation is parallel with the field (lower energy)
   Tissue magnetization builds up
   In thermal equilibrium ratio moment populations are given by the Boltzman distribution



#### Spin-lattice relaxation

- When B<sub>0</sub> = 0 proton moments are oriented randomly. Tissue magnetization = 0.
   When B<sub>0</sub> > 0 proton moments rearrange exchanging energy with the lattice.
   Preferential orientation is parallel with the field (lower energy)
   Tissue magnetization builds up
   In thermal equilibrium ratio moment populations are given by the Boltzman distribution
   At 1.5T for every 1 million protons there is excess of 5 "parallel" protons over the "anti-parallel"



#### Spin-lattice relaxation

T1 -- time scale for reaching equilibrium distribution of proton moments: spin-lattice relaxation time







		-





#### Magnetic Resonance

- The RF pulse applied at resonance tips equilibrium magnetization away from  $B_{\boldsymbol{\theta}}$ + "90° pulse" rotates magnetization by  $90^\circ$ 



\_

**Magnetic Resonance** 

- The RF pulse applied at resonance tips equilibrium magnetization away from  $B_{\boldsymbol{\theta}}$
- "90° pulse" rotates magnetization by 90° "180° pulse" rotates magnetization by 180° (populations of proton moments inverted)



#### Magnetic Resonance

- Immediately after the 90° pulse
   There is zero longitudinal magnetization. Entire magnetization is flipped to the transverse plane and becomes transverse magnetization

   The transverse magnetization precesses around B<sub>6</sub> with Larmor frequency

Precessing transverse magnetization is MRI signal













Transverse magnetization, T2\*, T2

- Following 90° pulse the signal induced by rotating transverse magnetization eventually disappears
- This relaxation has exponential envelope with a characteristic time-scale,  ${T_2}^{\star}$



#### Transverse magnetization, T2\*, T2

There are two major contributions to T2\* relaxation (besides T1 process): 1. Inhomogeneilies of the static magnetic field 2. T2 process

#### Transverse magnetization, T2\*, T2

- There are two major contributions to T2\* relaxation (besides T1 process): 1. Inhomogeneities of the static magnetic field 2. T2 process
- Magnetization is the **sum** of individual contributions built from proton magnetic moments



- There are two major contributions to T2\* relaxation (besides T1 process):
- 1. Inhomogeneities of the static magnetic field
   2. T2 process
   Magnetization is the sum of individual contributions built from proton magnetic moments
   Immediately after 90° pulse moments precess coherently or in-phase with one another



#### Transverse magnetization, T2\*, T2

- There are two major contributions to T2\* relaxation (besides T1 process):

  Inhomogeneities of the static magnetic field
  T2 process

  Dephasing due to field inhomogeneities can be eliminated using Spin Echo sequence.



#### Transverse magnetization, T2\*, T2

- There are two major contributions to T2\* relaxation (besides T1 process): 1. Inhomogeneities of the static magnetic field 2. T2 process
- Dephasing due to field inhomogeneities can be eliminated using Spin Echo sequence.



#### Transverse magnetization, T2\*, T2

#### There are two major contributions to T2\* relaxation (besides T1 process): 1. Inhomogeneities of the static magnetic field 2. T2 process

- - process Even if B<sub>y</sub> was perfectly homogeneous, transverse magnetization would still diminish Protons interact with randomly fluctuating magnetic fields from their neighbors Moments randomly de-phase and lose their coherence transverse magnetization decreases and ultimately disappears



TR (repetition time) - time between 90° pulses (excitation

# T1 and T2 contrast T1 and T2 are properties of tissue

#### T1 and T2 contrast



















