



Intra-Operative Radiation Therapy using Mobile Linear Accelerators

Sam Beddar, Ph.D., FCCPM, FAAPM
Professor & Chief of Research


Department of Radiation Physics
The University of Texas MD Anderson Cancer Center

THE UNIVERSITY OF TEXAS

**MD Anderson
Cancer Center**

Making Cancer History®

1905 – Carl Beck



History of IORT

- 1st report of IORT in 1907 by Carl Beck
- 1964 Abe in Japan
- 1st US IORT in 1975 @ Howard University
- MGH 1978, NCI 1979, Mayo Clinic 1981, Joint Center 1982, MCO 1983



- 1st Mobile LINAC: Novac7 in Italy in 1997
- 1st Mobile LINAC Mobetron prototype in the US: UCSF in 1998
- 1st Commercial Mobetron installed @ University Hospitals of Cleveland in 1999

Before the mobile systems


- The tremendous logistics to transport the patient under anesthesia from the OR to radiation oncology deterred many centers from implementing the procedure
- Many centers lost interest after only few cases



Novac 7



Mobetron 2000




Mobetron 2000



Mobetron 1000

Transport Mode



Treatment Mode



Commissioning – before & after



... It goes around corners



... to the OR floor



...it fits in elevators,
but
can't walk up stairs yet


... in the OR



Control console outside the OR



Applicator Placement




Almost ready to treat




Soft Docking – Laser Alignment System




Treatment Delivery





Multidisciplinary: Roles & Team Work







Advantages of IORT



Advantages of IORT


- **Provides a large dose of radiation to the tumor/ tumor bed at the time of surgery, while normal tissues can be displaced from the radiation field**
- **Potentially decreases side effects and complications of radiation therapy**
- **Shortens overall treatment time for the patients by decreasing the number of visits to the Radiation Oncology Department**
- **Social niche (breast cancer)**
 - **poor countries, with no access to breast conservation - an alternative**
 - **rich countries, where patients can't fit RT treatment into their schedules**

Advantages of IORT



Machine Characteristics

- Electron energies: 4, 6, 9 and 12 MeV with dose rates up to 10Gy/min
- Flat and beveled electron applicators (3.0 to 10.0 cm, 0.5 cm increments)
- Plastic Boluses (0.5 or 1.0 cm)
- Soft-docking system
- 50 cm nominal SSD, non-isocentric
- 6 degrees of freedom (3 transl., 3 rot.)
- Beam stopper tracks the beam (self-shielded)
- No need for shielding in the OR
- Plugs into normal 3-phase electrical outlet / single-phase input



Wootton et. al., 2016

Why No Shielding Requirements?

- Electron Beam Only– Low beam current greatly reduces radiation leakage.
- No Bending Magnet – The most significant source of leakage in a conventional accelerator is eliminated.
- X-ray Contamination – Extremely Low!



A self-shielded accelerator




Fig. 2. uGy/1000 Monitor Units at 2 meters in the Y Plane

Intraoperative radiation therapy using mobile electron linear accelerators: Report of AAPM Radiation Therapy Committee Task Group No. 72

TABLE OF CONTENTS

| | | | |
|----------------------------------------------------------------------------------------|------|---------------------------------------------------------------------------------------|------|
| I. INTRODUCTION..... | 1477 | C. Radiation survey..... | 1479 |
| II. IORT USING MOBILE VERSUS STATIONARY LINEAR ACCELERATORS..... | 1477 | D. Exposure rate measurements..... | 1479 |
| III. IMPLEMENTATION OF AN IORT PROGRAM WITHIN AN OPERATING ROOM ENVIRONMENT..... | 1478 | V. ACCEPTANCE TESTING AND COMMISSIONING..... | 1480 |
| IV. RADIATION PROTECTION..... | 1478 | A. Acceptance testing measurements..... | 1480 |
| A. Regulatory considerations..... | 1479 | B. Commissioning and dose measurements..... | 1482 |
| B. Radiation site plan..... | 1479 | VI. RECOMMENDED QUALITY ASSURANCE.... | 1484 |
| 1. Treatment operating rooms..... | 1479 | A. Previous quality assurance recommendations for medical linear accelerators..... | 1484 |
| 2. Commissioning and annual quality assurance location..... | 1479 | B. Quality assurance for mobile electron accelerators..... | 1484 |
| | | VII. CLINICAL ASPECTS OF IORT TREATMENT DELIVERY: DOSE SPECIFICATION..... | 1486 |
| | | VIII. RECOMMENDATIONS FOR FUTURE CONSIDERATIONS..... | 1487 |

References – Mobile Systems

- Beck C, “On external Roentgen treatment of internal structures eventration treatment,” N. Y. Med. Journal 89, 621–622 1909.
- Abe M, “Intraoperative radiotherapy – past, present and future,” Int. J. Radiat. Oncol., Biol., Phys. 10, 1987–1990 (1984).
- Beddar AS, Krishnan S. Intraoperative radiotherapy using a mobile electron LINAC: A retroperitoneal sarcoma case. J Appl Clin Med Phys 6(3):95-107, 2005.
- Beddar AS. Stability of a mobile electron linear accelerator system for intraoperative radiation therapy. Med Phys 32(10):3128-3131, 10/2005.
- Beddar AS, Biggs PJ, Chang S, Ezzell GA, Faddegon BA, Hensley FW, Mills MD. **Intraoperative radiation therapy using mobile electron linear accelerators: Report of AAPM Radiation Therapy Committee Task Group No. 72.** Med Phys 33(5):1476-1489, 5/2006.
- Mills MD, Fajardo LC, Wilson DL, Daves JL, Spanos WJ. Commissioning of a mobile electron accelerator for intraoperative radiotherapy. J Appl Clin Med Phys. 2:121–130, 2001.
- Beddar AS, Briere TM, Ouzidane M. Intraoperative radiation therapy using a mobile electron linear accelerator: Field matching for large-field electron irradiation. Phys Med Biol 51(18):N331-337, 9/2006.
- Wootton LS, Meyer J, Kim E, Phillips M. Commissioning, clinical implementation, and performance of the Mobetron 2000 for intraoperative radiation therapy, J Appl Clin Med Phys:18: 230–242, 2017.

Additional Slides for References

Disclaimer

I am only providing these additional slides for additional information and educational purposes to complement my personnel presentation and experience using the Mobetron 1000. The slides below were obtained from IntraOp Medical to show the new version of the Mobetron: Mobetron 2000.

The author has no direct affiliation with Intraop Medical, Inc. and received no financial support for the research reported in his work in the past/present or for giving this invited presentation.

The author thanks Derek Descioli, VP Global Sales. IntraOp and Dan Goer, Co-Founder of IntraOp in 1993 and serving currently as its Chief Scientist.



IntraOp Mobetron Treatment Applications

Treatment Application Characteristics





| Tumor | Energy [MeV] | | | Applicator diameter [cm] | Bevel [deg] | | | ½ cm sized applicators |
|-------------------|--------------|-----|-----|--------------------------|-------------|-----|-----|------------------------|
| | 6 | 9 | 12 | | 0 | 15 | 30 | |
| Breast | 26% | 48% | 24% | 4 – 6 (84%) | 53% | 31% | 15% | 35% |
| Colorectal | 35% | 50% | 10% | 4 – 7 (86%) | 8% | 4% | 88% | 21% |
| Pancreas | 24% | 26% | 49% | 5 – 7 (77%) | 73% | 12% | 16% | 22% |
| Sarcoma-Extremity | 45% | 46% | 8% | 10 (30%) | 39% | 13% | 48% | 12% |
| Sarcoma-RPS | 22% | 58% | 19% | 10 (23%) | 34% | 8% | 58% | 11% |

Source: IntraOp Mobetron User Database



IntraOp Mobetron Output Stability

Mobetron Output Variation for Various Energies




Source: Beddar, et al. 2005

Mobetron Long Term Stability

| Absolute Calibration Output (cGy/MU) | | | | |
|--------------------------------------|---------|---------|---------|---------|
| Electron Energy | 4 MeV | 6 MeV | 9 MeV | 12 MeV |
| Commissioning 1999 | 0.999 | 0.997 | 0.993 | 1.001 |
| Annual Calibration 2000 | 0.987 | 0.993 | 1.004 | 1.010 |
| Percent Change (%) | - 1.2 % | - 0.4 % | + 1.1 % | + 0.9 % |

Source: Beddar, Et al. 2005

Mobetron Energy Ratio Variation



Equivalent PDD Shift

| | |
|--------|--------|
| 4 MeV | 0.6 mm |
| 6 MeV | 0.8 mm |
| 9 MeV | 0.8 mm |
| 12 MeV | 0.4 mm |



Source: Beddar, et al. 2005



Technical and Clinical Review of Mobetron Usage

Mobetron Treatments Characteristics



- ▶ Worldwide, 56% of Mobetron treatments were for breast cancer (83% in Europe, 20% in the U.S, and 33% in Asia).
- ▶ For the non-breast data main sites were: Colorectal (21%), Pancreas (15%), Extremity Sarcomas (22%) and RPS (8%).



Breast data characteristics


► For the breast data received:

- 71% of breast patients were treated as a boost (75% EU, 72% U.S., 16% Asia).
- 29% of boost treatments followed by 3 weeks EBRT.
- 48% were treated at 9 MeV; 26% at 6 MeV; and 24% at 12 MeV (2% @ 4 MeV).




Breast data characteristics (continued)

- ▶ For the breast data received:
 - ▶ 84% of patients were treated with FS between 4-6 cm.
 - ▶ 35% of patients were treated with a ½ cm sized applicator.
 - ▶ 3% of patients were treated with a FS > 7 cm.



ISIORT Breast Treatments Characteristics

- ▶ 42 Collaborating Centers
- ▶ 8,075 Breast cancer patients
 - ▶ Median age 61 years (16 – 90)
 - ▶ 81.8% T1 and 16.1% T2
 - ▶ 96.5% Ductal carcinoma; 3.5% Lobular
 - ▶ 52.2% Surgery + IORT
 - ▶ 47.8% Surgery + IORT + EBRT
 - ▶ CT in 13.2% of cases





APBI Single Fraction Treatment Clinical Results

Single Fraction Clinical Results

| Reference | Median FU (yrs) | Total # pts | # LR | LR (%) | ASTRO Suitable | | | ESTRO Good | | |
|----------------------------------|-----------------|-------------|------------|--------|----------------------|----------|--------|------------------|-----------|--------|
| | | | | | # pts ASTRO Suitable | #LR | LR (%) | # pts ESTRO Good | #LR | LR (%) |
| ELIOT ⁽¹⁾ | 5.8 | 585 | 34 | 5.8% | 135 | 2 | 1.5% | - | - | - |
| ELIOT Out Trial ^(2,3) | 3.5* | 1822 | 75 | 4.1% | 294 | 3 | 1.0% | 573 | 7 | 1.2% |
| U. of Verona ^(4,5) | 5 | 226 | 4 | 1.8% | 128 | 1 | 0.8% | 160 | 3 | 1.9% |
| Brussels ⁽⁶⁾ | 2 | 204 | 1 | 0.5% | 87 | 1 | 1.1% | 151 | 1 | 0.7% |
| Japan ⁽⁷⁾ | 6 | 32 | 0 | 0% | 16 | 0 | 0% | 28 | 0 | 0% |
| Brazil ⁽⁸⁾ | 4.3 | 152 | 5 | 3.2% | 48 | 0 | 0% | 92 | 2 | 1.9% |
| Naples ⁽⁹⁾ | 4 | 13 | 0 | 0% | 7 | 0 | 0% | 13 | 0 | 0% |
| China ⁽¹¹⁾ | 4.3 | 36 | 2 | 5.6% | 2 | 0 | 0% | 6 | 0 | 0% |
| Stanford ⁽¹⁰⁾ | 6.8 | 64 | 1 | 1.6% | - | 0 | 0% | - | 0 | 0% |
| Total | | 3134 | 122 | | 717 | 7 | | 1023 | 13 | |

Source: Horst et al., Presented at 2016 ASCO Annual Meeting

5 yr LR of <2% for Low Risk Patients

Updated ASTRO APBI Guidelines Affirm the use of Electron IORT

- ▶ 2016 ASTRO APBI Consensus Guidelines update removes experimental status for electron IORT in suitable patients
- ▶ Citing Evidence from Multivariate analysis on randomized trial with median follow up of 5.8 years
- ▶ More than 3,000 patients receiving single treatment electron IORT have been studied in the literature
- ▶ Electron IORT is the only single treatment option recognized by the updated ASTRO guidelines

ASTRO Suitable Criteria

AGE: 50 years or older

SIZE: ≤ 2 cm

NODAL STATUS: pN0

HISTOLOGY: IDC/Favorable

MARGINS: Negative (≥ 2 mm)

ESTROGEN RECEPTOR: Positive

LVSI: Negative

DCIS : Low to intermediate ≤ 2.5 cm

DCIS MARGIN: Negative (≥ 3 mm)

SYSTEMIC THERAPY: No neoadjuvant

GRADE: Any
