

Rapid Modeling of HIFU Beam Propagation through Inhomogeneous Tissues: Validation and Applications

Douglas Christensen, Scott Almquist, Sara Johnson, Christopher Dillon, Dennis Parker and Allison Payne

Utah Center for Advanced Imaging Research
University of Utah
Salt Lake City, UT 84112

UCAIR
UTAH CENTER FOR ADVANCED IMAGING RESEARCH




Outline


- Overview of Hybrid Angular Spectrum (HAS) method
- Validation:
 - Comparison of HAS to k-Wave simulations
 - With experimental MRI temperature imaging
- Applications:
 - Phase aberration correction in UofU Breast System
 - Characterization of scattering by the skull

Overview of Hybrid Angular Spectrum Method*

- Extends traditional homogeneous angular spectrum method (in spatial-frequency domain) to include 3D heterogeneous media
- Leapfrogs between the space and spatial-frequency domains (next slide)
- Employs FFT commands, so very rapid
- Assumptions: steady state conditions, linearity and compressional waves only

*U. Vyas and D. A. Christensen, "Ultrasound beam simulations in inhomogeneous tissue geometries using the hybrid angular spectrum method," IEEE Trans UFFC 59 (6), 1093-1100, June 2012.






Validation

Validation 1 - Comparison of HAS to k-Wave Simulation


- k-Wave: a well-known open-source pseudospectral k-space method
- Both simulation methods modeled pressure from a phased-array transducer into a heterogeneous breast model
- HAS pattern was compared to steady-state result from k-Wave
- k-Wave: 8 hours, 32 min; HAS: 27 sec



309 x 181 x 308 MRI-segmented model
 0.25-mm resolution 1.0 MHz
 3 tissue types :
 breast fat 1480 m/s 0.75 dB/cm
 fibroglandular 1480 m/s 0.80 dB/cm
 cancer 1560 m/s 1.15 dB/cm


Results – Comparison of HAS to k-Wave Simulation

Pressure Profile through Focus



Normalized Root-Mean-Square Deviation (NRMSD) = 2.96%
 over 10 x 10 x 20-mm volume around focus


Validation 2 – With Experimental MRI Temperature Imaging




Two-fold Motivation

1. Compare simulations with experimental temperature profiles -
 Acoustic simulations: Hybrid Angular Spectrum method (HAS)
 Thermal simulations: Finite-Difference Time-Domain (FDTD) of Pennes's Bioheat Equation (PBHE)
2. Use Monte Carlo statistical analysis to calculate expected uncertainties in the simulated temperature rise given uncertainties of input parameters -
 Tissue-specific properties vary widely in literature

Flow Diagram




Measurements of Parameters and Their Uncertainties



Results – Peak Temperature


Temperature Rise (Spatial Peak Temporal Peak)
(ave of N = 3)



- Average Normalized Temperature Difference (N = 18) = 3.3%
- Error bars represent expected temperature range due to uncertainties in simulation input parameters (Monte Carlo)

Results – Spatial Profiles

one example: 50% milk at 25 W (ave of N = 3)







- Normalized Difference in FWHM (N = 18)
Transverse = 4.7%
Longitudinal = 12.1%
- Offset in location of peak (x,y,z) = (0,0,1.7 mm)

Applications


Application 1 – Phase Aberration Correction

What is phase aberration correction?



hydrophone scan
3D-printed aberrator
phased-array transducer


Concept : Correcting Phase Aberration with HAS Simulations




Propagate beam from each element with 0° phase
Find each phase at intended focus
Invert phase on respective elements
Achieve constructive interference at intended focus

Phase Aberration Correction - in UofU Breast System

UofU Breast System - cut-away view :



HAS simulations with phase aberration* :




laterally firing phased-array transducer

*A Farrer, et al. 2016 Med Phys 43(3): 1374-1384


Motivation: Experimentally demonstrate that HAS simulations can correct for phase aberrations in experimental breast-mimicking phantoms

Experimental HAS-based Phase Aberration Correction

- Breast mimicking phantoms
 - 250-bloom gelatin with 50% milk
 $c = 1555 \text{ m/s}$, $\alpha = 0.047 \text{ Np/cm MHz}$
 similar to fibroglandular tissue
 - canola oil (generates aberrations)
 $c = 1469 \text{ m/s}$, $\alpha = 0.0125 \text{ Np/cm MHz}$
 similar to breast fat
- MRgFUS heating
 - 20 -30 s sonications, 12 -24 acoustic watts
 - 3D MRTI in oblique plane
 - with and without phase aberration correction




Results



➔ Average of 23% increase in temperature rise (13 unique sonication locations)


Application 2 – Characterization of Scattering by the Skull



1024-element phased-array transducer

- Attenuation = absorption (heating) + scattering
- Acoustic scattering in the skull is a significant portion of attenuation in transcranial HIFU treatments
- Attenuation varies from patient to patient, and current clinical images (from CT Hounsfield Units) do not separate out the scattering portion of attenuation
- Modeling scattering requires finer resolution than clinical CT, thus MicroCT


Models – Developed from Human Skull Flap



Simulating Scattering from MicroCT Models with COMSOL

- FEM iterative solver with automatic meshing
- 1 MHz, 650 kHz, 220 kHz plane wave input
- No absorption or mode conversion included (scattering only)
- Acoustic properties:
 - Bone – 2900 m/s, 1900 kg/m³
 - Porous areas – 1500 m/s, 1000 kg/m³
- Simulated 3 x 13 MicroCT models, 3 each of 13 100-HU ranges

Results – Scattering Coefficient vs. Hounsfield Units



- Scattering increases at low HU (more small pores)
- Scattering increases at higher frequencies
- Next: Experimental validation with several small skull pieces covering various HU ranges

Conclusions

Validation:

- Comparison of HAS to k-Wave breast model pressure patterns was within 3%
- Comparison of simulated to experimental temperatures was within 3.3% for peak temperature, and within 12.1% for spatial FWHM

Applications:

- Phase aberration correction is beneficial for many sonication locations in the UofU Breast System
- Map being developed of scattering coefficient to clinical CT Hounsfield Units

Acknowledgements




Focused Ultrasound Lab

- Dennis Parker, Ph.D.
- Allison Payne, Ph.D.
- Rock Hadley, Ph.D.
- Doug Christensen, Ph.D.
- Lorne Hofstetter, M.S.
- Scott Almquist, B.S.
- Chris Dillon, Ph.D.
- Henrik Odeen, Ph.D.
- Sara Johnson, B.S.
- Bryant Svedin, Ph.D.
- Hailey McLean
- Emilee Minalga, M.S.
- Robb Merrill, M.S.
- Erik Dumont, Ph.D.
- Image Guided Therapy

This work was funded by NIH R01 CA172787, R03 EB023712 & R03 EB013433.

Thank you -
Any questions?


MicroCT Models - Registered to Clinical CT Voxels



~ 600 HU
 ~ 900 HU
 ~ 1300 HU

Form three samples of
 3D MicroCT models over
 13 HU ranges


Models - Developed from a Human Skull Flap



Motivation

- Acoustic scattering in the skull is a significant portion of attenuation in transcranial treatments
- Attenuation = absorption + scattering
- Attenuation varies from patient to patient, and current clinical images (from CT Hounsfield Units) do not separate out the scattering portion of attenuation.
- Modeling scattering requires finer resolution MicroCT

Results – Scattering Coefficient vs. HU



- Scattering increases at low HU
- absorption = attenuation - scattering
- At low HU, need to differentiate between clusters of small scatterers and partial volume effect

[1] S. Pichardo, V. W. Sin, and K. Hynynen, *Phys Med Biol* 56/1, Jan 2011.
