

Clinical Implementation of Pencil Beam Scanning (PBS) Proton Therapy

Commissioning PBS dose calculation algorithms and understanding their limitations



X. Ronald Zhu, PhD
University of Texas
MD Anderson Cancer Center
July 30, 2018



Particles contributing to doses

- Primary protons
 - Elastic interactions with electrons
 - Elastic proton-nucleus scattering
- Secondary particles
 - Non-elastic nuclear interactions
 - Secondary protons and other fragments (deuterons, tritons, alphas, neutrons, etc.)




Dose Algorithms

- Monte Carlo Simulation
 - Becoming available for clinical uses in commercial TPS
- Analytical calculation – pencil beam algorithms

$$D(x,y,z) = I[d(z)] \times LAT[x,y,d(z)]$$

■ $I(d)$ – integral depth dose

■ $LAT(x,y,d)$ – lateral dose profile




depth = d

$$I(d) = \int_0^d I(z) dz$$



Integral Depth Dose



Pedroni et al. PMB 2005

- (a) Total dose
- (b) Beam flux attenuation
- (c) Dose locally deposited
- (d) Dose from "beam halo" or low dose envelope
- (e) Dose from primary protons
- (f) Beam width σ broadening due to MCS in water




Lateral Dose Profile

- Multiple Coulomb Scattering (MCS)
 - In the patient
 - In the range shifter propagating through the air gap to the patient
 - Some devices in beamline
- Nuclear interaction
 - beam "halo" due to large angle inelastic nuclear fragments




Sawakuchi et al. PMB 2010



Low dose envelope

- High E – nuclear fragments
 - Depends on energy & depth
- Low E – MCS
 - In the devices and phantom/patient
- Small, but can be significant when thousands of spots are used



Sawakuchi et al. PMB 2010

Low dose envelope

- Spot isodoses for different energies at different depths

Lin et al. PMB 2013

Low dose envelope

- Methods of measurements:
 - Concentric square frames: Pedroni et al. PMB 2005 - protons
 - Field size factors: Sawakuchi et al. PMB 2010 - protons, Inaniwa et al. Med Phys 2009 - carbons
 - Concentric circles: Clisie et al. PMB, 2012 - protons

Sawakuchi et al PMB 2010

Commissioning - Example

Zhu et al. Med Phys 2013

- Total dose: Fluence x Beamlet dose

$$D(x, y, z) = \sum_{E_i} \left\{ \sum_{\text{Beamlet } j} [\Phi_{E_i}(x_j, y_j, z) D_{E_i}^{\text{Beamlet}}(x-x_j, y-y_j, d(z))] \right\}$$

- Beamlet dose: IDD's x Kernel

$$D_{E_i}^{\text{Beamlet}}(r, d(z)) = \frac{1}{\rho_{H,O}} [S(d) \times K(r, d)]$$

$$D_{E_i}^{\text{Beamlet}}(r, d(z)) = \frac{1}{\rho_{H,O}} [S_{pp}(d) K_{MCS,prim}(r, d) + S_{sp}(d) K_{MCS,sec}(r, d)]$$

pp - primary photons	$K_{MCS,prim}$ - MCS, Moliere theory, 2 Gaussians
sp - secondary particles	$K_{MCS,sec}$ - secondary particles, nuclear interaction

Input Data Requirements by Treatment Planning Systems

- In air profiles:
 - At 3 to 5 different positions from isocenter (e.g., ±200, ± 00, and ±0 mm) for every 10-20 MeV in both directions.
 - If a range shifting device is used, 2~3 complete data sets for 2~3 different thicknesses.



Input Data Requirements by the Treatment Planning System

- Integrated depth doses (IDDs):
 - Depth dose to be measured with a large p-p chamber with the radius R .

$$R \geq 3\sigma_{spot} = \sqrt{\sigma_{fluence}^2 + 2(0.0307 \times Range)^2}$$

- IDD are in unit of Gy•mm²/MU or Gy•mm²/Gp (Gp = 10⁹ ,giga protons)
- The p-p chamber might be too "small" and requires a correction for dose deposited outside of the p-p chamber

Claise et al PMB 2012



Correction factors for IDD

- MC simulation to determine CFs

$$CF(E, d = 2cm) = \frac{IDD_{MC}(E, d = 2cm; r_2)}{IDD_{MC}(E, d = 2cm; r_1)}$$

- Convert MC IDD to Gy•mm²/MU or Gp


$$IDD(E, d; r_2) = \frac{IDD_{MC}(E, d; r_2)}{IDD_{MC}(E, d = 2cm; r_2)} \cdot IDD_{meas}(E, d = 2cm; r_1) \cdot CF(E, d = 2cm)$$

Normalized MC IDD at 2 cm depth


Measured IDD at 2 cm depth

Zhu et al Med Phys 2013






Difference in field size factors



20 mono-energetic fields
 Avg ± Stdev = 0.2% ± 0.7% (-1.7% to 2.1%)




• In air lateral profiles: Comparison of Measurements, Monte Carlo and TPS calculated by single- and double Gaussian fluence model

Zhu et al Med Phys 2013


Fluence Model with Gaussians

- Colleagues at Mayo Clinic developed a in-house method to determine Gaussian fluence parameters



Shen et al Med Phys 2016


Dose at the center of SOBPs



The SOBP widths: 2 to 24 cm.
 Avg ± Stdev = 0.0% ± 0.6% (range, -1.9% to 1.2%)
 Zhu et al Med Phys 2013



Comparison of depth doses




Field size of 10 cm x 10 cm,
 SOBP width of 4 or 10 cm.

Zhu et al Med Phys 2013




Comparison lateral dose profiles





Zhu et al Med Phys 2013


Prostate Patient



Base of skull (Chordoma) patient







It is not perfect

- Commissioning will not exhaustively test all clinical scenarios
- Patient QA including dose measurements is desirable

Table 2: Summary of the gamma index passing percentages from the patient specific quality assurance of 2,187 treatment fields.

(a)	[2%, 2-mm]		[3%, 3-mm]		(b)	[2%, 2mm]		[3%, 3-mm]	
	[2%, 2-mm]	[3%, 3-mm]	[2%, 2mm]	[3%, 3-mm]		[2%, 2mm]	[3%, 3-mm]	[2%, 2mm]	[3%, 3-mm]
Overall	85.3±0.8%	96.2±0.3%	81.5±1.5%	95.2±0.8%	SFD	81.5±1.5%	95.2±0.8%	81.5±1.5%	95.2±0.8%
CNS	85.9±1.0%	95.0±1.2%	83.3±1.1%	95.6±0.6%	MPO	83.3±1.1%	95.6±0.6%	83.3±1.1%	95.6±0.6%
HN	82.7±1.2%	94.9±0.7%	81.9±1.0%	94.8±0.6%	RS	81.9±1.0%	94.8±0.6%	81.9±1.0%	94.8±0.6%
Prostate	100.0%	100.0%	86.1±2.0%	99.0±0.0%	NRS	86.1±2.0%	99.0±0.0%	86.1±2.0%	99.0±0.0%
Therapeutic	80.1±1.0%	97.2±0.6%							

Mackin et al. IJPT 2014

Summary

- Accurately modeling the low dose envelope is one of the most important elements during PBS commissioning
- Analytical dose models have limitations
- Patient specific QA should include dose measurements to continue validating the dose model.
- Better dose calculation methods such as Monte Carlo simulation are becoming available in commercial TPS
- Methods presented here could be used for the basic validations of Monte Carlo dose calculation as well
- Inhomogeneity phantoms (e.g., IROC phantoms) should be used for end-to-end validations of imaging, planning and delivery.

Acknowledgements

Gabriel Sawakuchi, PhD
 Falk Poenisch, PhD
 Uwe Titt, PhD
 Xiaodong Zhang, PhD
 Narayan Sahoo, PhD
 Yupeng Li, M.S.
 Heng Li, PhD
 Jim, Lii, MS
 Radhe Mohan, PhD
 Michael Gillin, PhD
 Many others



Thank you!



Core, Halo, Aura and Spray

- Gottschalk et al have a more precise definition for the dose distribution of a proton beam stopping in water - Basic physics component and others from beam contaminations
 - Core – for the primary beam
 - Halo – for the low dose region from charged secondaries
 - Aura – for the low dose region from neutrals
 - Spray – for beam contamination

*Gottschalk et al PMB 2015


Basic Information about Bragg Peak Chamber

- Nominal sensitive volume: 10.5 cm³.
- Sensitive volume: r = 40.8 mm, t = 2 mm.
- Nominal response: 325 nC/Gy.
- Reference point 3.5 mm front chamber surface.
- Entrance window: 3.47 mm PMMA.
- WET window: 4 mm.
- $N_{D,wk_p} = (3.181 \pm 0.023) \times 10^6$ Gy/C*
 - Average 3 inter-comparison




*Gillln et al, Med Phys 2010

Brain - SFIB



- 11 year old boy
- Glioma of tectal plate
- Single field (Simultaneous) integrated boost (SFIB)



SOBP does not have to be flat

