




Functional Image-guided Thoracic and Hepatic Radiation Therapy: synergies between nuclear medicine and radiation oncology

Stephen R. Bowen, PhD, DABR
Associate Professor, Radiation Oncology & Radiology (joint)
July 15, 2019


UW Medicine

Disclosures




- No conflicts of interest to disclose
- Funding: NIH / NCI R01CA204301


Synergy: Nuc Med Imaging for Precision Rad Onc




- Key premise is that therapeutic ratio (efficacy / toxicity) can be individually optimized
- Nuc med imaging as a biomarker to risk stratify rad onc patients
 - Select therapies, dosing, fractionation
 - Define target / organs-at-risk volumes + dose objectives
- Nuc med imaging as a biomarker to evaluate individual RT response
 - Adapt therapies, dosing, fractionation
 - Adapt target / organs-at-risk volumes + dose objectives
- Nuc med imaging to spatially optimize RT dose distribution
- Applications of nuc med imaging to radiation therapy for thoracic / hepatic cancers











FDG PET Predicts RT Failure Patterns

- FDG PET as spatial map of local treatment failure risk distribution
Bougjar Radiother Oncol 2010

NM Imaging to Optimize RT Spatial Dose Distribution

- Dose painting by contours / Subvolume boosting
 - Biological target volumes (BTV) for uniform dose escalation
(Lung 2009, Tome 2003, Madani 2009)
- Dose painting by numbers accounts for intratumoral variations in response to therapy
(Brahme 2003, Bentzen 2005, Le 2007, Stewart 2007)

Meijer et al. 2011 Radiother Oncol

FDG PET-guided Lung Dose Painting Clinical Trials

- NKI-Maastricht: average dose increase to PET avid areas up to 85+ Gy in 24 fractions
van Elmpt Radiother Oncol 2012
- RTOG 1106: dose escalation based on mid Tx PET (up to 80.4 Gy in 30 fractions)
Kang et al. NRG



Nuc Med Lung Functional Imaging

- Perfusion
 - ^{99m}Tc -MAA SPECT/CT
 - ^{68}Ga -MAA PET/CT
- Ventilation
 - ^{99m}Tc -DTPA SPECT/CT
 - ^{99m}Tc -Technegas SPECT/CT
 - ^{68}Ga -Galligas PET/CT
- Inflammation
 - FDG PET/CT

Siva et al RAO 115(2): 157-162

Petit, DeRuyscher IJROBP 2011

Functional Lung Dosimetry for Toxicity Prediction

- Patients with similar clinical characteristics & anatomic dosimetry (MLD)
 - Upper lobe primary tumors
 - Conventionally fractionated chemoRT
 - Anatomic mean lung dose 16.6 Gy (top) vs. 16.4 Gy (bottom)
- Different functional lung dosimetry (pMLD)
 - Top (pneumonitis): perfused mean lung dose 24.2 Gy
 - Bottom (no pneumonitis): perfused mean lung dose 8.6 Gy
- Combined MLD and pMLD best predict for Grade 2+ pneumonitis in initial (AUC = 0.92) and expanded patient cohorts (AUC = 0.94)

Dahmi et al. *Strahlenther Onkol* 2017
Lee et al. *IJROBP* 2018

Functional Lung RT Dose-response Modeling

- Perfusion decline spatially correlates with increased inflammation in high dose regions
- Baseline high perfusion regions have steeper (more sensitive) dose-response curves

Thomas et al. *BJR* (under review)
Owen et al. *IJROBP* 2018

Nuc Med Liver Functional Imaging

<ul style="list-style-type: none"> PET/CT <ul style="list-style-type: none"> - FDG - FDGalactose SPECT/CT <ul style="list-style-type: none"> - GSA - HIDA - Sulfur colloid 	^{99m} Tc-Sulfur colloid Homogeneous Regional Uptake	High Global Uptake Untreated CP-A6	Low Global Uptake Untreated CP-B7
	Heterogeneous Regional Uptake	Treated CP-A5 Chemoembolization	Untreated CP-B7

Bowen et al. EJNMI Res 2016

Functional Liver Dosimetry for Outcome Stratification

Group	Median OS (days)	1 year OS	Log rank p
CP A	-	50%	-
CP B/C	202	39%	0.002

Group	Median OS (days)	1 year OS	Log rank p
Low Risk	-	100%	-
Intermediate Risk	525	68%	0.0002
High Risk	202	30%	-


Schaub et al. UROBP 2018

Functional Liver RT Dose-response Modeling

Price et al. UROBP 2018




Functional Liver Avoidance Clinical Trials



Differential Hepatic Avoidance Radiation Therapy (DHART) Planning (U Wash)


Functional liver image-guided hepatic therapy (FLIGHT) trial (Indiana)

- HIDA SPECT/CT
- Maximize functional liver volume receiving < 15 Gy
- Inter-patient variability in benefit




Bowen et al. *Radiother Oncol* 2015

Functional Lung Avoidance Clinical Trials






Perfused Lung Sparing Trial (Peter Mac)



Siva et al *R&O* 115(2): 157-162

- IMRT plans adapted to avoid perfused lung on 4D MAA PET/CT
- Ongoing functional lung avoidance trials using NM imaging
 - 4D MAA PET/CT (Peter MacCallum)
 - MAA SPECT/CT (U. Washington, U. Michigan)

Functional Lung Avoidance & Response-adaptive Escalation (FLARE) RT

Leo et al. *Med Phys* 2017



FLARE-RT Phase II Trial (NCT02773238, R01CA204301, PI: Bowen, Zeng)

- All patients get functional lung avoidance RT: potential quality of life benefit
- Only high local failure risk patients get FDG PET-guided dose escalation: potential survival benefit
- 1st endpoint: 2 yr overall survival vs. RTOG 0617 60 Gy arm
- 2nd endpoints: 1 yr local control, grade 2+ pneumonitis incidence vs historical rates

FLARE-RT Mid-Tx PET Response Assessment

Median 10 days, Range < 1 month

Median 24 Gy in 12 fx
 Range 20-26 Gy in 10-13 fx

Baseline Pre RT

Week 3 Mid RT

FLARE008 Responder

FLARE010 Non-Responder

- PET Responders: Δ SUV_{peak} decrease 40% (32-48%), Δ MTV decrease 40% (32-62%)
- PET Non-responders: Δ SUV_{peak} decrease 13% (8-23%), Δ MTV decrease 9% (-5-21%)

FLARE-RT Plan Adaptation Decisions

PET Responder

PET Non-responder

Pre-RT

3wk Mid-RT

FDG PET/CT

- PET Responder: required anatomic adaptation to 60 Gy in 30 fx
- PET Non-responder: required functional adaptation to 74 Gy in 30 fx



FLARE-RT Boost to PET Non-responders

- 13 / 34 (38%) have received FLARE RT boost
 - 3 week PET non-responders & no evidence distant mets
 - 74 Gy to PTV_{max} [90+ Gy to SUV_{peak}]
- Well tolerated
 - 0 Grade 3+ esophagitis
 - 2 Grade 2+ pneumonitis

Future: Multiscale Imaging Response Prediction

- Multiscale imaging response prediction uses information from individual image voxels combined with regional, tumor, and patient factors
- Voxel Forecast Tool:** custom generalized least squares (GLS) regression to predict tumor voxel response on mid-RT PET
 - Matérn model variogram to account for spatially correlated voxel data
 - Jack-knife bias-corrected estimator validated on simulations of known voxel response patterns

Hippe AAPM 2018 Best in Physics Bowen et al. Clin Cancer Res 2019

Voxel Forecast Tool: PET Voxel Response Prediction for Rad Onc Decision Support

PET Responder
 NED 491 days
 14% under-responding tumor voxels
 MAE = 1.1 SUV

PET Non-responder
 Died 323 days
 94% under-responding tumor voxels
 MAE = 3.3 SUV

Hippe AAPM 2018 Best in Physics Bowen et al. Clin Cancer Res 2019



Future: Personalizing NM+RO for HCC management

Favorable planned dosimetry **Unfavorable planned dosimetry**

- **Favorable planned dosimetry:** personalize ^{90}Y -microsphere prescriptions to achieve isotoxicity in liver → verify on ^{90}Y -SPECT/PET dosimetry
- **Unfavorable planned dosimetry:** modify ^{90}Y -microsphere injection site → boost any remaining cold ^{90}Y -SPECT/PET dosimetry regions with SBRT

Summary: Synergies between Nuc Med and Rad Onc

- Nuc Med imaging to personalize Rad Onc
 - Risk stratify patients
 - Define targets and functional normal tissues
 - Optimize prescriptions and planned radiation dose distributions
 - Assess early response for adaptive therapy
- **Future:** combined NM image-guided RT with NM therapies
 - RT + targeted radionuclide therapy (TRT) + immunotherapy (IO)
 - NM imaging for targeting / avoidance / dosimetry / verification / response
 - Radiomics / machine learning of NM imaging to personalize NM + RO Tx

Acknowledgments

Financial Support NIH/NCI R01CA204301	Imaging Research Lab Paul Kinahan Robert Myczka Adam Alessio Larry Pierce Darrin Byrd Rebecca Cristofel Nathan Bell	SCCA Proton Therapy Tony Wong Jay Sani Dominic Maes Alex Egan Rajesh Reamji Erick Leuro Zander Nevitt Dosi, Sim, RT staff
Radiation Oncology Jing Zeng Jim Ajayaramhanarax Ramesh Rengan Nina Mayr Shilpa Patel George Sandison Matt Nyflot Charles Bloch Juergen Meyer Mark Phillips Rob Stewart Eric Ford Priya Vissamraju Dosi, Sim, RT staff	Radiology Hubert Vesselle Manuela Mettessan William Yuh Dan Hippie Vivi Wu Missy Warner Fam Pham NM Tech staff	Trainees Kevin Horn Toby Chapman Stephanie Schaub Ryan Price Eunsin Lee Nina Dharmi Howard Lee Sally Duan Hannah Thomas



