

Joint AAPM/SNMMI Symposium
Nuclear Medicine Theranostics and Functional Image-Guided
Radiation Therapy for Precision Oncology

Image-guided ⁹⁰Y-radionuclide treatment planning, delivery, and verification for hepatic cancers

S. Cheenu Kappadath, PhD
Associate Professor, Department of Imaging Physics
<http://www.mdanderson.org/kappadath-lab>



1

Disclosures

- Research grants: BTG International, GE Healthcare
- Consultant: BTG International, Terumo Medical Systems, ABK Biomedical

S. Cheenu Kappadath, PhD


AAPM 2019 2

⁹⁰Y-microsphere Radioembolization, or Selective Internal Radiation Therapy (SIRT)

- Trans-arterial delivery of radioactive ⁹⁰Y-labeled microspheres via a catheter directly at disease sites (targeted infusion)
- Microspheres (20-30 μm) trapped in tumor capillary vessels due to their embolic size and targeted delivery



© SIRTEX



© TheraSphere



- β emissions from trapped ⁹⁰Y-microspheres deliver radiation dose to proximal tissue (tumors) while sparing distal (normal liver) tissue → max range of 10 mm


S. Cheenu Kappadath, PhD

(Murthy et al, Sem IR, 2008; Sarfaraz et al, Med Phys, 2004)

AAPM 2019 3

Rationale for liver-directed therapy

- Normal liver blood flow
 - 75% portal vein
 - 25% hepatic artery
- Hepatic neoplasm, >3mm metastases
 - 80- 100% supply from hepatic artery
 - Greater vascular density in neoplasm
- Primary site of disease in HCC and ICC
- Dominant organ of metastases in CRC and NET
- Resection improves survival HCC, CRC, NET
 - Control of liver disease should increase survival





Indications

Non surgical candidate
 Not amenable to ablative therapy
 Bridge to transplant or resection
 Palliative for liver-only or liver-dominant disease

S. Cheenu Kappadath, PhD

AAPM 2019 4

90Y-SIRT Workflow





- Tumors fed by hepatic artery

SIR-Spheres, Sirtex
 TheraSphere, BTG
 QuiremSpheres, Terumo*
 Eye90, ABK Biomed*
 *not FDA-approved

S. Cheenu Kappadath, PhD

AAPM 2019 5

90Y-SIRT Workflow




99mTc-MAA


S. Cheenu Kappadath, PhD

AAPM 2019 6

90Y-SIRT Workflow




S. Cheenu Kappadath, PhD




AAPM 2019 7

90Y-SIRT Workflow



S. Cheenu Kappadath, PhD



AAPM 2019 8

SIRT is "IGRT" but focused on Safety


- Lung Doses:
 - Arterio-venous shunting in neoplastic vasculature
 - Prevention of radiation pneumonitis
 - Mean lung dose <30 Gy per treatment
- Liver Doses:
 - Maintain upper limit to mean to dose total liver
 - SirSpheres < 80 Gy & TheraSphere < 80-150 Gy
 - Assume uniform uptake in tumor and normal liver
- Major Challenge for SIRT: Current therapy planning not designed to deliver specific dose to target lesions
 - Accurate dosimetry models not routinely used
 - Tumor dose-response and toxicities not well established

S. Cheenu Kappadath, PhD

AAPM 2019 9

90Y-SIRT should be based on dosimetry

- Radiation is the actuator of therapeutic effect not embolization
- Intent: Curative vs Palliative ← disease stage, prior treatments
- Organs at Risk (OAR) in SIRT: Lung and Normal Liver
- Aim to increase therapeutic ratio ← max tumor dose yet acceptable OAR dose



S. Cheenu Kappadath, PhD


AAPM 2019 10

The SIRT Dosimetry Conundrum

Efficacy is predicated by good match between planned and actual radiation dose distribution

PROSPECTIVE


- PLANNING: MAA is not a consistently reliable predictor of microsphere distribution (dose)
- TARGET: Doses necessary for tumor response not fully established (recent results are promising)



S. Cheenu Kappadath, PhD

AAPM 2019 11


Does planning ^{99m}Tc-MAA represent ⁹⁰Y microspheres distribution after therapy?



S. Cheenu Kappadath, PhD

AAPM 2019 12


Package Insert: Tumor Dose Response



S. Cheenu Kappadath, PhD

AAPM 2019 13


Spatial Representation of SIRT Dosimetry Models



S. Cheenu Kappadath, PhD

AAPM 2019 14


Voxel Dosimetry: Cumulative DVH



S. Cheenu Kappadath, PhD

AAPM 2019 15

Different Dosimetry Models On the Same Patients With Matched VOIs Result in Large Differences in Absorbed Dose Estimates



S. Cheenu Kappadath, PhD


(Mikell et al, Int J Rad Onc Phys Bio, 2016)

AAPM 2019 16



Post-therapy ⁹⁰Y-PET/CT

- ⁹⁰Y also emits β+ (E_{max} ~ 800 keV) with BR = 32 x 10⁻⁶
 - Internal pair-production in the 0⁺-0⁻ transition of ⁹⁰Zr from ⁹⁰Y decay (first works circa 1955; Selwyn et al, App Rad Iso 65, 2007)
- First clinical ⁹⁰Y PET image published in 2009 (30 min/bed)
- Quantitative accuracy depends on coincidence counts, system hardware, acquisition & reconstruction parameters
 - Background activity errors ~10%
 - Spheres < 2.5 cm underestimated by ~20%




S. Cheenu Kappadath, PhD

AAPM 2019 17




Post-therapy ⁹⁰Y-SPECT/CT

- Standardized acquisition & reconstruction yet to be established
 - Monte-Carlo based techniques excellent image quality
 - Practical approaches can also provides clinically meaningful evaluation of in vivo ⁹⁰Y distribution
 - Partial volume errors for tumors < 3-4 cm
- Quantitative ⁹⁰Y SPECT/CT
 - Self-calibration but VOI choice is critical
 - Calibration errors vary 25%-70%, therefore
 - Consistent acquisition & reconstruction parameters is paramount




S. Cheenu Kappadath, PhD

AAPM 2019 18



90Y-SIRT Voxel Dosimetry

- Start with quantitative 90Y SPECT/CT or 90Y PET/CT
- Voxel dosimetry calculations
 - Monte Carlo transport = Local Deposition = Soft-tissue kernel (liver only)



(Mikell et al, EINMMI Physics 2, 2015)

S. Cheenu Kappadath, PhD

AAPM 2019 19

Tumor Dose Response Study

- Single-institutional retrospective study (n=34) (Kappadath et al, Int J Rad Onc Phys Bio, 2018)
 - 53 HCC tumors from 34 90Y glass-microsphere treatments
- Tumors and liver lobes segmented by Interventional Radiologist
 - Diagnostic CT or MR images co-registered with 90Y SPECT/CT
 - Tumors diameters > 2.5 cm; Maximum of 3 tumors per patient
- Calculate voxel-level absorbed dose (AD) and biological effective dose (BED)
 - Activity & Tissue distributions from quantitative 90Y SPECT/CT
 - Local dose deposition
- Association of tumor response with AD and BED evaluated

$$BED_{ijk} = D_{ijk} + \frac{D_{ijk}^2}{(\alpha/\beta)_{ijk}} \left(\frac{\lambda}{\mu_{ijk} + \lambda} \right)$$

(Dale et al, Can Bio Radiopharm 20, 2005)

	HCC tumor	normal liver
λ eff. decay rate (h ⁻¹)	0.0108	0.0108
α/β radiosensitivity (Gy)	15	2.5
T _{1/2} repair (h)	1.5	2.5

(Tai et al, IJROBP 70, 2008)

S. Cheenu Kappadath, PhD

AAPM 2019 20

Response Metric?

- Tumor response evaluated by IR on follow-up CT or MR at 3 and 6 months using WHO, RECIST, and mRECIST
 - Appropriate choice of response metric is essential

WHO	• No Correlation
RECIST	• No Correlation
mRECIST	• Significant Correlation (p<0.01) • Most Significant: AD20-70 & BED30-70


mRECIST
 ORR = 57%
 DCR = 96%

(Kappadath et al, Int J Rad Onc Phys Bio, 2018)

S. Cheenu Kappadath, PhD

AAPM 2019 21

TCP curves for HCC following SIRT



S. Cheenu Kappadath, PhD

(Kappadath et al, Int J Rad Onc Phys Bio, 2018)

AAPM 2019 22

HCC Tumor Response Dose Thresholds

Study	No. of Patients, Tumors	Device Used	Voxel Dose Image	Dosimetry Model	Threshold Dose
Strigari 2010 ¹	73 Patients >73 Tumors	SIR-Spheres	⁹⁰ Y SPECT	Voxel	Dmean > 97 Gy TCP 50%
Chan 2018 ²	27 Patients 38 Tumors	TheraSphere	⁹⁰ Y PET/CT	Voxel	Dmean > 200 Gy
Kappadath 2018 ³	34 Patients 53 Tumors	TheraSphere	⁹⁰ Y SPECT/CT	Voxel	Dmean > 160 Gy TCP 50%
Garin 2013 ⁴	71 Patients >71 Tumors	TheraSphere	^{99m} Tc-MAA SPECT/CT	Partition	Dmean > 205 Gy
Chiesa 2015 ⁵	52 Patients 60 Tumors	TheraSphere	^{99m} Tc-MAA SPECT	Voxel	Dmean > 390 Gy TCP 50%


Patient selection (BCLC stage), treatment volume (whole liver vs lobar vs segmental) have large effects on patient response

1. Strigari et al. *INM* 51, 2010;
2. Chan et al. *UROBP* 101, 2018;
3. Kappadath et al. *UROBP* 102, 2018;
4. Garin et al. *EJNMMI* 40, 2013;
5. Chiesa et al. *EJNMMI* 42, 2015.

S. Cheenu Kappadath, PhD

AAPM 2019 23

Many Confounding Factors for Dosimetry



S. Cheenu Kappadath, PhD

AAPM 2019 24

Reporting of dose and dose-response

- Radioembolization community needs to be more specific when reporting dosimetry
 - Dose (e.g., 160 Gy, 60 Gy)
 - Methodology (e.g., Voxel dosimetry, Partition model)
 - Device (SIR-Spheres, TheraSphere)
 - Disease (e.g., HCC, mCRC, mNET)
- Estimate of dose deposited depends on model used
 - 100 Gy MIRD \neq 100 Gy Partition \neq 100 Gy Voxel**
- Biological effect of dose depends on device properties
 - 100 Gy SIR-Spheres \neq 100 Gy TheraSphere**


S. Cheenu Kappadath, PhD

AAPM 2019 25



MDACC HCC OS Study (n=181)

- Median OS: 13.4 months (95% CI 9.7-17.2)
- Stratification: Tumor burden (<50% or >50%) & Aggressive disease features (Y or N)





S. Cheenu Kappadath, PhD

(Teyateeti et al, SNMMI 2019)

AAPM 2019 26



90Y-SIRT + SBRT



S. Cheenu Kappadath, PhD

AAPM 2019 27



Looking Forward

- Time is right to focus on the personalized treatment plan
 - Better knowledge on tumor dose and dose response
 - Better understanding of errors in dose quantification
- Opportunities for improvements in SIRT
 - 1. Need for standardization and consistency in practice
 - 2. Need to be more descriptive when reporting dosimetry
 - 3. Improved dose response models for OAR are needed
- Improve understanding on how to incorporate SIRT as part of combination treatments
 - SBRT, proton therapy, systemic, immunotherapy

Growing Value

Cultivating Safety

Building Bridges
