Myocardial Elastography and 3D Electromechanical Wave Imaging of the Heart

Julien Grondin Department of Radiology, Columbia University, New York, NY

61st AAPM Annual Meeting and Exhibition, July 14-18, San Antonio, TX Session: Advanced Ultrasound Imaging in Clinical Applications

Coronary artery disease (CAD)

- Atherosclerotic plaques build up in the coronary arteries have
- > Can lead to ischemia and infarction

• Methods to diagnose CAD in clinic: Conventional coronary angiography: Identify and measure the degree of stenosis in a coronary artery

- Catheter insertion
 Fluoroscopy

d d

Nuclear stress test: Assess myocardial perfusion

Radioactive tracer injection
 Treadmill or pharmacological agent

Treadmill or pharmacological agent
 Qualitative

Echo stress test: Assess wall motion abnormalities

ly et al, *J Nucl Med*, 2008

¹Roth et al., N Engl J Med, 2015



Radiofrequency (RF)-based speckle tracking performs better than envelope-based tracking for small tissue deformation as they contain signal phase information^{1,2} ¹Ma et al., Ultrasonics, 2013 ²Li et al., Ultrasound Med Biol, 2015

2-D Myocardial Elastography (ME)

- Tracking of Radiofrequency (RF) signals¹ in axial and lateral direction²
- Transmural myocardial strain imaging
- Angle-independent
- Multi-sector and ecg-gating
- Validation against tagged MRI (tMRI) in humans in vivo³



 Associated with reduced blood flow in Left Anterior Descending artery (LAD) in canines in vivo⁴



Can be performed during free breathing in a single heart cycle
 Kondagou et al. Ultrasound Med Biol, 2002
 Lee et al., Pfrys Med Biol, 2011
 Lee et al., Pfrys Med Biol, 2011

Objectives

- Assess the performance of single heartbeat Myocardial Elastography to distinguish normal vs ischemic patients
- Investigate the sensitivity of Myocardial Elastography to the territories perfused by the coronary arteries



Myocardial Elastography: Patient population







RF signals at time t





Konofagou et al., Ultrasound Med Biol, 1998

RF signals at time t+1

· Motion estimation rate: 500 Hz

· 10:1 linear interpolation of RF lines in the lateral direction to improve lateral displacement

· 1-D cross-correlation1 (window size: 5.9 mm, 90% overlap)

1-D kernel in a 2-D search²

¹Luo et al., IEEE TUFFC, 2010



¹Luo et al., IEEE TUFFC, 2008

²Luo et al., IEEE TUFFC, 2004

Mean radial strain in each territory

• Computation of mean end-systolic radial strain ($\vec{\epsilon}$) in all the territories



LV: Left ventricle RV: Right ventricle LAD: Left Anterior Descending artery LCX: Left Circumflex artery RCA: Right Coronary Artery T: Total cross-section

 Computation of mean end-systolic radial strain in each territory perfused by a coronary artery











Grondin et al., Ultrasound Med Biol, 2017





Conclusion

- Myocardial Elastography can differentiate normal subjects from CAD patients determined by angiography or nuclear stress test
- Radial strains are lower in territories perfused by occluded arteries or with perfusion defect than in normal territories

Ongoing work:

- Improve quality of strain estimation using coherent compounding^{1,2}
- Validate ME against nuclear imaging (SPECT and PET) in a large cohort of patients



300

¹Sayseng et al., IEEE TUFFC, 2018









Objectives

<u>Hypothesis:</u> Local onset of myocardial shortening imaged by ultrasound is caused by local electrical activation

- Investigate the relationship between the electromechanical and the electrical activation of the canine heart *in vivo*
- Demonstrate that 3D-rendered EWI can predict arrhythmia origin location















Methods: Validation against electroanatomical mapping



















3D-rendered EWI can localized accessory pathways

Double-blinded study in a cohort of 15 WPW pediatric patients for the first time (1 patient excluded for poor echocardiographic windows)



• EWI accurately predicted the AP location before ablation in 100% cases

• 12-lead ECG analysis correctly predicted 78.6% (11/14) of AP locations

Melki et al., JACC Clinical EP, 2019

Conclusions

- ٠ Electrical and electromechanical activation are well correlated $(R^2 = 0.64 - 0.82)$
- EWI can identify PVC and quantify response to CRT •
- · 3D-rendered EWI accurately predicted the AP location before ablation with higher accuracy (100%) than 12-lead ECG (78.6%)

Ongoing and future works:

- Real-time implementation of EWI

- Treatment planning and assessment Single-heartbeat 3D EWI of the full heart



2D matrix array : 32x32 elements

Acknowledgements

National Institutes of Health:

- ≻ R01-EB006042
- > R01-HL114358
- > R01-HL140646
- > Wallace H. Coulter Foundation

Contributing authors:

Ethan Bunting, PhD Alexandre Costet, PhD Alexandre Costet, PhD Jose Dizon, MD Alok Gambhir, MD, PhD Hasan Garan, MD Christopher S. Grubb, BS Elisa Konofagou, PhD Leonardo Liberman, MD Lea Melki, MS

Pierre Nauleau, PhD James Peacock, MD Vincent Sayseng, MS Eric Silver, MD Marc Waase, MD, PhD Elaine Wan, MD Daniel Wang, MD Rachel Weber, RDCS, RVT







