

Commissioning and clinical implementation of an MRI brachytherapy suite

Ananth Ravi, PhD, MCCPM
Medical Physicist, Sunnybrook Odette Cancer Centre
Clinical Operations Lead, Brachytherapy
Associate Professor, Department of Radiation Oncology, UofT



Medicine
UNIVERSITY OF TORONTO



Sunnybrook
HEALTH SCIENCES CENTRE



Disclosure

I am a consultant to MOLLI Surgical



Outline

- Value of MR imaging in Brachytherapy
- MR brachytherapy workflows and considerations
 - MR imaging + registration
 - Intraoperative MR imaging
- Resource requirements
- Commissioning/QA





MRI in Brachytherapy




MRI-based Brachytherapy

- MRI allows us to see applicators, OAR and tumor
 - Safe dose escalation and OAR avoidance is now possible
 - Dose is patient specific and not applicator specific



CT



MRI



Impact of including MR in brachytherapy

Gyne

- Overall Survival

- Benefit of 10% OS

- retroEMBRACE compared with historical controls

- Toxicities

- Reduction in 3 – 6 % per organ

Prostate

- Volume

- CT volumes are 16% larger than MR

- TRUS volumes are 10% smaller than MR

- Dosimetry

- MR limits dose spill at Apex and Base

- Can control dose to urinary sphincter

Potter/Tanderup et al, EMBRACE Review, ctRO 75 (2018)

Sturdza et al, retroEMBRACE, Radiother Oncol 120 (2016) 428-433

Smith et al, IJROBP 64 (2007) 1238-1247

Takiar et al, Brachytherapy, 13 (2014) 68-74




Medicine
UNIVERSITY OF TORONTO


 **Sunnybrook**
HEALTH SCIENCES CENTRE




MRI in clinical workflow



MR + CT hybrid process with rigid registration





CT



MR



Contour based deformable registration




Structure	DIL			
Metric	ROs	MR2US	Rigid	B-Spline
DSC	0.80 ± 0.10	0.80 ± 0.13	0.65 ± 0.20	0.51 ± 0.30
MDA (mm)	1.24 ± 0.73	1.30 ± 0.53	1.71 ± 0.80	3.10 ± 2.00
Distance between centroids (mm)	6 ± 2	5 ± 2	7 ± 5	18 ± 11
Registration time (sec)	227 ± 27	11 ± 2	7 ± 1	199 ± 38
Volume (cc)	3.52 ± 2.00	3.31 ± 2.00	2.83 ± 1.74	2.30 ± 1.64
Difference between volumes (cc) ^a	0.86 ± 0.50	1.10 ± 0.50	1.50 ± 1.00	2.10 ± 1.20

Quantitative measures of image registration accuracy – TG132

Shaaer, A., et al. (2018). *Brachytherapy* 18.1 (2019): 95-102.


Brock, Kristy K., et al. TG132. *Medical physics* 44.7 (2017): e43-e76.

MRI in clinical workflow




CT-Based Technique and Catheter Displacement

Planned



Delivered



Patient Transfer System


- Q-fix Symphony



- Diacor – Zephyr HDR





MRI in clinical workflow




Medical University of Vienna (Low Field - Open Bore)

- Low Field (0.35T MRI)
- Open bore
- Improved field homogeneity
- Low susceptibility artifact
- Reduced chemical shift artefact




Princess Margaret Cancer Centre (MR on rails)

- 1.5T MR on rails
 - MR sim
 - Linac bunker
 - Brachy suite



OCC MR image guided brachytherapy suite




Design Considerations:

- Room Dimensions
- Weight of the MR unit
- Shielding
- Penetrations
- MR and RT safety systems




Shielding for RF and radiation



- **RT safety systems:**
 - Interlocks
 - Emergency Stops
 - Indicators
 - Radiation Monitoring & Source stuck kit
- **MR safety:**
 - Appropriate Zoning
 - Ferromagnetic detectors at the door
 - Quench button
 - Evacuation fans





MRI in clinical workflow



Cost-Utility of MR guided brachy vs. CT or 2D


- Public health single payer perspective
- For all stages MRgBT provides systemic savings over CT and 2D
- Driven by cost of recurrence



Perdrizet, Johnna, et al. *IJROBP* (2020).

2D Brachytherapy at OCC

2D Planning
1 hour




20 Consecutive 2D cervix cases


Total Duration									
1:40:48									
1:26:24									
1:12:00									
0:57:36									
0:43:12									
0:28:48									
0:14:24									
0:00:00									



Start of the OCC program




20 consecutive MR guided cases



Intraoperative Brachytherapy Experience

MR-only Planning
2.5 hours



20 consecutive intraoperative cases

Total Duration								
4:19:12								
3:50:24								
3:21:36								
2:52:48								
2:24:00								
1:55:12								
1:26:24								
0:57:36								
0:28:48								
0:00:00								




MR safe brachytherapy equipment

- **MR safe anesthetic cart**
 - Compatible with Anesthesia supplies
- **MR safe patient monitor**
 - Remote display capability
- **MR Conditional Afterloader**
 - Plastic cables
 - RF shielded afterloader
 - Shielded data cable




Selection of MR Safe Applicators

- **Plastic applicators create signal voids**
 - No local distortions
 - MR markers/Model based reconstruction may aid reproducibility
- **Metal Applicators**
 - Distortion and artefacts need to be quantified
 - Vendor must provide MR safety information




MRI in clinical workflow




Commissioning

- Geometric fidelity checks of MR sequences
- Data transfer integrity
- Source path characterization
- Applicator model validation




Applicator Commissioning

- Index position of first Dwell position
- Distance of dwell to outer surface of applicator
- Spot checks of source path for curved applicators
- Imaging artefacts introduced by metallic applicators




Reconstruction aids



MR sequences for catheter delineation


T2W 3D coarse – (1 mm isotropic)

- Catheters appear as voids
- Negative contrast
- ~3 mins



T1W 3D coarse – (1 mm isotropic)

- Catheters appear as bright
- Positive Contrast
- ~ 3 mins



QA for MR-guided brachy program


Patient specific QA

- Real-time peer review of contours
- Free length verification
- Automated second checks

Programmatic QA

- CPQR QA program
- Monthly QA of MR markers

Channel number	Last Dwell	Tandem (Y/N)	Measured Free Length/mm	Estimated Last dwell	Error
1				-2	2
2				-2	2
3				-2	2
4				-2	2
5				-2	2
6				-2	2
7				-2	2
8				-2	2
9				-2	2
10				-2	2
11				-2	2
12				-2	2
13				-2	2
14				-2	2
15				-2	2
16				-2	2
17				-2	2
18				-2	2
19				-2	2
20				-2	2
21				-2	2
22				-2	2
23				-2	2



HFN: 3357022
 Fraction number: 2
 Plan type: SN Transperineal - WITH tandem
 Plan name: fx2
 Plan date: 20200210
 Plan time: 104257

QA STATUS: SEE WARNINGS IN REPORT

QA REPORT

PASS - Channel length check: All channels have correct indexer length and 1 tandem in channel 5 has a 1300 mm indexer length
 PASS - Plan approval status
 PASS - Flexitron HDR 192-Ir afterloader is being used

WARNING - Treatment position check: Most distal dwell in channel 1 is 21 mm from the channel end

WARNING - Treatment position check: Most distal dwell in channel 4 is 28 mm from the channel end

WARNING - Treatment position check: Most distal dwell in channel 7 is 31 mm from the channel end

WARNING - Treatment position check: Most distal dwell in channel 12 is 51 mm from the channel end!

WARNING - Same treatment position warning as above affects other channels... Please check!

PASS - Channel directionality check: All channels have proper directionality assuming head first scan orientation.

PASS - Channel offset check: All channels have zero offset and 1 tandem in channel 5 has a -6 mm offset

PASS - Transperineal channel free length checks

PASS - Dose calculation based on TG-43 formalism

PASS - Channel mapping check: All channels have correct mapping.

PASS - Activity check

All checks completed. Treatment plan pdf and QA map copied to UN IMPAC HDR



Thank you



Medicine
UNIVERSITY OF TORONTO



Sunnybrook
HEALTH SCIENCES CENTRE

