

# Task Group 182: Electronic Intracavitary Brachytherapy QM based on Risk Analysis

Session #3:

ZEISS INTRABEAM® System and Simple Design of Quality Management Program

presented by

Susha Pillai, M.Sc,M.S,DABR

Medical Physicist

Oregon Health & Science University

1

Disclosures: NONE

2



### Scope of this session:

- Model a risk based Quality Management (QM) program
- Risk analysis model: TG 100 based FMEA tool
- Intracavitary Breast IORT for early stage breast cancer
- ZEISS INTRABEAM® System
- Importance of QC & QA activities for a robust QM program

3

### Outline

- ZEISS INTRABEAM® system
- Intracavitary Breast IORT procedure
- Design risk based QM program




4



## ZEISS INTRABEAM® system

- X-ray probe (XRS source): 50kV, 40μA
- Electron Gun
- Accelerator Section
- Beam Deflector
- Drift tube: 10cm length, 3.2mm diameter
- Gold target: (0.5 μm)
- Internal Radiation Monitor(IRM)





4. INTRABEAM SYSTEM PRS 500 with XRS 4 software version 2.2

5

## ZEISS INTRABEAM® system cont'd....

- Superficial kV X-rays (HVL in mm of Al)
- Nearly spherical dose distribution
- Mean energy in water:  
28keV at the surface of the applicator
- Dose fall off at the rate of  $1/r^3$  in tissue



6. Ebert, M. A., & Carruthers, B. (2003) Medical Physics, 30(9), 2424–243

6



## ZEISS INTRABEAM® system (cont'd)

- Spherical Applicator set
  - Solid, hollow applicators
  - Water equivalent plastic: Polyetherimide material
  - Size: 1.5cm – 5cm in 0.5cm increments



Image courtesy Carl Zeiss Meditec.

7

## ZEISS INTRABEAM® system cont'd....

- Surgical Arm (Floor Stand)
  - 6 degrees of motion
- Treatment Console
- QA devices
  - PAICH (Probe Adjuster & Ion Chamber Holder)
    - XRS Probe straightening tool
    - Output constancy check
  - Soft x-ray ionization chamber (PTW 23342)
  - PDA (Photo Diode Array)
    - Deflection/Isotropy check




courtesy Frank Weigand of Zeiss, Inc.

8




## Intracavitary Breast IORT procedure




9

## QM program for Intracavitary Breast Therapy




10





11



12



## FMEA for "Consultation &amp; Surgery


Process Steps	Failure Modes	Cause	Effect of Failure	AVG (O)	AVG (S)	AVG (D)	RPN= O*S*D
Radiation Oncologist Schedule	RadOnc is not available	Oncologists not alerted	Tx delayed or aborted	6	3.5	1.25	24.75

13

## Process Map:

## Pre-Treatment procedure

- All applicators are available
- Surgical stand check
- Console & XRS probe check
- Isotropy check
- Source output check



14



## FMEA Analysis for “Pre-Tx Procedure”


Processes	Process steps	Failure Modes	Cause	Effects of failure	AVG (O)	AVG (S)	AVG (D)	RPN = O*S*D
Applicator availability	Appl Present	Appl not present	Applicator misplaced	Tx aborted	5.5	3.5	1.5	30.8

15

## Process Map:

## Surgery &amp; IORT procedure

- Applicator selection
- Time calculation
- Treatment delivery



16




## FMEA Analysis for "Surgery & IORT"


Process	Process Steps	Failure Modes	Cause	Effects of failure	AVG (O)	AVG (S)	AVG (D)	RPN = O*S*D
Applicator Selection	Fit lumpectomy cavity	Applicator too small	Surgeon error in fit	Wrong/heterogenous dose	5.8	6.8	5.8	266.3
			Correct size not available	Wrong/heterogenous dose	5.8	6.8	3.8	185

17

*Histogram: RPN*



*Histogram: S values*



# of failure modes = 61

18



# Fault Tree Analysis(FTA)




- Address the high risk (higher RPN value) failure modes
- Graphical representation to identify the cause & effect relationship
  - Use logical diagram (OR gates & AND gates)
- Easy to identify the potential locations of quality measures

19

## FTA cont'd..

- Highest RPN score: Wrong Dose
- Causes ?
  - Applicator size failure
  - Tx Time calculation error
  - Delivery failure




20




## FTA cont'd..

- Highest RPN score: Wrong Dose
  - Applicator size failure ←
  - Tx Time calculation error
  - Delivery failure



21

## Fault Tree :Wrong Dose - Applicator size failure



22




## FTA cont'd..

- Highest RPN score: **Wrong Dose**

- Applicator size failure


- **Tx Time Calculation error** ←

- Delivery Failure



23

## Fault Tree :Wrong Dose-Treatment Time Calculation Failure




24



FTA cont'd..


Highest RPN score: **Wrong Dose**

- Applicator size failure
- Tx Time Calculation error
- Delivery Failure** ←



25

Fault Tree:Wrong Dose - delivery failure



26



## Quality Management System

- Concentrate on common causes
  - Inadequate training
  - Performance failure
- Eliminate common causes and interrupt failure propagation
  - Re-design the procedure
  - Standardize procedures & checklists
  - Adequate staffing & staff training
  - Routine maintenance of software & hardware
  - Clear line of communications

27

## SUMMARY

- Model QM model program & Quality checklist
- FMEA sorting based on RPN vs S values
- FTA should incorporate all required QC/QA measures recommended by the vendor
- FTA helps to identify and incorporate additional QA & QC measures

28



## SUMMARY

- Process Map: Team effort & involvement of multi-disciplinary team members
- Training, facilitator, core team members
- Risk-assessment tools: Task Group 100 & Quality & Safety resources
- Periodic review of QM program:
  - Evaluate the existing QC measures
  - Incorporate institutional ILS

29

Thank you!

30



## Acknowledgements:

### Task Group 182 Members

Bruce R. Thomadsen, PhD, FAAPM(Chair)	Jonathan N. Law, MS
Barrett S. Caldwell, PhD	Mark J. Rivard, PhD, FAAPM, FASTRO FABS
Chris Stacey	Peter J. Biggs, PhD, FAAPM
Habib Safigholi, PhD	Robert A. Cormack, PhD, FAAPM
Gene A. Cardarelli, PhD, FAAPM, FACMP	Sujatha Pai, MS, FAAPM
Habib Safigholi, PhD	Sushakumari Pillai, MS
Henry T. Heaton II, MS, FAAPM	Timothy J. Waldron, MS
James C. H. Chu, PhD, FAAPM, FACR	Tina L Pike, PhD
Jessica R Hiatt, PhD	Wenzheng Feng, MS
Jeffrey P. Limmer, MS, FAAPM	Zoubir Ouhib, MS, FAAPM, FACR, FABS

31

## References:

1. Huq, M. S., Fraass, B. A., Dunscombe, P. B., Gibbons, J. P., Ibbott, G. S., Mundt, A. J., Mutic, S., Palta, J. R., Rath, F., Thomadsen, B. R., Williamson, J. F., & Yorke, E. D. (2016). The report of Task Group 100 of the AAPM: Application of risk analysis methods to radiation therapy quality management. *Medical Physics*, 43(7), 4209–4262.
2. 2013 Summer school session - Failure Mode and Effect Analysis by Frank J. Rath, M.S., University of Wisconsin, Madison, WI
3. AAPM. (2020). *Quality & Safety Resources*. <https://www.aapm.org/QualitySafety/default.as>
4. INTRABEAM SYSTEM Manual : INTRABEAM SYSTEM PRS 500 with XRS 4 software version 2.2
5. Shaikh, M. Y., Nalichowski, A., Joiner, M. C., & Burmeister, J. (2020). Dosimetric evaluation of the INTRABEAM system for breast intraoperative radiotherapy: A single-institution experience. *Medical Dosimetry*, 45(2), e1–e6.
6. Ebert, M. A., & Carruthers, B. (2003). Dosimetric characteristics of a low-kV intra-operative x-ray source: Implications for use in a clinical trial for treatment of low-risk breast cancer. *Medical Physics*, 30(9), 2424–243
7. Clausen, S., Schneider, F., Jahnke, L., Fleckenstein, J., Hesser, J., Glatting, G., & Wenz, F. (2012). A Monte Carlo based source model for dose calculation of endovaginal TARGIT brachytherapy with INTRABEAM and a cylindrical applicator. *Zeitschrift Für Medizinische Physik*, 22(3), 197–204.
8. Biggs, D. S., & Thomson, E. S. (1996). Radiation properties of a miniature X-ray device for radiosurgery. *The British Journal of Radiology*, 69(822), 544–547.

32