

PET/CT QC/QA

Magnus Dahlbom, Ph.D.
Division of Nuclear Medicine
Ahmanson Biological Imaging Center&
Ahmanson Translational Imaging Division
David Geffen School of Medicine at UCLA
Los Angeles



David Geffen
School of Medicine

1

Learning Objectives

- Summarize the signal processing steps for coincidence detection
- Understand the components of a daily QA procedure
- Identify and trouble shoot possible sources of failure in daily QA
- List the recommended frequency of QA/QC tests
- Describe the process of scanner calibration
- Name the different components of the NEMA test used for PET acceptance testing
- Understand the meaning of the results of the NEMA tests.



David Geffen
School of Medicine

2

Outline

- PET Basics
- PET Detectors and data flow
- Daily QC Procedures
- QC Evaluation Techniques
- Examples of Common Problems
- Quarterly and Annual QC Procedures



David Geffen
School of Medicine

3

Quality Control in PET


- Verify the operational integrity of the system
 - Ensure Stability of the Detectors
 - Ensure Stability of the Acquisition Electronics
- Maintain consistent and high image quality
- Minimize chances for image artifacts
- Catch potential problems early
- Maintain quantitative accuracy
- Eliminate unnecessary repeat scans



David Geffen
School of Medicine


4

Coincidence Detection




David Geffen
School of Medicine

Coincidence Detection




David Geffen
School of Medicine

Coincidence Detection



David Geffen
School of Medicine


Data Flow in a PET System



David Geffen
School of Medicine

8

Time-of-flight PET



$$s = v \cdot t$$

$$R + x = vt_1$$

$$R - x = vt_2$$


$$2x = v(t_2 - t_1) \Rightarrow x = \frac{c\Delta t}{2}$$

Time Resolution: ~300 ps



David Geffen
School of Medicine

PET Detectors




David Geffen
School of Medicine

PET Detectors

- To identify the detector elements, X- and Y- positions are calculated for each event:




$$X = \frac{(PMT_A + PMT_B) - (PMT_C + PMT_D)}{PMT_A + PMT_B + PMT_C + PMT_D}$$

$$Y = \frac{(PMT_A + PMT_C) - (PMT_B + PMT_D)}{PMT_A + PMT_B + PMT_C + PMT_D}$$



David Geffen
School of Medicine

Pixelated Detector System




PET Detector Calibration

- Balancing of PMTs
- Identification of Peaks/Detectors in Flood Histogram
- Acquire Energy Spectra / Set Energy Thresholds
- Time Alignment
- Challenge:
 - Large number of detector blocks (~200)
 - Large number of detector elements (>30,000)



David Geffen
School of Medicine

13



Daily QC

- Detect small drifts in the system
 - Detector Gain Drifts (Efficiency, Scatter Fraction, Randoms)
 - Timing
 - Quantification
- Challenges:
 - Large number of detector blocks (~200)
 - Large number of detector elements (>30,000)
 - Time constraint (< 1 hr)




David Geffen
School of Medicine


15

Daily QC Scan

- Transmission/Rotating Rod Sources
- Uniform cylinder phantom (20 cm Ø)
- Also used for calibration (Counts \rightarrow Bq/ml)



Water Phantom
Calibrated w. Dose Calibrator




Solid ^{68}Ge Phantom
Pre-calibrated



David Geffen
School of Medicine

16

Sinogram




- All coincidence lines that are parallel at a given angle form a projection in the sinogram.



David Geffen
School of Medicine

Sinogram




- All coincidence lines (or lines of response) for a given detector form a diagonal trace in the sinogram.





David Geffen
School of Medicine


- Visually inspect the sinograms for:
 - Apparent streaks (hot or cold)
 - Consistency – Comparison to a reference scan




Item	Low Limit	High Limit	Current Reading	N/V/C
[X] Consistency View	13.9255	13.9305	13.9278	Green
[X] Consistency Vary	13.9255	13.9313	13.9278	Green
[X] Straps View	24.945	24.955	24.949	Green
[X] Straps Vary	24.945	24.955	24.949	Green
[X] View Oscillation	-0.1853	0.1843	0.1203	Green
[X] Timing View	-0.05	0.04	0.02	Green
[X] Energy View	-8.0	8.0	0.174095	Green









When things go wrong...




Normal



Detector Failure



Detector Controller Failure




Sorting Memory Failure




David Geffen
School of Medicine

Detector Block Failure





David Geffen
School of Medicine


Detector Block Failure




David Geffen
School of Medicine










System Calibration

- Calculates Conversion Factor : Image A.U. \rightarrow Bq/ml
- Scan uniform 20 cm \varnothing ^{68}Ge or ^{18}F Cylinder
 - Calibrated Activity, Time & Volume
 - Emission & Transmission / CT
- Reconstruct:
 - All corrections applied
 - Standard reconstruction parameters
- Calibration Factor :
Calibrated Act. Conc. / ROI of Image (A.U.)



David Geffen School of Medicine

32

Daily QC Scan / Siemens mCT

Acquire ^{68}Ge Cylinder Scan (100-200Mcts) **Centered!**

Generation of Normalization Factors

Evaluation of Daily QC Components

- Block Noise
- Block Efficiency
- Measured Randoms
- Scanner Efficiency
- Scatter Ratio
- Calibration Factor) & Image Plane Efficiency
- Timing Offset, Width and Time Alignment Fit



David Geffen School of Medicine

33

Daily QC Scan / Siemens mCT

Item	Upper Bound	Lower Bound	Value	Decision
Block Noise	3 [crystal]	0 [crystal]	0 Blocks out of range	Passed
Block Efficiency	120 [%]	80 [%]	0 Blocks out of range	Passed
Measured Randoms	115 [%]	85 [%]	101.9 [%]	Passed
Scanner Efficiency	27.04 [cps/Bq/cc]	14.56 [cps/Bq/cc]	20.5 [cps/Bq/cc]	Passed
Scatter Ratio	36.3 [%]	29.7 [%]	31.7 [%]	Passed
Scanner efficiency correction factor (ECF)	4.1e+007 [Bq*s/ECAT counts]	2.21e+007 [Bq*s/ECAT counts]	3.257e+007 [Bq*s/ECAT counts]	Passed
Image Plane Efficiency	5 [%]	-5 [%]	0 Planes out of range	Passed
Block Timing Offset	0.5 [bin]	0 [bin]	0 Blocks out of range	Passed
Block Timing Width	5 [bin]	0 [bin]	0 Blocks out of range	Passed
Time Alignment Residual	2 [mm]	0 [mm]	1.1 [mm]	Passed
Time Alignment Fit (x / y)	2 [mm]	0 [mm]	0.3 [mm] / 0.3 [mm]	Passed

Scatter Fraction / Energy Threshold




CT Daily QA Scan

- Normal operations include the following 3 tasks (in order):
 - **Tube Warmup**- A built-in prep scan that gradually increases heat loading in the X-ray tube in order to prevent thermal cracking and eliminate the potential for an arc to occur. It includes a series of exposures made at incrementing kVp
 - **Daily Air Cals**- A built-in prep scan that performs a series of exposures at varying techniques in order to normalize the detector response using *air* as the attenuating media. These scans essentially adjust the detector gains to achieve a uniform response
 - **Daily QA Phantom scan**- Provides data for 3 areas of concern in daily quality assurance

© Medtronic MDS&C

Quality Control CT daily regimen



- Scan water layer
 - Measure water HU
 - CTAC
 - Check for artifacts
 - Ring artifacts
 - Redo conditioning
 - Redo Air Cal



Austin Health

Courtesy: Stefan Eberl, RPA

Quarterly QC Procedures

- Quantification Check
- Gantry alignment (for PET/CT)
- Other cross calibrations (well counter, etc)



David Geffen
School of Medicine

39

Quantification Scan

- Scan uniform 20 cm Ø ^{68}Ge or ^{18}F Cylinder

- Emission & Transmission / CT


- Reconstruct:

- All corrections applied
- Standard reconstruction parameters

- Visual inspection


- Compare image ROI activity to calibrated activity

- Always perform after any service or detector adjustments



David Geffen
School of Medicine

40



David Geffen
School of Medicine

Quantification Scan


^{18}F Cylinder
0.0299 $\mu\text{Ci/ml}$

	FBP		OSEM	
Calculated Atten. Corr.	0.0296	-1.1%	0.0295	-1.5%
Meas. Atten. Corr. CT	0.0285	-4.7%	0.0284	-4.0%
Meas. Atten. Corr. Rods	0.0242	-19.2%	0.0239	-20.2%



David Geffen
School of Medicine

PET/CT Gantry Alignment



```
ERROR: undefined
OFFENDING COMMAND: f'~

STACK:
```