

Monte Carlo treatment planning in the clinic - successes and challenges

Part II-electron beams

Joanna E. Cygler, Ph.D.

The Ottawa Hospital Cancer Centre, Ottawa, Canada

Carleton University Dept. of Physics, Ottawa, Canada

University of Ottawa, Dept. of Radiology. Ottawa, Canada



The Ottawa
Hospital | L'Hopital
Cancer Centre d'Ottawa



Objectives – electron beams

- Currently available commercial MC-based treatment planning systems for electron beams.
- Commissioning of such systems in terms of beam models and dose calculation modules.
- Factors associated with MC dose calculation within the patient-specific geometry, such as statistical uncertainties, CT-number to material density assignments, and reporting of dose-to-medium versus dose-to-water.
- Possible clinical impact of MC-based electron beam dose calculations



Rationale for Monte Carlo dose calculation for electron beams

- Difficulties of commercial pencil beam based algorithms
 - Monitor unit calculations for arbitrary SSD values
 - large errors*
 - Dose distribution in inhomogeneous media has large errors for complex geometries

** can be circumvented by entering separate virtual machines for each SSD - labour consuming*



Monte Carlo based Treatment Planning Systems

M C dose calculations give in general the right answer

- There are no significant approximations
 - no approximate scaling of kernels is needed
 - electron transport is fully modelled
 - geometry can be modelled as exactly as we know it
 - all types of heterogeneities can be properly handled
- There are many experimental benchmarks showing M C calculations can be very accurate (see the references)



Components of Monte Carlo based dose calculation system


There are **two** basic components of MC dose calculations, see the next slide:

1. Particle transport through the accelerator head
 - Explicit transport (e.g. BEAM code)
 - Accelerator head model (parameterization of primary and scattered beam components)
2. Dose calculation in the patient



Clinac 2100C

9 MeV electrons



Particle transport through the machine head - beam models

- Direct MC simulation of the accelerator head
 - beam simulations can be done accurately if all the parameters are known - but they often are not
- Beam models provide a solution to the above problem
 - is any algorithm that delivers the location, direction and energy of particles to the patient dose-calculating algorithm.



Example of a beam model

Sub-sources

1 - the main diverging source of electrons and photons;


2 - edge source of electrons;

3 - transmission source of photons;

4 - line source of electrons and photons.

Beam model:
Multiple source model

Dose calculation
in patient



M.K. Fix et al,

Phys. Med. Biol. 58 (2013) 2841-2859




Commercial implementations

- **MDS Nordion** (Nucletron - now Elekta) **2001**
 - First commercial Monte Carlo treatment planning for electron beams
 - Kawrakow's VMC++ Monte Carlo dose calculation algorithm (2000)
 - Handles electron beams from all clinical linacs
- **Varian Eclipse eMC** **2004**
 - Neuenschwander's MMC dose calculation algorithm (1992)
 - Handles electron beams from Varian linacs only (23EX)
 - work in progress to include beam models for linacs from other vendors (*M.K. Fix et al, Phys. Med. Biol. 58 (2013) 2841-2859*)
- **CMS (now Elekta) XiO eMC for electron beams** **2010**
 - Based on VMC (Kawrakow, Fippel, Friedrich, 1996)
 - Handles electron beams from all clinical linacs



Nucletron Electron Monte Carlo Dose Calculation Module



510(k) clearance (June 2002)


- Originally released as part of Theraplan Plus
- Currently sold as part of Oncentra Master Plan
- Fixed applicators with optional, arbitrary inserts, or variable size fields defined by the applicator like DEVA
- Calculates absolute dose per monitor unit (Gy/MU)
- User can change the number of particle histories used in calculation (in terms of particle $\#/cm^2$)
- Data base of 22 materials
- Dose-to-water is calculated in Oncentra
- Dose-to-water or dose-to-medium can be calculated in Theraplan Plus MC DCM
- Nucletron performs beam modeling



Varian Macro Monte Carlo transport model in Eclipse

- An implementation of Local-to-Global (LTG) Monte Carlo:
 - Local: Conventional MC simulations of electron transport performed in well defined local geometries ("kugels" or spheres).
 - Monte Carlo with EGSnrc Code System - PDF for "kugels"
 - 5 sphere sizes (0.5-3.0 mm)
 - 5 materials (air, lung, water, Lucite and solid bone)
 - 30 incident energy values (0.2-25 MeV)
 - PDF table look-up for "kugels"

The above step is performed off-line.




- Global: Particle transport through patient modeled as a series of macroscopic steps, each consisting of one local geometry ("kugel")



Varian Macro Monte Carlo transport model in Eclipse

- Global geometry calculations
 - CT images are pre-processed to user defined calculation grid
 - HU in CT image are converted to mass density
 - The maximum sphere radius and material at the center of each voxel is determined
 - Homogenous areas → large spheres
 - In/near heterogeneous areas → small spheres



Varian Eclipse Monte Carlo

- User can control
 - Total number of particles per simulation
 - Required statistical uncertainty
 - Random number generator seed
 - Calculation voxel size (several sizes available)
 - Isodose smoothing on / off
 - Methods: 2-D Median, 3-D Gaussian
 - Levels: Low, Medium, Strong
- Dose-to-medium is calculated



CMS XiO Monte Carlo system

- XiO eMC module is based on the early VMC* code
 - simulates electron (or photon) transport through voxelized media
- The beam model and electron air scatter functions were developed by CMS
- CMS performs the beam modeling
- The user can specify
 - voxel size
 - dose-to-medium or dose-to-water
 - random seed
 - total number of particle histories per simulation
 - or the goal Mean Relative Statistical Uncertainty (MRSU)
 - minimum value of dose voxel for MRSU specification

*Kawrakow, Fippel, Friedrich, Med. Phys. 23 (1996) 445-457;

*Fippel, Med. Phys. 26 (1999) 1466-1475



User input data for MC based TPS

Treatment unit specifications:

- Position and thickness of jaw collimators and MLC
- For each applicator scraper layer:
 - Thickness
 - Position
 - Shape (perimeter *and edge*)
 - Composition
- For inserts:
 - Thickness
 - Shape
 - Composition

No head geometry details required for Eclipse, since at this time it only works for Varian linac configuration



User input data for MC TPS *cont*

Dosimetric data for beam characterization (beam model), as specified in User Manual, for example:

Beam profiles without applicators:

- in-air profiles for various field sizes
- in-water profiles
 - central axis depth dose for various field sizes
 - some lateral profiles

• Beam profiles with applicators:

- Central axis depth dose and profiles in water
- Absolute dose at the calibration point

Dosimetric data for verification

- Central axis depth doses and profiles for various field sizes



Clinical implementation of MC treatment planning software

- Beam data acquisition and fitting
- Software commissioning tests*
 - Beam model verification
 - Dose profiles and MU calculations in a homogeneous water tank
 - In-patient dose calculations
- Clinical implementation
 - procedures for clinical use
 - possible restrictions
 - staff training

**should include tests specific to Monte Carlo*

A physicist responsible for TPS implementation should have a thorough understanding of how the system works.




Software commissioning tests: goals

- Setting user control parameters in the TPS to achieve optimum results (acceptable statistical noise, accuracy vs. speed of calculations)
 - Number of particle histories
 - Required statistical uncertainty
 - Voxel size
 - Smoothing
- Understand differences between water tank and real patient anatomy based monitor unit values





XiO: 9 MeV - Trachea and spine importance of high quality data



SU-E-T-669

Example of beam model verification

CMS eMC: cutout factors




Vandervoort and Cygler, COMP 56th Annual Scientific Meeting, Ottawa, June 2010




Monte Carlo Settings: Noise in the dose distributions


Varying MRSU, voxel size= $2.5 \times 2.5 \times 2.5 \text{ mm}^3$, dose-to-medium, 6 MeV beam, $10 \times 10 \text{ cm}^2$ applicator



MRSU=10%



MRSU=5%





MRSU=2%




Eclipse eMC

Effect of voxel size and smoothing



Dose-to-water vs. dose-to-medium

Small volume
of water





Voxel of medium

D_m - energy absorbed in a medium voxel divided by the mass of the medium element.


D_w - energy absorbed in a small cavity of water divided by the mass of that cavity.

$$D_w = D_m \left(\frac{S}{\rho} \right)_m^w$$



Dose-to-water vs. Dose-to-medium

Dose-to-water vs. dose-to-medium,
MRSU=2%, voxel size=4×4×4 mm³, 6 MeV
beam, 15×15 cm² applicator, both 602 MU



MU MC vs. hand calculations

Monte Carlo

Real physical dose
calculated on a patient
anatomy

Inhomogeneity
correction included

Arbitrary beam angle

Hand Calculations

Rectangular water
tank


No inhomogeneity
correction

Perpendicular beam
incidence only



9 MeV, full scatter phantom (water tank)

RDR=1 cGy/MU



Point Dose Data


File Edit

	9MeV full scatter
Weight	N/A
MU or min	200
Iso: 9MeV full scatter	N/A
dMax: 9MeV full scatter	N/A
DMax: 9MeV full scatter (0.00,0.00,-2.10)	200.000



Lateral scatter missing

Point Dose Data		
File	Edit	
	9MeV_X2cm	Total Rel. Dose
Weight	N/A	N/A
MU or min	234	N/A
Iso: 9MeV_X2cm	N/A	N/A
dMax: 9MeV_X2cm	N/A	N/A
DMax: 9MeV_X2cm (P) (0.00,0.00,-2.10)	200.000	N/A
Plan Data - Norm: 100.00		
Number of Fractions	1	
MU or min / Fraction	234	



Real contour / Water tank =
 $= 234\text{MU} / 200\text{MU} = 1.17$


Reason for more MU: % isodose at the nominal (reference)
 d_{max} depth < 100%




MU real patient vs. water tank

Patient - 9MEV

Source
RT Eye
Bolus




0.00 cm (# 17733)



$$\text{MC} / \text{Water tank} = 292 / 256 = 1.14$$




Internal mammary nodes





$$\text{MC / Water tank} = 210 / 206 = 1.019$$



MU-real patient vs. water tank: Impact on DVH



Posterior cervical lymph node irradiation - impact on DVH




How long does it take?

- MC gives **entire dose distribution** in the irradiated volume, not just a few points
- time for N beams **is the same as** for 1 beam
- timing is a complex question since it depends on
 - **statistical uncertainty** and how defined
 - **voxel size**
 - **field size**
 - **beam energy** and whether photons or electron
 - speed of **CPU** and **optimization** of compiler
 - complexity of patient specific **beam modifiers**



Monte-Carlo Settings: Effect on computation time



Timing Results XiO TPS:

For 9 and 17 MeV beams, 10x10 cm² applicator and the trachea and spine phantom, timing tests were performed for a clinical XiO Linux workstation, which employs 8 processors, 3 GHz each, with 8.29 GB of RAM.



Summary - electron beams

- Commercial MC based TP systems are available
 - fairly easy to implement and use
 - MC specific testing required
- Fast (minutes) and accurate 3-D dose calculations
- Single virtual machine for all SSDs
- Large impact on clinical practice
 - Accuracy of dose calculation improved
 - More attention to technical issues needed
 - Dose-to-medium is calculated, although some systems calculate dose-to-water as well
 - MU based on real patient anatomy (including contour irregularities and tissue heterogeneities)
- Requirement for well educated physics staff



Acknowledgements

George X. Ding

George Daskalov

Gordon Chan

Robert Zohr

Ekaterina Tchistiakova

Junior Akunzi

Indrin Chetty

Margarida Fragoso

Charlie Ma

Eric Vandervoort

Neelam Tyagi

David W.O. Rogers

I have received research support from Nucletron and Varian.

TOHCC has a research agreement with Elekta.


I hold a research grant from Elekta



Thank You



Rationale for Monte Carlo dose calculation for electron beams



98-10-21

Timing - Nucletron TPS Oncentra 4.0

Anatomy - 201 CT slices
Voxels 3 mm³
10x10 cm² applicator
50k histories/cm²

4 MeV Timer Results:

Init = 0.321443 seconds

Calc = 42.188 seconds

Fini = 0.00158201 seconds

Sum = 42.5111 seconds

20 MeV Timer Results:


Init = 0.311014 seconds

Calc = 110.492 seconds

Fini = 0.00122603 seconds

Sum = 110.805 seconds

System

Manufacturer:	Hewlett-Packard Company		
Model:	HP Z800 Workstation		
Rating:	 Windows Experience Index		
Processor:	Intel(R) Xeon(R) CPU	E5520	@ 2.27GHz 2.26 GHz
Installed memory (RAM):	12.0 GB		
System type:	64-bit Operating System		

Faster than pencil beam!



Timing – Varian Eclipse

Eclipse MMC, Varian single CPU Pentium IV

XEON, 2.4 GHz

10x10 cm², applicator, water phantom,

cubic voxels of 5.0 mm sides

6, 12, 18 MeV electrons,

3, 4, 4 minutes, respectively

Chetty et al.: AAPM Task Group Report No. 105: Monte Carlo-based treatment planning, Med. Phys. 34, 4818-4853, 2007

