



Monte Carlo treatment planning for electron beams

Joanna E. Cygler, PhD, FCCPM, FAAPM

Senior Medical Physicist

Professor, Radiology, University of Ottawa

Adjunct Research Professor, Physics, Carleton University

The Ottawa Hospital Cancer Centre, Ottawa, Canada

**Carleton University Dept. of Physics, Ottawa,
Canada**

**University of Ottawa, Dept. of Radiology. Ottawa,
Canada**



AAPM Spring Meeting, March 15 – 18, 2014, Denver, Colorado


Objectives

- To discuss currently available commercial MC-based treatment planning systems for electron beams.
- To describe commissioning of such systems in terms of beam models and dose calculation modules.
- To discuss the factors associated with MC dose calculation within the patient-specific geometry, such as statistical uncertainties, CT-number to material density assignments, and reporting of dose-to-medium versus dose-to-water.
- To discuss possible clinical impact of MC-based electron beam dose calculations

Rationale for Monte Carlo dose calculation for electron beams

- Difficulties of commercial pencil beam based algorithms
 - Monitor unit calculations for arbitrary SSD
 - values – large errors*
 - Dose distributions in heterogeneous media have large errors for complex geometries

**can be circumvented by entering separate virtual machines for each SSD – labor consuming*



Ding, G. X., et al, *Int. J. Rad. Onc. Biol. Phys.*

(2005) 63:622-633



Components of Monte Carlo based dose calculation system

There are **two** basic components of MC dose calculation, see the next slide:

1. Particle transport through the accelerator head
 - explicit transport (e.g. BEAM code)
 - accelerator head model (parameterization of primary and scattered beam components)
2. Dose calculation in the patient

Clinac 2100C

9 MeV electrons




Example of a beam model

Sub-sources

- 1 - the main diverging source of electrons and photons;
- 2 - edge source of electrons;
- 3 - transmission source of photons;
- 4 - line source of electrons and photons.


M.K. Fix et al, Phys. Med. Biol. 58 (2013) 2841–2859



Commercial implementations

- **MDS Nordion** (Nucletron→Elekta) 2001
 - First commercial Monte Carlo treatment planning for electron beams
 - Kawrakow's VMC++ Monte Carlo dose calculation algorithm (2000)
 - Handles electron beams from all clinical linacs
- **Varian Eclipse eMC** 2004
 - Neuenschwander's MMC dose calculation algorithm (1992)
 - Handles electron beams from Varian linacs only (23EX)
 - work in progress to include beam models for linacs from other vendors (*M.K. Fix et al, Phys. Med. Biol. 58 (2013) 2841–2859*)
- **Elekta-CMS XiO eMC** for electron beams 2010
 - Based on VMC (Kawrakow, Fippel, Friedrich, 1996)
 - Handles electron beams from all clinical linacs

Nucletron Electron Monte Carlo Dose Calculation Module



510(k) clearance (June 2002)


- Originally released as part of Theraplan Plus
- Currently sold as part of Oncentra Master Plan
- Fixed applicators with optional, arbitrary inserts, or variable size fields defined by the applicator like DEVA
- Calculates absolute dose per monitor unit (Gy/MU)
- User can change the number of particle histories used in calculation (in terms of particle $\#/cm^2$)
- Data base of 22 materials
- Dose-to-water is calculated in Oncentra
- Dose-to-water or dose-to-medium can be calculated in Theraplan Plus MC DCM
- Nucletron performs beam modeling

Varian Macro Monte Carlo transport model in Eclipse

- An implementation of Local-to-Global (LTG) Monte Carlo:
 - Local: Conventional MC simulations of electron transport performed in well defined local geometries (“kugels” or spheres).
 - Monte Carlo with EGSnrc Code System - PDF for “kugels”
 - 5 sphere sizes (0.5-3.0 mm)
 - 5 materials (air, lung, water, Lucite and solid bone)
 - 30 incident energy values (0.2-25 MeV)
 - PDF table look-up for “kugels”


This step is performed off-line.

- Global: Particle transport through patient modeled as a series of macroscopic steps, each consisting of one local geometry (“kugel”)



Varian Macro Monte Carlo transport model in Eclipse

- Global geometry calculations
 - CT images are pre-processed to user defined calculation grid
 - HU in CT image are converted to mass density
 - The maximum sphere radius and material at the center of each voxel is determined
 - Homogenous areas → large spheres
 - In/near heterogeneous areas → small spheres



from C. Zankowski et al "Fast Electron Monte Carlo for Eclipse"

Varian Eclipse Monte Carlo

- User can control
 - Total number of particles per simulation
 - Required statistical uncertainty
 - Random number generator seed
 - Calculation voxel size (several sizes available)
 - Isodose smoothing on / off
 - Methods: 2-D Median, 3-D Gaussian
 - Levels: Low, Medium, Strong
- Dose-to-medium is calculated

Elekta - CMS XiO Monte Carlo system

- XiO eMC module is based on the early VMC* code
 - simulates electron (or photon) transport through voxelized media
- The beam model and electron air scatter functions were developed by CMS
- The user can specify
 - voxel size
 - dose-to-medium or dose-to-water
 - random seed
 - total number of particle histories per simulation
 - or the goal Mean Relative Statistical Uncertainty (MRSU)
 - minimum value of dose voxel for MRSU specification
- CMS performs the beam modeling

**Kawrakow, Fippel, Friedrich, Med. Phys. 23 (1996) 445-457;
Fippel, Med. Phys. 26 (1999) 1466–1475

Clinical implementation of MC treatment planning software

- Beam data acquisition and fitting
- Software commissioning tests*
 - Beam model verification
 - Dose profiles and MU calculations in a homogeneous water tank
 - In-patient dose calculations
- Clinical implementation
 - procedures for clinical use
 - possible restrictions
 - staff training

**should include tests specific to Monte Carlo*

A physicist responsible for TPS implementation should have a thorough understanding of how the system works.

User input data for MC based TPS

Treatment unit specifications:

- Position and thickness of jaw collimators and MLC
- For each applicator scraper layer:
 - Thickness
 - Position
 - Shape (perimeter *and edge*)
 - Composition
- For inserts:
 - Thickness
 - Shape
 - Composition

No head geometry details required for Eclipse, since at this time it only works for Varian linac configuration

User input data for MC TPS *cont.*

Dosimetric data for beam characterization (beam model), as specified in User Manual, for example:


- Beam profiles without applicators:
 - in-air profiles for various field sizes
 - in-water profiles
 - central axis depth dose for various field sizes
 - some lateral profiles
- Beam profiles with applicators:
 - Central axis depth dose and profiles in water
 - Absolute dose at the calibration point
- Dosimetric data for verification
 - Central axis depth doses and profiles for various field sizes

Software commissioning tests: goals


- Setting user control parameters in the TPS to achieve optimum results (acceptable statistical noise, accuracy vs. speed of calculations)
 - Number of particle histories
 - Required statistical uncertainty
 - Voxel size
 - Smoothing
- Understand differences between water tank and real patient anatomy based monitor unit values

XiO: 9 MeV - Trachea and spine

importance of high quality data





In high-dose-gradient region sub-millimetre accuracy is required



Example of beam model verification


CMS eMC: cutout factors



Monte Carlo settings: noise in the dose distributions

Effect of Mean Relative Statistical Uncertainty (MRSU):

6 MeV beam, 10x10 cm² applicator, voxel size=2.5 × 2.5 × 2.5 mm³, dose-to-medium





MRSU=10%

MRSU=5%


MRSU=2%

Eclipse eMC

Effect of voxel size and smoothing





Dose-to-water vs. dose-to-medium



D_m - energy absorbed in a medium voxel divided by the mass of the medium element.


D_w - energy absorbed in a small cavity of water divided by the mass of that cavity.

$$D_w = D_m \left(\frac{S}{\rho} \right)_m^w$$



Dose-to-water vs. Dose-to-medium

6 MeV beam, 15x15 cm² applicator, both 602 MU
MRSU=2%, voxel size=4 × 4 × 4 mm³



MU - MC vs. hand calculations

Monte Carlo

Real physical dose
calculated on a patient
anatomy

Heterogeneity correction
included

Arbitrary beam angle

Hand calculations


Rectangular water tank

No inhomogeneity
correction

Perpendicular beam
incidence only

9 MeV, full scatter phantom (water tank)

RDR=1 cGy/MU




Point Dose Data	
File	Edit
	9MeV full scatter
Weight	N/A
MU or min	200
Iso: 9MeV full scatter	N/A
dMax: 9MeV full scatter	N/A
DMax: 9MeV full scatter (0.00,0.00,-2.10)	200.000

100% isodose at the nominal (reference) dmax depth

Lateral scatter missing

Point Dose Data		
File	Edit	
	9MeV_X2cm	Total Rel. Dose
Weight	N/A	N/A
MU or min	234	N/A
Iso: 9MeV_X2cm	N/A	N/A
dMax: 9MeV_X2cm	N/A	N/A
DMax: 9MeV_X2cm (P) (0.00,0.00,-2.10)	200.000	N/A



Plan Data - Norm: 100.00	
Number of Fractions	1
MU or min / Fraction	234



Real contour / Water tank =
 $=234\text{MU} / 200\text{MU}=1.17$


Reason for more MU: % isodose at the nominal (reference) d_{\max} depth is less than 100%

MU real patient vs. water tank




$$\text{MC / Water tank} = 292 / 256 = 1.14$$

Internal mammary nodes





$$\text{MC / Water tank} = 210 / 206 = 1.019$$

MU-real patient vs. water tank: Impact on DVH



Posterior cervical lymph node irradiation - impact on DVH



MU eMC vs. pencil beam


<i>Patient</i>	<i>Smoothing</i>				<i>PB</i>
	<i>None</i>	<i>Low</i>	<i>Medium</i>	<i>Strong</i>	
	<i>eMC</i>				
1 (neck)	120	128	134	138	147
2 (head)	192	198	204	210	224
3 (breast boost)	204	208	215	221	224
4 (breast boost)	192	199	200	202	201
5 (breast boost)	189	192	195	198	199
6 (breast boost)	187	192	196	199	200
7 (breast boost)	176	179	183	188	201
8 (breast boost)	212	219	225	229	222
9 (chest wall)	195	203	204	213	220
10 (melanoma in sternum)	421	429	436	440	448
11 (breast boost)	193	199	205	208	212
12 (internal mammary node)	56	57	59	59	61
13 (breast boost)	186	189	192	194	199
14 (breast boost)	204	206	214	216	224
15 (breast boost)	192	200	207	214	222

How long does it take?

- MC gives entire dose distribution in the irradiated volume, not just a few points
- time for N beams is the same as for 1 beam
- timing is a complex question since it depends on
 - statistical uncertainty and how it is defined
 - voxel size
 - field size
 - beam energy and whether photons or electron
 - speed of CPU and optimization of compiler
 - complexity of patient specific beam modifiers



Monte-Carlo Settings: Effect on computation time



Timing Results XiO TPS:

For 9 and 17 MeV beams, 10x10 cm² applicator and the trachea and spine phantom, timing tests were performed for a clinical XiO Linux workstation, which employs 8 processors, 3 GHz each, with 8.29 GB of RAM.

Timing – Nucletron TPS Oncentra 4.0

Anatomy - 201 CT slices
Voxels 3 mm³
10x10 cm² applicator
50k histories/cm²

4 MeV Timer Results:

Init = 0.321443 seconds

Calc = 42.188 seconds

Fini = 0.00158201 seconds

Sum = 42.5111 seconds

System

Manufacturer:	Hewlett-Packard Company		
Model:	HP Z800 Workstation		
Rating:		6.1 Windows Experience Index	
Processor:	Intel(R) Xeon(R) CPU	E5520	@ 2.27GHz 2.26 GHz
Installed memory (RAM):	12.0 GB		
System type:	64-bit Operating System		

20 MeV Timer Results:

Init = 0.311014 seconds

Calc = 110.492 seconds

Fini = 0.00122603 seconds

Sum = 110.805 seconds

Faster than pencil beam!

Timing – Varian Eclipse

Eclipse MMC, Varian single CPU Pentium IV

XEON, 2.4 GHz

10x10 cm², applicator, water phantom,

cubic voxels of 5.0 mm sides

6, 12, 18 MeV electrons,

3, 4, 4 minutes, respectively

Chetty et al.: AAPM Task Group Report No. 105: Monte Carlo-based treatment planning, Med. Phys. 34, 4818-4853, 2007

Summary

- Commercial MC based TP systems are available
 - fairly easy to implement and use
 - MC specific testing required
- Fast (minutes) and accurate 3-D dose calculations
- Single virtual machine for all SSDs
- Large impact on clinical practice
 - Accuracy of dose calculation improved
 - More attention to technical issues needed
 - Dose-to-medium is calculated, although some systems calculate dose-to-water as well
 - MU based on real patient anatomy (including contour irregularities and tissue heterogeneities)
- Requirement for well educated physics staff

Selected references

1. *Kawrakow, I., M. Fippel, and K. Friedrich. (1996), 3D electron dose calculation using a Voxel based Monte Carlo algorithm (VMC), Med Phys 23 (4):445-57.*
2. *Kawrakow, I., VMC++ electron and photon Monte Carlo calculations optimized for radiation treatment planning, Proceedings of the Monte Carlo 2000 Meeting, (Springer, Berlin, 2001) pp229-236.*
3. *Neuenschwander H and Born E J 1992 A Macro Monte Carlo method for electron beam dose calculations, Phys. Med. Biol. 37 107 - 125.*
4. *Neuenschwander H, Mackie T R and Reckwerdt P J 1995 MMC—a high-performance Monte Carlo code for electron beam treatment planning, Phys. Med. Biol. 40 543-74.*
5. *Janssen, J. J., E. W. Korevaar, L. J. van Battum, P. R. Storchi, and H. Huizenga. (2001). "A model to determine the initial phase-space of a clinical electron beam from measured beam data." Phys Med Biol 46:269-286.*

Selected references cont.

6. *Traneus, E., A. Ahnesjö, M. Åsell. (2001) "Application and Verification of a Coupled Multi-Source Electron Beam Model for Monte Carlo Based Treatment Planning," Radiotherapy and Oncology, 61, Suppl. 1, S102.*
7. *Cygler, J. E., G. M. Daskalov, and G. H. Chan, G.X. Ding. (2004). "Evaluation of the first commercial Monte Carlo dose calculation engine for electron beam treatment planning." Med Phys 31:142-153.*
8. *Ding, G. X., D. M. Duggan, C. W. Coffey, P. Shokrani, and J. E. Cygler. (2006). "First Macro Monte Carlo based commercial dose calculation module for electron beam treatment planning-new issues for clinical consideration." Phys. Med. Biol. 51 (2006) 2781-2799.*
9. *Popple, RA., Weinberg, R., Antolak, J., (2006) "Comprehensive evaluation of a commercial macro Monte Carlo electron dose calculation implementation using a standard verification data set". Med Phys 33:1540-1551.*

Selected references cont.

10. *Faddegon, B.A. and Cygler, J.E., Use of Monte Carlo Method in Accelerator Head Simulation and Modelling for Electron Beams, Integrating New Technologies into Clinic: Monte Carlo and Image-Guided Radiation Therapy, AAPM Monograph No. 32, edited by B.H. Curran, J.M. Balter, I.J. Chetty, Medical Physics Publishing (Madison, WI, 2006) P.51-69.*
11. *Cygler, J.E., Heath, E., Ding, G.X., Seuntjens, J.P., Monte Carlo Systems in Preclinical and Clinical Treatment Planning: Pitfalls and Triumphs, Integrating New Technologies into Clinic: Monte Carlo and Image-Guided Radiation Therapy Monograph No. 32, edited by B.H. Curran, J.M. Balter, I.J. Chetty, Medical Physics Publishing (Madison WI, 2006) p.199-232.*
12. *Chetty, I., Curran, B., Cygler, J.E., et al.,(2007) Report of the AAPM Task Group No. 105: Issues associated with clinical implementation of Monte Carlo-based photon and electron external beam treatment planning. Med Phys 34, 4818-4853.*
13. *Reynaert, N., S. C. van der Marck, D. R. Schaart, et al. 2007. Monte Carlo treatment planning for photon and electron beams . Radiat Phys Chem 76: 643-86.. Radiat Phys Chem 76: 643-86.*

Selected references cont.

14. *Fragoso, M., Pillai, S., Solberg, T.D., Chetty, I., (2008) "Experimental verification and clinical implementation of a commercial Monte Carlo electron beam dose calculation algorithm". Med Phys 35:1028-1038.*
15. *Edimo, P., et al., (2009) Evaluation of a commercial VMC++ Monte Carlo based treatment planning system for electron beams using EGSnrc/BEAMnrc simulations and measurements. Phys Med, 25(3): 111-21.*
16. *Cygler, J.E., and Ding, G.X., "Electrons: Clinical Considerations and Applications " in Monte Carlo Techniques in Radiation Therapy, ISBN-10: 1466507926, Taylor & Francis (CRC Press INC) Boca Raton 2013, p 155-166*
17. *Fix, M. K., Cygler, J. E., Frei, D., Volken, W., Neuenschwander, H., Born, E.J., and Manser, P., (2013), Generalized eMC implementation for Monte Carlo dose calculation of electron beams from different machine types, Phys. Med. Biol. 58, 2841–2859,*
18. *Zhang, A., Wen, N., Nurushev., T., Burmeister, J., Chetty, I., J., (2013), J. .Appl. Clin. Med. Phys., 14, (2), 127-145*

Thank you