

Dose Calculation Algorithms and Commissioning

The Status of Intensity Modulated Proton and Ion Therapy

AAPM Symposium 2014

X. Ronald Zhu, PhD
July 24, 2014



Particles contributing to dose

- Primary protons
 - Elastic interactions with electrons
 - Elastic proton-nucleus scattering

- Secondary particles
 - Non-elastic nuclear interactions
 - Secondary protons and other fragments (deuterons, tritons, alphas, neutrons, etc.)




Dose Algorithms

- Monte Carlo Simulation
 - Not routinely used in the clinics yet
- Analytical calculation – pencil beam algorithms

$$D(x, y, z) = I(d) \times LAT(x, y, d)$$

- $I(d)$ – integral depth dose
- $LAT(x, y, d)$ – lateral dose profile




$$I(d) \propto \int_{-\infty}^{\infty} \int_{-\infty}^{\infty}$$

Effect of low dose envelope

■ Methods of measurements:

• Field size factors:

Sawakuchi *et al.* PMB 2010 - protons
 Inaniwa *et al.* Med Phys 2009 - carbons





Sawakuchi *et al.* PMB 2010

• Concentric circles:



Clasie *et al.* PMB, 2012 - protons

• Concentric square frames:

Pedroni *et al.* PMB 2005 - protons




Effect of low dose envelope



Improved agreement with nuclear interaction term when modeling the integral depth dose


Zhang *et al.* PMB 2011

Effect of low dose envelope



Modified Cauchy-Lorentz function is better choice than Gaussian for lateral profile modeling

Li *et al.* PMB 2012



Spot isodoses for different energies at different depths

Lin *et al.* PMB 2013



Commissioning - Example

Zhu et al. Med Phys 2013

- Total dose: Fluence x Beamlet dose

$$D(x, y, z) = \sum_{E_k} \left\{ \sum_{\text{Beamlet } j} [\Phi_{E_k}(x_j, y_j, z) D_{E_k}^{\text{Beamlet}}(x - x_j, y - y_j, d(z)) \right\}$$

- Beamlet dose: IDD x Kernel

$$D_{E_k}^{\text{Beamlet}}(r, d(z)) = \frac{1}{\rho_{H_2O}} [S(d) \times K(r, d)]$$

$$D_{E_k}^{\text{Beamlet}}(r, d(z)) = \frac{1}{\rho_{H_2O}} [S_{pp}(d) K_{lat,prim}(r, d) + S_{sp}(d) K_{lat,sec}(r, d)]$$

pp - primary photons	$K_{lat,prim}$ - MCS, Moliere theory, 2 Gaussians
sp - secondary particles	$K_{lat,sec}$ - secondary particles, nuclear interaction

Input Data Requirements by the Treatment Planning System

- In air profiles:
 - At 3 to 5 different positions from isocenter (e.g., ±200, ±100, and ±0 mm) for every 10-20 MeV in both directions.
 - If a range shifting device is used, 2~3 complete data sets for 2~3 different thicknesses.



Input Data Requirements by the Treatment Planning System

- Integrated depth doses (IDDs):
 - Depth dose to be measured with a large p-p chamber and suggests

$$R = 3\sigma_{spot} = \sqrt{\sigma_{fluence}^2 + 2(0.0307 \times Range)^2}$$

- IDD is in unit of Gy•mm²/MU.



Fluence Model with Gaussians

- Fluence for individual spot:

$$\phi_{E_k}(x, y; x_m, y_m, z) = \phi_{E_k}^m \sum_i \left[\frac{w_i(E_k)}{2\pi\sigma_i^2(E_k, z)} \exp\left(-\frac{(x-x_m)^2 + (y-y_m)^2}{2\sigma_i^2(E_k, z)}\right) \right]$$

Gaussians


$$w_1(E_k) + w_2(E_k) = 1$$

$$\sigma_i(E_k, z) = \sqrt{\left(\frac{A_i(E_k)}{2} + B_i(E_k)z + \frac{C_i(E_k)}{2}z^2\right)}$$


A, B & C phase space parameters

- Parameters were initially determined fitting input data to analytical formula





- Adjusted based on field size factors




Difference in field size factors




20 mono-energetic fields
 Avg ± Stdev = 0.2% ± 0.7% (-1.7% to 2.1%)




Brain - SFIB



- 11 year old boy
- Glioma of tectal plate
- Single field (Simultaneous) integrated boost (SFIB)



Head & Neck – MFO/IMPT



- 67 yo male
- Squamous cell carcinoma
- Right base of tongue
- CTV66, CTV60 & CTV54
- 3 fields: G280°/C15°, G80°/C345° & G180°/C0°

Acknowledgements

Xiaodong Zhang, PhD
 Falk Poenisch, PhD
 Narayan Sahoo, PhD
 Gabriel Sawakuchi, PhD
 Uwe Titt, PhD
 Yupeng Li, M.S.
 Heng Li, PhD
 Jim, Lii, MS
 Richard Wu, MS
 Kazumichi Suzuki, PhD
 Radhe Mohan, PhD
 Michael Gillin, PhD
 Many others



Thank you!



Basic Information about Bragg Peak Chamber

- Nominal sensitive volume: 10.5 cm³.
- Sensitive volume: r = 40.8 mm, t = 2 mm.
- Nominal response: 325 nC/Gy.
- Reference point 3.5 mm front chamber surface.
- Entrance window: 3.47 mm PMMA.
- WET window: 4 mm.
- $N_{D,wk_p} = (3.181 \pm 0.023) \times 10^6 \text{ Gy/C}^*$
 - Average 3 inter-comparison



*Gillin et al. Med Phys 2010
