




- Respiratory and cardiac motion are inherent problems in medical imaging
- Limits scan quality
 - Resolution
 - Quantification
 - Lesion detectability
 - AC artifacts



2015 AAPM spring clinical meeting







- Respiratory motion considered the resolution limiting factor in thorax imaging
- Future perspective

Daou D. Respiratory motion handling is mandatory to a
 2015 AAPM spring clinical meeting

"Respiratory motion handling is mandatory to accomplish the high-resolution PET destiny"










- State of respiratory gating technology in nuclear medicine:
 - o 10+ years of research
 - Major vendors sell integrated systems
 Many clinics own necessary equipment
- Respiratory gating rarely used in routine imaging
- (my) question: why is respiratory motion correction .


- failing its transition into the clinic?
- (my) answer: cost / benefit
- (my) solution: stick around for the talk!

• 2015 AAPM spring clinical meeting














Respiratory gating in PET

- Hardware driven gating is the field standard
- In recent years several software-based methods have been presented to extract respiratory signal to be used for gating
- Software driven methods appeal
 - o Ease of use
 - Operator independent
 - None of the errors in the application of hardware
 - If integrated properly, their implementation would be a software add-in, and require no change to current clinical protocols


• 2015 AAPM spring clinical meeting

<section-header><section-header><text><text><text><text><text>




• 15


013

































Discussion

- Software gating appeared to work as well as (and in some cases better) than hardware gating

 Limitations not yet seen
- Software gating has obvious advantages:
 Uses existing information that is prematurely thrown out
 Requires no changes to current clinical procedures
- Fits within "doing more with less" framework


 All existing PET scanners are (theoretically) capable of software based gating – require a software patch
- The "low cost" implementation of software gating can reasonably support a PET field where motion corrected images are ubiquitously available for review.


•23


SUV max in images displayed increased an average of
 2015 AAPM spring clinical meeting



 Separating available statistics into phase-bins results in decreased image quality – less statistics per bin










Methods

- Non-linear image morphing has been proposed for recombining gated data
- We present an alternative strategy for utilizing the additional information provided by motion characterization - "gating+"
 - Basic precept: Movement of signal in space is expressed in intensity fluctuations in individual voxels over the gated frames
 - Our methods are based on isolating the fluctuations in voxels, and modulating them according to their reliability










				1		_	1	La								
					X	e	sul	τs								
 Simulat 	ion								Relative co	unt stal	istics					
c							1	4		13		55		200	0	
				Ungal	ied		- 100			-					-	
Lesion/background ratio = 3 Upper disphragm/background = 1.5 Lower disphragm/background = 3.0				Gated				1.0		•					1	
				Gating+			10	*				•			1	
				Motion Map				84		100		-20	8	1	÷¢	
Relative Counts	1			4				13	13		55			2000		
	ungated	gated	gating+	ungated	gated	gating+	ungated	gated	gating+	ungated	gated	gating+	ungated	gated	gati	
Max	76%	174%	76%	58%	112%	87%	58%	94%	90%	55%	79%	79%	55%	76%	75	
Volume (70% max)	29%	4%	29%	154%	4%	29%	171%	46%	71%	188%	96%	92%	183%	100%	96	
	63%	187%	63%	49%	112%	77%	49%	73%	72%	47%	66%	67%	47%	65%	65	
SUV (mean/bckgrnd)	03%															














- Discussion
 Algorithm works with effective signal
 Irrespective of reconstruction algorithm, smoothing, etc
 Irrespective of quality of signal
 - Areas not benefiting from gating, or entire scans not benefiting from gating, will be returned to their ungated embodiment


• 34

0.35

- Algorithm utilizes available information and optimizes its transformation into Cartesian space.
 Does not preclude the use of non-linear morphing algorithms
- Potential applications:
 - o Support use of routine gating thorax imaging
 - Respiratory, Cardiac imaging
 Human, small animal
 - PET, SPECT, CT (low dose 4D CT), MRI...

• 2015 AAPM spring clinical meeting

2015 AAPM spring clinical meeting



- When information is optimized in frequency space, during the gating+ processing, there is an opportunity to shift the phase of the signal by rotating the frequency components in real and imaginary space. This allows a user to extract a voxel value at any and all phases of the cycle.
 values adhere to the optimized frequency information
- By repeating process for all voxels, can reconstruct phase shifted images

 -0.02 seconds processing per slice
- With this process, we can reconstruct continuous motion image (CMI) sequences

• 37

• 2015 AAPM spring clinical meeting


<section-header><text><list-item><list-item><list-item><list-item>




Workflow results

- All images created using standard FDG PET acquisitions
- Animations created with 90 frames/cycle, displayed with 30 frames/second


 2015 AAPM spring clinical meeting 	• 41



































• Acknowledgements

Current collaborators Phillip Koo (UC Denver) Kate Petrova (UC Denver) Jonathan Chung (NJH) Vicki LaRue (NJH)

<u>Siemens support</u> Jim Hamill Judson Jones

● 49

Mentors Daniel Silverman (UCLA) Stig Palm (IAEA) Nanette Freedman (Hebrew University)

<u>Co-(data driven) pioneers</u> Paul Schleyer (Kings College) Florian Büther (University of Münster)

• 2015 AAPM spring clinical meeting