NIST Air-Kerma Standard for Electronic Brachytherapy Calibrations

Michael G. Mitch, Ph.D.

Leader, Dosimetry Group

Radiation Physics Division

Physical Measurement Laboratory

National Institute of Standards and Technology



Disclaimers

Certain commercial equipment, manufacturers, instruments, or materials are identified in this presentation in order to specify the experimental procedure adequately. Such identification is for informational purposes only and is not intended to imply recommendation or endorsement by the National Institute of Standards and Technology, nor is it intended to imply that the manufacturer, materials, or equipment are necessarily the best available for the purpose.

Xoft, Inc. provided funding for the development of the NIST electronic brachytherapy facility and supplied systems for source control.

NIST Dosimetry Group Strategic Element

Develop dosimetric standards for x rays, gamma rays, and electrons

based on the SI unit, the gray, $1 \text{ Gy} \equiv 1 \text{ J} / \text{kg}$





free-air chamber



calorimetry





cavity chamber

ultrasonic/optical

NIST Standards for Radiation Therapy



- External beam (⁶⁰Co, orthovoltage and MV x rays, electrons, protons)
- Brachytherapy

Low-Energy, Low-Dose-Rate (¹²⁵I, ¹⁰³Pd, ¹³¹Cs seeds) High-Energy, Low-Dose-Rate (¹⁹²Ir seeds, ¹³⁷Cs sources) *High-Energy, High-Dose-Rate (¹⁹²Ir sources)* Low-Energy, High-Dose-Rate (miniature x-ray sources)

NIST Standards for Radiation Therapy



Safety and efficacy requires accurate treatment planning

Dosimetry traceable to primary standards

Dosimetry of X Rays (*E* < 300 keV)



$$K = \frac{\mathrm{d}E_{\mathrm{tr}}}{\mathrm{d}m} = \boldsymbol{\Phi} \cdot \boldsymbol{E} \cdot \left(\frac{\mu_{\mathrm{tr}}}{\rho}\right)$$

KERMA = <u>Kinetic Energy Released per unit MA</u>ss

transferred to electrons by x rays

$$\frac{\mu_{\rm tr}}{\rho} = \frac{f_{\rm pe}\sigma_{\rm pe} + f_{\rm incoh}\sigma_{\rm incoh}}{\int uA}$$
photoelectric Compton

Photon and Charged-Particle Data Center



XCOM: Photon Cross Sections Database http://www.nist.gov/pml/data/xcom/index.cfm

Photon Cross Sections

Bibliography

http://www.nist.gov/pml/data/photon_cs/index.cfm



http://physics.nist.gov/PhysRefData/Star/Text/ESTAR.html

Dosimetry of X Rays (*E* < 300 keV)



$$K_{\rm air} = Q_{\rm air} \left(\frac{\overline{W}_{\rm air}}{e}\right) \frac{1}{\rho_{\rm air} V}$$

KERMA = <u>Kinetic Energy Released per unit MA</u>ss

$$\mathbf{C} \cdot \left[33.97 \, \frac{\mathbf{J}}{\mathbf{C}} \right] \cdot \frac{1}{\mathbf{kg}} = \frac{\mathbf{J}}{\mathbf{kg}} = \mathbf{G}\mathbf{y}$$

Air kerma can be measured absolutely with a free-air ionization chamber

Free-Air Ionization Chamber (*E* < 300 keV)







NIST Free-Air Chambers





Lamperti

Ritz

Chamber	X-ray tube potential (kV)	Plate separation (mm)	Plate height (mm)	Collector length (mm)	Aperture diameter (mm)	Air absorption length (mm)	Electric field strength (V / cm)
Lamperti	10 to 60	40	50	10	5	39	750
Ritz	20 to 100	90	90	70	10	127	55





Maze entry (leaded glass)



Control area

- The Xoft x-ray source can not be continuously rotated (like a brachytherapy seed)
- Lamperti free-air chamber and HPGe spectrometer **rotate around the source**





source in water cooling catheter



leaded glass shield



Comparison of Lamperti and Ritz Free-Air Chambers



NIST Electronic Brachytherapy Calibration Facility, v. 1 PROBLEM – Alignment not reproducible



NIST Electronic Brachytherapy Calibration Facility, v. 2 SOLUTION – Optical table for rigid mounting of instruments



NIST Electronic Brachytherapy Calibration Facility, v. 2 SOLUTION – Larger lead-glass surround







Pulse Height Distribution - Xoft Source at 50 kV



Fluorescence peaks at 14.9 keV and 16.7 keV are from Y

Peaks from 8 keV to 12 keV are from the W anode

Spectrometry of X-Ray Sources

For a photon detector, the measured pulse-height distribution, H(h), is given by

$$H(h) = \int S(E)R(E,h) dE$$

S(E) is the incident photon spectrum

R(*E*,*h*), the *response function*, is the probability per pulse height that a photon incident with energy *E* will produce a pulse of height *h*

Spectrometry of X-Ray Sources

The response function can be written

$$R(E,h) = T(E) \int D(E,\varepsilon) G(\varepsilon,h) d\varepsilon$$

T(*E*) is the *window-attenuation factor*

 $D(E,\varepsilon)$, the energy-deposition spectrum, is the probability per deposited energy that a photon incident with energy E deposits an energy ε in the detector

 $G(\varepsilon,h)$, the *intrinsic resolution function*, is the probability per pulse height that the deposition of energy ε will give rise to a pulse of height h

Spectrometry of X-Ray Sources

The *energy-deposition spectrum* $D(E,\varepsilon)$ depends on the detector dimensions: cylinder of radius *r* and height *z*

For *E* < 300 keV:

 $D(E,\varepsilon) = P_0(E,\varepsilon) \,\,\delta(\varepsilon - E)$

Photopeak (complete absorption)

+ $P_{x\alpha}(E,\varepsilon) \,\delta(\varepsilon - E + E_{\alpha}) + P_{x\beta}(E,\varepsilon) \,\delta(\varepsilon - E + E_{\beta})$

Ge K_{α} and K_{β} fluorescence x-ray escape

 $+ C(E,\varepsilon)$

Compton continuum

Accurately calculated by Monte Carlo

Seltzer, S.M., "Calculated response of intrinsic germanium detectors to narrow x-ray beams with energies up to 300 keV," *Nucl. Instr. Meth.* **188**, 133-151 (1981).



Unfolded Spectrum: Xoft source at 50 kV

 $H(h) = \int S(E) R(E,h) dE$



Spectrum of Xoft Source at 50 kV



Free-Air Chamber Correction Factors for Xoft Source at 50 kV

$$\dot{K}_{air} = I_{air} \left(\frac{\overline{W}_{air}}{e} \right) \frac{1}{\rho_{air} V} \prod_{i} k_{i}$$

Air-kerma rate at 50 cm

Fa	octor	For:	Lamperti	Ritz
1	k _{ion}	ion recombination	≈1.0000	≈1.0000
2	k _{humidity}	humidity of air	0.998	0.998
3	k _{att}	attenuation	1.0087	1.0283
4	k _{el}	electron loss	1.0008	1.0000
5	k _{sc}	photon scatter	0.9987	0.9970
6	k _{fl}	fluorescence reabsorption	0.9979	0.9969
7	k _{br} /(1-g)	effects of bremsstrahlung	1.0	1.0
8	k _{ii}	initial ion	1.0	1.0
9	k _{dia}	diaphragm scatter	1.0	1.0
П	k ₁₋₉		1.0041	1.0200

Uncertainty Budget for Xoft Source at 50 kV

		Relative standard	
Component	For:	Туре А	Туре В
	net charge or current	s _Q a, s _l a	0.06
Q _{net} , I _{net}	typical value	0.14 ^b	
W/e	mean energy per ion pair	-	0.15
ρ ₀	air density	0.01	0.07
V _{eff}	effective volume	0.04	0.01
k _{ion}	ion recombination	0.03	
k _{humidity}	humidity of air		0.04
k _{att}	attenuation	-	0.11
k _{el}	electron loss	-	0.06
k _{sc}	photon scatter	-	0.03
k _{fl}	fluorescence reabsorption	-	0.05
k _{br} /(1-g)	effects of bremsstrahlung	-	0.02
k _{ii}	initial ion	-	0.04
k _{dia}	diaphragm scatter	-	0.10
k _d	electric field distortion	-	0.20
	polarity difference	0.02	
Combined	air kerma	0.15	0.321

^a Determined as the standard deviation of the mean of the measurement.

^b Typical value for sources measured in 2013/2014

U = 0.71 % (k = 2)

Air-Kerma Rate vs. Air-Kerma Strength



Air-kerma strength

Factor		Lamperti
П k ₁₋₉		1.0041
K _{vac} /K _{air}	conversion to air-kerma strength	1.12

Measurement Traceability for Brachytherapy Sources



AAPM TG-43 Formalism for Brachytherapy Dose Calculations



Dose rate in water

 $\dot{D}(r,\theta) = S_K \cdot \Lambda \cdot \frac{G_L(r,\theta)}{G_L(r_0,\theta_0)} \cdot g_L(r) \cdot F(r,\theta)$

Geometry Function

$$G_L(r,\theta) = \frac{\beta}{Lr\sin\theta} \quad G_L(r,0) = (r^2 - L^2/4)^{-1}$$

Dose Rate Constant (NIST-traceable S_{κ})

Radial Dose Function

2D Anisotropy Function

$$\Lambda = \frac{\dot{D}(r_0, \theta_0)}{S_K} \qquad \begin{array}{c} r_0 = 1 \text{ cm} \\ \theta_0 = \pi / 2 \end{array} \qquad g_X(r) = \frac{\dot{D}(r, \theta_0)}{\dot{D}(r_0, \theta_0)} \frac{G_X(r_0, \theta_0)}{G_X(r, \theta_0)} \quad F(r, \theta) = \frac{\dot{D}(r, \theta)}{\dot{D}(r, \theta_0)} \frac{G_L(r, \theta_0)}{G_L(r, \theta)} \frac{G_L(r, \theta_0)}{G_L(r, \theta)}$$

Modified Formalism for Electronic Brachytherapy Sources

Dose-rate conversion coefficient χ

DeWerd, Culberson, Micka, and Simiele, Brachytherapy (2015)

- TG-43 point-source approximation
- 2D Anisotropy Function applicable due to polar anisotropy
- *i* subscript denotes applicator

Dose rate in water

$$\dot{D}_i(r,\theta) = \dot{K}_{50cm} \cdot \chi_i(r_0,\theta_0) \cdot G_P(r,\theta) \cdot g_i(r) \cdot F_i(r,\theta)$$

Dose Rate Conversion Coefficient (NIST-traceable
$$\dot{K}_{50cm}$$
)

Radial Dose Function

2D Anisotropy Function

$$\chi_{i} = \frac{\dot{D}_{i}(r_{0},\theta_{0})}{\dot{K}_{50cm}} \qquad \begin{array}{c} r_{0} = 1 \text{ cm} \\ \theta_{0} = \pi / 2 \end{array} \qquad g_{i}(r) = \frac{\dot{D}_{i}(r,\theta_{0})}{\dot{D}_{i}(r_{0},\theta_{0})} \frac{G_{P}(r_{0},\theta_{0})}{G_{P}(r,\theta_{0})} \qquad F$$

$$F_i(r,\theta) = \frac{\dot{D}_i(r,\theta)}{\dot{D}_i(r,\theta_0)}$$

Geometry Function

$$G_P(r,\theta) = \frac{1}{r^2}$$

AAPM Dosimetric Prerequisites

LDR Brachytherapy

- Air-kerma strength calibrations traceable to NIST
- TG-43 parameters published (experimental and Monte Carlo)
- NIST standard transferred to the ADCLs
- Annual comparisons between NIST and ADCLs

AAPM Dosimetric Prerequisites

LDR Brachytherapy

- Air-kerma strength calibrations traceable to NIST
- TG-43 parameters published (experimental and Monte Carlo)
- NIST standard transferred to the ADCLs
- Annual comparisons between NIST and ADCLs

Electronic Brachytherapy

- AAPM Task Group proposed

Summary

- NIST air-kerma standard for electronic brachytherapy realized
- Standard transferred to AAPM ADCL using a well chamber
- Proficiency test with AAPM ADCL completed
- New calibration service pending: "Well Ionization Chamber Calibration with Electronic Brachytherapy Sources"
- Clinical implementation of new standard in progress

Acknowledgements

- Xoft, Inc. funding for the development of the NIST electronic brachytherapy facility; sources and associated equipment
- Mel McClelland and Dave Eardley design, fabrication, and installation of mechanical, electrical, and electronic systems
- Ron Tosh computer control code
- Jason Walia spectrometer calibration