

Beam steering and monitoring photon and electron beam energies with ion chamber array

Song Gao, Ph.D

Department of Radiation Physics

MD Anderson Cancer Center



Disclosure

Receives research funding from Sun Nuclear Corporation.



Learning Objectives

- 1. To understand the use of Ion Chamber Array (ICA) to steer beams with the real-time feedback.**
- 2. Monitoring of photon beam energy with ICA .**
- 3. Using ICA combine with double wedge phantom measuring electron beam energy.**



Beam Steering using the IC PROFILER

Device Setup:

100 cm source to the surface of device distance

Photons:

Field: 30 x 30 cm² (maximum size on ICP)

Effective Depth: 2 cm solid water + 0.9 cm intrinsic (eliminates many of the contamination electrons)

Electrons:

Cone: 25 x 25 cm²


Effective Depth: Energy dependent (1-3 cm)




Beam Steering using IC PROFILER

Software Setup:

- CAX Point difference (CPD) Symmetry metric
- Instant Rate mode
- Invert: Real time



Real time: profile, mirror profile, the deviation between profile and its mirror




Beam Steering using IC PROFILER

Ensure machine was steered and working correctly.

Servos:
zeroed,
enabled.

Symmetry:
 $\leq 0.5\%$
is
achievable.



Evaluation of beam steering with ion chamber array

After beams have been steered with ICA


Comparison of Profiles: ICA vs. 3D water scan:

- a) Compare the overlaid profile shapes
- b) calculation of point-by-point differences.
- c) examining the symmetry differences in Varian's “**Point Difference Symmetry**” metric



Profiles: ICP vs. 3D SCANNER

- a) Qualitative compare profiles: ICA vs 3DS
- b) point-by-point symmetry differences: ICA vs. 3DS




Profile Comparison Results

c) Difference of point symmetry between ICA and 3DS:

- 90 photon and 160 electron beams
- All beams were steered with the ICA.

95% beam \leq 0.7%; 100% beam \leq 1.0%



Beam Steering using IC PROFILER

- **Solid water was added during steering to ensure measurements were beyond the buildup region.**
- **Beam was steered in the radial and transverse directions with servos open.**
- **Make sure using “Instant Rate mode” and “CAX Point Difference Symmetry” metric.**
- **With steering complete, servos were zeroed and enabled. Symmetry $\leq 0.5\%$ is achievable.**



Beam Steering using the IC PROFILER

- **2D Ion Chamber Array can be used to steer linear accelerator photon and electron beams**
- **IC Array is able to achieve beam symmetry matches that of 3D water scanning system.**



Monitoring beam energy constancy with IC Array



- 1. Monitoring changes in photon energy (off-axis ratio metric) with the ICA**
- 2. Monitoring changes in electron energy (R_{50}) with the combination of ICA and an aluminum double-wedge plate**



Monitoring photon beam energy with ICA


Changes in photon energy (change bending magnet current)

- Percent depth dose (PDD)
- **Off-axis ratio (OAR) from beam profile**



Off-axis ratio vs. PDD

PDD vs. Off-axis ratio (OAR)



Changes in OAR is ~ 3 times than changes in PDD



Monitoring photon beam energy with ICA

PDD vs. OAR

Minimum detectable change in energy (%)

- PDD metric: $\pm 2.5\%$ changes in energy
- OAR metric: $\pm 0.5\%$ changes in energy

OAR-based metric is more sensitive to energy change than PDD metric




Monitoring photon beam energy with ICA

Long term Results:

- Monthly of 4 years
- 260 measurements
- Differences from the baseline:



$98.4\% \leq \pm 0.5\%$;

$100\% \leq \pm 0.7\%$



Electron beam energy constancy

IC Profiler/double wedge plate combination





Gao, et al., Monitoring linear accelerators electron beam energy constancy with a 2D ionization chamber array and double-wedge phantom, JACMP (2020) 21(1): 18-25.





Electron beam energy constancy

- In DW Profile: D_{50} is directly related to the amount of Al in the beam, related to energy
- (R_{50} FWHM) of the diagonal profile as its energy metric



Electron beam energy constancy

- DW profiles:
25x25 cm² cone at 100cm SSD with DW Plate
- R₅₀ values of known energy beams



Electron beam energy constancy


Calibration:

Establish a linear relationship:

R_{50} vs. FWHM (x)

Measurement:

- R_{50} can be determined from FWHM of the DW profile.
- IC Profiler software reports R_{50} .



Evaluate R_{50} measured with ICP/DWP

R_{50} : ICP/DWP vs. in water scans

For 38 beams of known and unknown energies for the same **Type** of linac as the one for calibration.

R_{50} measured with ICP/DWP agreed with R_{50} in Water:

- mean $\pm \sigma$: 0.1 ± 0.3 mm
- maximum discrepancy: 0.7 mm.



Electron beam energy constancy

R_{50} : Long-term reproducibility with ICP/WD


5 electron beams monthly over a 3 year

R_{50} : Baseline vs. ICP/DW:

99.4% $\leq \pm 0.5$ mm

100% $\leq \pm 1.0$ mm

185 measurements



Energy	6e	9e	12e	16e	20e
$\delta R_{50} \pm \sigma$ (cm)	0.040 \pm 0.008	0.020 \pm 0.011	0.010 \pm 0.011	0.020 \pm 0.010	0.020 \pm 0.013



Electron beam energy


For different Type of Linacs: **Individual vs. universal calibration**

Individual Calibration:

- FWHM, DW profile normalized to CAX.
- Ignore off-axis effect


Universal Calibration:

- FWHM, DW profile normalized to open field profile
- Taking into account off-axis corrections



Electron beam energy


Individual calibrations vs. Universal calibration



Electron beam energy measurement

R_{50} : Individual vs. universal calibration

Total:



$$\text{Maximum } \delta R_{50_max} = R_{50} \text{FWHM} - R_{50} \text{Water}$$

2.5 mm vs. 0.6 mm



Conclusions

- 1. Use of the ICA greatly speeds up the beam steering process because of its the real-time feedback and reduces effort by eliminating the need to setup a 3D water scanning tank.**
- 2. Monitoring photon and electron beam energy with IC Profiler could be incorporated in routine QA procedures saving time while maintaining quality.**



Acknowledgements

- **Peter Balter, Ph.D**
- **Mark Rose, Ph.D**
- **Bill Simon, M.S**
- **Mikhail Chetvertkov, Ph.D**
- **Amir Sadeghi, Ph.D**

