

ARTIFICIAL INTELLIGENCE IN LUNG CANCER: USING RADIOMICS TO PREDICT TUMOR RECURRENCE

Sarah A. Mattonen, Ph.D.


Assistant Professor, Depts. of Medical Biophysics, Oncology, and Biomedical Engineering

Baines Research Chair in Translational Cancer Imaging

Western University, London, Ontario, Canada

Non-Small Cell Lung Cancer

- ▶ Survival rates remain quite poor despite advances in diagnosis and treatment.
- ▶ 5-year survival rates:
 - > Stage I: 55-75%
 - > Stage II: 40-50%
 - > Stage III: 5-35%
 - > Stage IV: <5%




Lang-Lazdunski, L. et al. Eur Resp Soc. (2013) 382-404.
SEER Stat Fact Sheets: Lung and Bronchus Cancer.

The Role of Imaging

Computed Tomography (CT)




Positron Emission Tomography (PET)




The Role of Imaging

Computed Tomography (CT)




Positron Emission Tomography (PET)



Radiomics

- Radiomics aims to extract more complex quantitative information (e.g., texture) from standard medical images.



Positron Emission Tomography (PET)


- ▶ SUVmax has been shown to predict a higher risk of recurrence or death in NSCLC.



Travis WD, et al. J of Thor Onc. (2011) 6:244-285.
Kadota K, et al. J of Thor Onc. (2015) 10:806-814.
Lee et al., European Radiology (2017) 27:1912-1921.

Positron Emission Tomography (PET)

- ▶ Tumour invasion from the main tumour mass.
- ▶ Dissemination of disease throughout the body.




Travis WD, et al. J of Thor Onc. (2011) 6:244-285.
Kadota K, et al. J of Thor Onc. (2015) 10:806-814.
Lee et al., European Radiology (2017) 27:1912-1921.

Objective

To develop a software system integrating PET imaging and non-imaging biomarkers to improve lung cancer prognosis and risk stratification.

Computer-Aided Risk Stratification



Materials

► **Training Cohort (n = 145):**


- > Selected from two local medical centers.
- > All patients had pre-operative PET/CT performed prior to surgery.
 - Feature selection and model development

► **Testing Cohort (n = 146):**


- > Selected from three local medical centers.
- > Underwent PET/CT imaging prior to definitive treatment as part of observational biomarker study.
 - Model evaluation

Segmentation: Tumour

- ▶ The **metabolic tumour volume (MTV)** was segmented on the PET image.
- ▶ A 3-dimensional **penumbra** region was also generated surrounding the MTV to sample surrounding uptake.
- ▶ Three regions were evaluated:
 - > MTV only
 - > Penumbra only (excluding the MTV)
 - > MTV plus penumbra




Methods: Bone Marrow Segmentation



MIM (MIM Software Inc., Cleveland, OH).

Mattonen, et al. Radiology: 293(2), 451-459, 2019.

Radiomic Feature Extraction




- A total of **668 radiomics features** were extracted from the volumes of interest.

Echegaray, S., et al. (2017) J Digit Imaging, 1-12 .

GitHub: riip/3d_qifp

Methods: Model Training


- ▶ Top predictive features were selected using randomizations of 4-fold cross-validation of LASSO Cox regression.



Clinical


Methods: Model Training

- Top predictive features were selected using randomizations of 4-fold cross-validation of LASSO Cox regression.



Clinical


+



Tumour Plus Penumbra


Methods: Model Training

- Top predictive features were selected using randomizations of 4-fold cross-validation of LASSO Cox regression.




Clinical

+



Tumour Plus Penumbra

+





Bone Marrow

Results: Multivariate Model

Clinical + Tumour + Bone Marrow


Feature Type	Feature	HR [95% CI]	p-value
Clinical	Stage	1.98 [1.45-2.70]	p<0.001*
Blood	WBC (1000/uL)	0.99 [0.88-1.11]	p=0.81
	Hemoglobin (g/dL)	0.99 [0.82-1.20]	p=0.93
	Platelets (1000/uL)	1.00 [1.00-1.01]	p=0.93
Tumor	MTV Plus Penumbra GLCM Energy (MAD)	0.69 [0.40-1.19]	p=0.18
	Penumbra GLCM Entropy (IQR)	1.35 [0.97-1.86]	p=0.07
	Penumbra GLCM Cluster Shade (Max)	1.17 [0.84-1.63]	p=0.36
Bone Marrow	GLCM Sum Mean (Skewness)	0.52 [0.32-0.84]	p=0.008*
	GLCM Cluster Tendency (Skewness)	1.62 [1.02-2.59]	p=0.04*

Results: Risk Stratification




Qualitative Results

(A)



$SUV_{max} = 10.3$
Stage I
High-Risk Radiomics
Recurrence


(B)



$SUV_{max} = 10.1$
Stage II
Low-Risk Radiomics
No Recurrence

Qualitative Results

(A)




Stage I
High-Risk Radiomics
Recurrence

(B)




Stage III
Low-Risk Radiomics
No Recurrence


Results: Summary




Concordance = 0.69 [0.60-0.77]




+



Concordance = 0.75 [0.67-0.82]



+



+



Concordance = 0.78 [0.70-0.85]


CAN ADDING CT FEATURES IMPROVE PERFORMANCE?

Jaryd Christie, CAMPEP PhD Candidate



Tumor and Peri-tumoral CT Segmentation

- ▶ MATLAB based-GUI for semi-automatic tumor segmentation on CT




Christie, Jaryd R., et al. Vol. 12036. SPIE, 2022.


<https://github.com/baines-imaging-mattonen-lab/CT-Lung-Tumour-Segmentation>

Qualitative Features

- ▶ Tumour features that describe the location and geometry
- ▶ Features which characterize the lung tissue, bronchi, and lumen



Spiculated



Severe Emphysema

Results: Feature Selection


- ▶ **Seven** selected features:
 - > **One** clinical feature
 - **Stage**
 - > **Six** radiomic features (3 texture, 3 first-order)
 - **Three CT (2 Tumour, 1 Peritumoural)**
 - **Three PET (2 Peritumoural, 1 Bone Marrow)**

Results: Model Evaluation

- ▶ Training: **Stage** vs **Radiomics + Stage**
 - > **Concordance: 0.67** [95% CI: 0.58 – 0.76] vs **0.78** [95% CI: 0.70-0.86]
 - > **p < 0.005**
- ▶ Testing: **Stage** vs **Radiomics + Stage**
 - > **Concordance: 0.60** [95% CI: 0.48 – 0.74] vs **0.76** [95% CI: 0.59-0.87]
 - > **p = 0.008**
- ▶ Radiomics model significantly stratified patients into high- and low-risk of recurrence


Results: Risk Stratification

Training, n=101



Concordance = 0.78

Testing, n=34



Concordance = 0.76

Conclusions

- ▶ These radiomics based tools have the potential to identify NSCLC patients at a higher risk of recurrence and may add clinical utility for risk stratification.
- ▶ This assist physicians in distinguishing patients who may benefit from adjuvant or more aggressive personalized treatment options.

Next Steps

- ▶ Collaborations for external validation of models.
- ▶ Implementation of standardized radiomics features and open-source software


Translational Cancer Imaging



Computer-Aided Decision Support



Improve Patient Outcomes



AFFILIATIONS



 **BAINES
IMAGING**



FUNDING AGENCIES



Acknowledgements



Department of Oncology, Western University

Pencilla Lang, MD, PhD
Andrew Arifin, MD
David Palma, MD, PhD
Jon Snir, PhD
Kathleen Surry, PhD

Department of Medical Imaging, Western University

Mohamed Abdelrazek, MD
Narinder Paul, MD

Division of Thoracic Surgery, Western University

Richard Malthaner, MD
Mehdi Qiabi, MD
Rahul Nayak, MD
Deb Lewis

Department of Pathology, Western University

Matt Cecchini, MD, PhD

Department of Radiology, Stanford University

Sandy Napel, PhD
Kristen Yeom, MD
Guido Davidzon, MD, MSc
Ann Leung, MD
Daniel Rubin, MD, MS
Minal Vasanawala, MD

Fred Hutchinson Cancer Center, University of Washington

Vish Nair, MD
Perrin Romine, MD



@MattonenLab

<http://www.bainesimaging.com/>

Translational Cancer Imaging



Computer-Aided Decision Support



Improve Patient Outcomes

